

## Cleaning FlexiPort™ Blood Pressure Cuffs

EQUIPMENT	PREPARATION	CLEANING
<p>FlexiPort™ Reusable Cuffs, tubing and port fittings</p>	<ul style="list-style-type: none"> <li>Remove FlexiPort™ fitting and tubes from cuff, seal open port with accessory plug. Close the hook and loop.</li> </ul>	<p>Use one or more of the following methods and allow to air dry:</p> <ul style="list-style-type: none"> <li>Wipe with mild detergent and water solution (1:9 solution). Rinse.</li> <li>Wipe with Enzol® per manufacturer's instructions. Rinse.</li> <li>Wipe with 0.5% bleach and water solution. Rinse.</li> <li>Wipe with 70% isopropyl alcohol.</li> <li>Launder with mild detergent in warm water (60° C max), normal wash cycle. Cuff is compatible with 5 wash cycles.</li> </ul> <p><b>Low-Level Disinfection Procedure</b> Prepare Enzol® enzymatic detergent according to manufacturer's instructions. Apply port-cap to cuff. Spray detergent solution liberally onto cuff and use a sterile brush to agitate the detergent solution over entire cuff surface for five minutes. Rinse continuously with distilled water for five minutes. To disinfect, first follow the cleaning steps above, then spray cuff with 10% bleach solution until saturated, agitate with a sterile brush over entire cuff surface for five minutes. Rinse continuously with distilled water for five minutes. Wipe off excess water with sterile cloth and allow cuff to air dry.</p>
<p>FlexiPort™ Disposable Cuffs, tubing and port fittings</p>		<p>Use one or more of the following methods and allow to air dry:</p> <ul style="list-style-type: none"> <li>Wipe with mild detergent and water solution (1:9 solution). Rinse.</li> <li>Wipe with Enzol per manufacturer's instructions. Rinse.</li> <li>Wipe with 0.5% bleach and water solution. Rinse.</li> <li>Wipe with 70% isopropyl alcohol.</li> </ul>