# **SARS-CoV-2 Rapid Antigen Test**

REF	$\nabla$	SYSTEM
9901-NCOV-01G	25	visual reading

English Intended use The SARS-CoV-2 Rapid Antigen Test is a rapid chromatographic immunoassay for the qualitative detection of specific antigens of SARS-CoV-2 present in the human nasopharynx. This test is intended to detect antigen from the SARS-CoV-2 virus in individuals suspected of COVID-19. This product is strictly intended for professional use in laboratory and Point of Care

**Summary**Coronaviruses can cause a variety of acute and chronic diseases. Common signs of a person Coronavruses can cause a vanety of acute and crironic diseases. Common signs of a person infected with a coronavirus include respiratory symptoms, fever, cough, shortness of breath, and dyspnea. In more severe cases, infection can cause pneumonia, severe acute respiratory syndrome, kidney failure, and even death. The 2019 new coronavirus, or SARS-CoV-2, was discovered due to Wuhan viral pneumonia cases in 2019 and a pandemic was declared by the World Health Organization on March 11, 2020. WHO confirmed that COVID-19 can cause colds and more serious diseases such as severe acute respiratory syndrome (SARS).

Test principle
The SARS-CoV-2 Rapid Antigen Test has two pre-coated lines: A "C" Control line and a "T" Test line on the surface of the nitrocellulose membrane. Both the control line and test line in the result window are not visible before applying any samples. Mouse monoclonal anti-SARS-CoV-2 antibody is coated on the test line region and mouse monoclonal anti-Chicken IgY antibody is coated on the control line region. Mouse monoclonal anti-Chicken IgY antibody conjugated with color particles are used as detectors for the SARS-CoV-2 antitigen device. During the test, the SARS-CoV-2 antitigen in the sample interacts with monoclonal anti-SARS-CoV-2 antibody conjugated with color particles making an antipen-antibody configuration. SARS-CoV-2 antigen in the sample interacts with monoclonal anti-SARS-CoV-2 antibody conjugated with color particles making an antigen-antibody color particle complex. This complex migrates on the membrane via capillary action to the test line, where it is captured by the mouse monoclonal anti-SARS-CoV-2 antigens are present in the sample. The intensity of the colored test line varies depending upon the amount of anti-SARS-CoV-2 antigen present in the sample.

Note: Even if the test line is very faint or not uniform the test result should be interpreted as a positive result. If anti-SARS-CoV-2 antigens are not present in the sample, no color appears in the test line. The control line is used for procedural control, and always appears if the test result should be interpreted as a supplication.

#### Reagents

- mAb anti-COVID19 antibody
- mAb anti-Chicken IgY
- mAb anti-COVID-19 antibody-gold conjugate
- Purified chicken IgY-gold conjugate

### Precautions and warnings Do not re-use the test kit

- Do not use the test kit if the pouch is damaged or the seal is broken.
- Do not use the extraction buffer tube of a different lot.
- Do not smoke, drink or eat while handling sample.
- Wear personal protective equipment, such as gloves and lab coats when handling kit reagents. Wash hands thoroughly after the tests are done.

is valid. If no control line is visible the test result should be considered as invalid.

- Clean up spills thoroughly using an appropriate disinfectant.
- Handle all samples as if they contain infectious agents.
- Observe established precautions against microbiological hazards throughout testing
- Dispose all samples and materials used to perform the test as biohazardous waste. Laboratory chemical and biohazardous wastes must be handled and discarded in accordance
- with all local, state, and national regulations.

  Desiccant in foil pouch is to absorb moisture and keep humidity from affecting products. If the desiccant status indicator changes from yellow to green, the test device in the pouch should

## be discarded.

Storage and stability
Store the kit at 2-30 °C / 36-86 °F out of direct sunlight. Kit materials are stable until the expiry date printed on the outer box. Do not freeze the kit

## Materials provided

- Test device (individually in a foil pouch with desiccant)
- Extraction buffer tube Nozzle cap
- Sterile swab
- Film (can be attached to the test device when performing outdoor testing)
- Instructions for use
- Quick Reference Guide

## Materials required (but not provided)

## Timer

# Test preparation and sample collection Carefully read the instructions for using the SARS-CoV-2 Rapid Antigen Test. Please also see the enclosed Quick Reterence Guide (with illustrations) before performing a test.

- 1. Check the expiry date on the back of the foil pouch. Do not use the test, if the expiry date has
- 2 Open the foil pouch and remove the test device and the desiccant package. Use the test immediatedly after opening the pouch.
- 3. Ensure that the test device is undamaged and that the desiccant status indicator shows valid

### 4. Perform a QC as required according to the Instructions for Use of the QC material. Collecting a sample (Nasopharyngeal swab)

- 1. To collect a nasopharyngeal swab sample, insert a sterile swab into the nostril of the patient,
- reaching the surface of the posterior nasopharynx.

  2. Using gentle rotation, push the swab until resistance is met at the level of the turbinate.
- 3. Rotate the swab a few times against the nasopharyngeal wall.
- 4. Remove the swab from the nostril carefully.
- 5. Insert the swab into the provided extraction buffer tube. While squeezing the buffer tube, stir the swah more than 5 times
- 6. Remove the swab while squeezing the sides of the tube to extract the liquid from the swab.
- 7. Press the nozzle cap tightly onto the tube.
- 8. The sample should be tested as soon as possible after collection.
- 9. Samples may be stored at room temperature for up to 1 hour or at 2-8 °C/ 36-46 °F for up to 4 hours prior to testing.

## Preparing a sample from viral transport media

Prepare a sample from a viral transport medium as snown	repare a sample from a viral transport medium as snown in the QHG illustration.						
Viral transport medium (VTM)	Recommended s	torage conditio					
Viral transport medium (VTM)	2 °C to 8 °C	<b>25</b> ℃					
Recommended VTMsa)	12 hours	8 hours					

Viral transport medium (VTM)	Recommended storage condition			
That transport moduli (TTM)	2 °C to 8 °C	<b>25</b> ℃		

(VTM), it is important to ensure that the VTM containing the sample is warmed to room temperature. Cold samples will not flow correctly and can lead to erroneous or invalid results. Several minutes will be required to bring a cold

a) Only use the following VTMs: Copan UTM™ Universal Transport Media, BD™ Universal Viral Transport, STANDARD™ Transport Medium.

#### Test procedure

- 1. Apply 3 drops of extracted sample to the specimen well of the test device.
- 2. Read the test result at 15-30 minutes ♠Do not read test results after 30 minutes. It may give false results.

# Reading and interpreting results:

- A colored line appears in the top section of the result window to show that the test is working properly. This line is the control line (C). Even if the control line is faint or not uniform, the test should be considered to be performed properly. If no control line is visible the test result should be considered as invalid.
- In case of a positive result, a colored line appears in the lower section of the result window. This line is the test line of the SARS-CoV-2 antigen (T). Even if the test line is very faint or not uniform the test result should be interpreted as a positive result.

QC
A control kit including positive and negative quality control is available separately from Roche (STANDARD™ COVID-19 Ag Control, SD Biosensor).

- The test procedure, precautions and interpretation of results for this test must be followed
- strictly when testing. The test should be used for the detection of SARS-CoV-2 antigen in human nasopharyngeal
- swab samples.

  This is a qualitative test, therefore quantitative values of SARS-CoV-2 antigen concentration cannot be determined.
- The immune response cannot be assessed with this test and needs other testing methods.
- The test result should not be used as a sole basis for treatment or patient management decisions, and should be considered in the context of the patient's recent exposures, history and the presence of clinical signs and symptoms consistent with COVID-19.
- A negative result may occur if the concentration of antigen in a sample is below the detection limit of the test or if the sample was collected or transported improperly. Therefore a negative test result does not eliminate the possibility of SARS-CoV-2 infection, and should be confirmed by viral culture or a molecular assay or ELISA, if necessary for patient management.
- Positive test results do not rule out co-infections with other pathogens.
   Positive test results do not differentiate between SARS-CoV-2 and SARS-CoV.
- Negative test results are not intended to rule in or rule out other coronavirus infection.

## Specific performance data

Specific performance data

Clinical evaluation

The sensitivity of the SARS-CoV-2 Rapid Antigen Test for rapid detection of SARS-CoV-2 antigen was established in prospective, randomized, single blinided studies conducted during the SARS-CoV-2 pandemic in Brazil and India. A total of 115 positive samples were tested using the SARS-CoV-2 Rapid Antigen Test. These samples consisted of nasopharyngeal swabs from symptomatic and asymptomatic patients. The specificity of the SARS-CoV-2 Rapid Antigen Test was tested using 311 negative samples. The sensitivity and specificity of the SARS-CoV-2 Rapid Antigen Test was compared to commercialized molecular assays.

Test sensitivity & specificity

The SARS-CoV-2 Rapid Antigen Test showed 96.52 % of sensitivity and 99.68 % of specificity.

			PCR		
		Positive	Negative	Total	
SARS-CoV-2	Positive	111	1	112	
Rapid Antigen Test	Negative	4	310	314	
	Total	115	311	426	
Sensitivity		96.52 % (111/115, 95 % CI 91.33-99.04 %)			
Specificity 99.68 % (310/311, 95 % CI 98.22-9		22-99.99 %)			

## Analytical performance

1. Limit of detection (LoD):
The study used the "SARS-CoV-2 (2019-nCOV) NCCP 43326/2020/Korea" strain. The titer of cultured virus was confirmed by PCR. The cell is inactivated and spiked into a nasopharyngeal swab sample. The LoD is 3.12 x 10²² TCID<sub>20</sub>/ml.

	2019-nCoV Strain Tested: NCCP 43326/2020 / Korea									
	Stock 2019-nCoV Titer: 1 X 10 <sup>6.2</sup> TCID <sub>50</sub> /ml									
Dilution	1/ 10	1/ 100	1/ 200	1/ 400	1/ 800	1/ 1600	1/ 3200	1/ 6400	1/ 1280- 0	1/ 2560- 0
Concen- tration <sup>b)</sup>	1 X 10 <sup>5.2</sup>	1 X 10 <sup>4.2</sup>	5 X 10 <sup>3.2</sup>	2.5 X 10 <sup>3.2</sup>	1.25 X 10 <sup>3.2</sup>	6.25 X 10 <sup>2.2</sup>	3.12 X 10 <sup>2.2</sup>	1.56 X 10 <sup>2.2</sup>	7.8 X 10 <sup>1.2</sup>	3.9 X 10 <sup>1.2</sup>
Call rate (5)c)	100- % (5/5)	100- % (5/5)	100- % (5/5)	100% (5/5)	100% (5/5)	100% (5/5)	100% (5/5)	0% (0/5)	0% (0/5)	0% (0/5)
Call rate (20) <sup>d)</sup>	NA	NA	NA	NA	NA	100% (20/2- 0)	100% (20/2- 0)	0% (0/20)	NA	NA
Lowest concentration with uniform positivity per parameter: 3.12 X 10 <sup>22</sup> TCID <sub>50</sub> /r								3.12 X 10	22 TCID	<sub>50</sub> /ml

Limit of Detection (LoD) per virus strain: 3.12 X 10<sup>22</sup> TCID<sub>50</sub>/ml

b) in dilution tested TCID<sub>so</sub>/ml

c) of 5 replicates

d) of 20 replicates near cut-of

There was no cross-reaction with potential cross-reactive substances except SARS coronavirus.

ross-reactivity testing with SARS-CoV-2 negative samples:						
Virus/ Bacteria/ Parasite	Strain	Concentration	Results			
SARS Coronavirus	Urbani <sup>e)</sup>	3.5 μg/ml	POS			
MERS Coronavirus	Jeddah_1_2013 <sup>f)</sup>	10 μg/ml	NEG			

Virus/ Bacteria/ Parasite	Strain	Concentration	Results
	Type 1g)	3 X 10 <sup>5</sup> TCID <sub>50</sub> /ml	NEG
	Type 3g)	1.5 X 10 <sup>6</sup> TCID <sub>50</sub> /ml	NEG
	Type 5 <sup>g)</sup>	4 X 10 <sup>5</sup> TCID <sub>50</sub> /ml	NEG
	Type 7 <sup>g)</sup>	1.5 X 10 <sup>6</sup> TCID <sub>50</sub> /ml	NEG
Adenovirus	Type 8g)	4 X 10 <sup>5</sup> TCID <sub>50</sub> /ml	NEG
	Type 11 <sup>g)</sup>	4 X 10 <sup>5</sup> TCID <sub>50</sub> /ml	NEG
	Type 18 <sup>g)</sup>	4 X 10 <sup>5</sup> TCID <sub>50</sub> /ml	NEG
	Type 23 <sup>g)</sup>	4 X 10 <sup>5</sup> TCID <sub>50</sub> /ml	NEG
	Type 55 <sup>g)</sup>	4 X 10 <sup>5</sup> TCID <sub>50</sub> /ml	NEG
	H1N1 Denverh)	3 X 10 <sup>5</sup> TCID <sub>50</sub> /ml	NEG
	H1N1 WS/33h)	3 X 10 <sup>5</sup> TCID <sub>50</sub> /ml	NEG
Influenza A	H1N1 Pdm-09h)	3 X 105 TCID <sub>50</sub> /ml	NEG
	H1N1 New Caledoniah)	3 X 10 <sup>5</sup> TCID <sub>50</sub> /ml	NEG
	H1N1 New Jerseyh)	3 X 10 <sup>5</sup> TCID <sub>50</sub> /ml	NEG
	Nevada/03/2011h)	3 X 10 <sup>5</sup> TCID <sub>50</sub> /ml	NEG
Influenza B	B/Lee/40h)	2.5 X 10 <sup>4</sup> TCID <sub>50</sub> /ml	NEG
	B/Taiwan/2/62h)	3 X 10 <sup>5</sup> TCID <sub>50</sub> /ml	NEG
Respiratory syncytial virus	Type A <sup>h)</sup>	3 X 10 <sup>5</sup> TCID <sub>50</sub> /ml	NEG
Respiratory syncytial virus	Type b <sup>h)</sup>	3 X 10 <sup>5</sup> TCID <sub>50</sub> /ml	NEG
Lastandla	Bloomington-2h)	5 X 10 <sup>4</sup> cells/ml	NEG
Legionella pneumophila	Los Angeles-1h)	5 X 10 <sup>4</sup> cells/ml	NEG
	82A3105h)	5 X 10 <sup>4</sup> cells/ml	NEG
	K <sub>i)</sub>	5 X 10 <sup>4</sup> cells/ml	NEG
	Erdman <sup>i)</sup>	5 X 10 <sup>4</sup> cells/ml	NEG
Mycobacterium tuberculosis	HN878 <sup>i)</sup>	5 X 10 <sup>4</sup> cells/ml	NEG
	CDC1551 <sup>i)</sup>	5 X 10 <sup>4</sup> cells/ml	NEG
	H37Rv <sup>i)</sup>	5 X 10 <sup>4</sup> cells/ml	NEG
	4752-98 [Maryland (D1)6B-17] <sup>h)</sup>	5 X 10 <sup>4</sup> cells/ml	NEG
Streptococcus pneumonia	178 [Poland 23F-16]h)	5 X 10 <sup>4</sup> cells/ml	NEG
priodinonia	262 [CIP 104340] <sup>h)</sup>	5 X 10 <sup>4</sup> cells/ml	NEG
	Slovakia 14-10h) [29055]	5 X 10 <sup>4</sup> cells/ml	NEG
Streptococcus pyrogens	Typing strain T1 [NCIB 11841, SF 130] <sup>h)</sup>	5 X 10 <sup>4</sup> cells/ml	NEG
	Mutant 22h)	5 X 10 <sup>4</sup> cells/ml	NEG
Mycoplasma pneumoniae	FH strain of Eaton Agent [NCTC 10119] <sup>h)</sup>	5 X 10 <sup>4</sup> cells/ml	NEG
	M129-B7 <sup>h)</sup>	5 X 10 <sup>4</sup> cells/ml	NEG
Pooled human nasal wash <sup>jj</sup>	NA <sup>k)</sup>	NA	NEG
	229E <sup>1)</sup>	1 X 10 <sup>4.5</sup> cells/ml	NEG
Coronavirus	OC43 <sup>()</sup>	1 X 10 <sup>5</sup> cells/ml	NEG
	NL63 <sup>()</sup>	1 X 10 <sup>4</sup> cells/ml	NEG
MERS Coronavirus	Florida /USA-2_Saudi Arabia_2014 <sup>()</sup>	4 X 10 <sup>4</sup> TCID <sub>50</sub> /ml	NEG
Human Meta- pneumo virus 3 (Type B1)	Peru2-2002 <sup>()</sup>	1 X 10 <sup>5</sup> cells/ml	NEG
Human Meta- pneumovirus 16 (Type A1)	IA10-2003 <sup>()</sup>	1 X 10 <sup>5</sup> cells/ml	NEG
	Type 1 <sup>I)</sup>	1 X 10 <sup>5</sup> cells/ml	NEG
Parainfluenzavirus	Type 2 <sup>()</sup>	1 X 10 <sup>5</sup> cells/ml	NEG
ı aralılılucı izavilus	Type 3 <sup>()</sup>	1 X 10 <sup>5</sup> cells/ml	NEG
	Type 4/A <sup>I)</sup>	1 X 10 <sup>5</sup> cells/ml	NEG
		1 X 10 <sup>5</sup> cells/ml	NEG

		Type 4/A <sup>I)</sup>	1 X 10 <sup>5</sup> cells/ml	١		
	Rhinovirus A16	N/A <sup>I)</sup>	1 X 10 <sup>5</sup> cells/ml	١		
e) BEI / inactivated virus						
	f) Bionote / recombinant	protein				
	g) Korea Bank for Patho	genic Viruses / live				
	h) ATCC / live virus					

j) to represent diverse microbial flora in the human respiratory tract

k) Bionote / Normal pooled human nasal wash from healthy employees SD biosensor / Normal pooled human nasal wash from healthy employees

I) Zentomatrix / inactivated

Note: Human coronavirus HKU1 has not been tested. The % identity of the nucleocapsid protein sequence between HKU1 and SARS-CoV-2 is below 35 %

# 3. Endogenous / exogenous interference substances studies: There was no interference on the test result from potentially interfering substances listed below SARS\_CAV\_2 onetitive and negrative samples were tested.

Potential interfering substance	Concentration	Results	
Respiratory s	samples	'	
Mucin: bovine submaxillary gland, type I-S	100 μg/ml	NEG	
Blood (human), EDTA anticoagulated	5 % (v/v)	NEG	
Biotin	100 μg/ml	NEG	

Potential interfering substance		Concentrati	ion	Res	ults
Nasal s	prays or dro	ps			
Neo-Synephrine (Phenylephrine)		10 % (v/v)		NE	G
Afrin Nasal Spray (Oxymetazoline)		10 % (v/v)		NE	3
Saline Nasal Spray		10 % (v/v)		NE	3
Homeopathic a	allergy relief	medicine			
Homeopathic Zicam Allergy Relief Nasal Ge	el	5 % (v/v)		NE	G
Sodium Cromoglycate		20 mg/ml		NE	3
Olopatadine Hydrochloride		10 mg/ml		NE(	G
Anti-	viral drugs				
Zanamivir (Influenza)		5 mg/ml		NE(	G
Oseltamivir (Influenza) 10 mg		10 mg/ml		NE	3
Artemether-lumefantrine (Malaria)	temether-lumefantrine (Malaria) 50 µM			NEG	
Doxycycline hyclate (Malaria)		70 μM		NEG	
Quinine (Malaria)		150 µM		NEG	
Lamivudine (Retroviral medication)		1 mg/ml		NEG	
Ribavirin (HCV)		1 mg/ml		NE	G
Daclatasvir (HCV)		1 mg/ml		NEG	
Anti-inflam	matory medi	cation			
Acetaminophen		200 μΜ		NE	G
Acetylsalicylic acid		3.7 mM		NE	G
Ibuprofen		2.5 mM		NEG	
А	ntibiotic				
Mupirocin		10 mg/ml		NE	G
Tobramycin		5 μg/ml		NE	G
Erythromycin		81.6 μM		NE	G
Ciprofloxacin		30.2 μM		NE	G
) Results from interference testing with SAF	8S-CaV-2 no	itive samples			
Potential interfering substance	Concer		Viral strai	in	Res-

Ant	ibiotic				
Mupirocin		10 mg/ml		NE	G
Tobramycin		5 μg/ml		NE	G
Erythromycin		81.6 μM		NE	G
Ciprofloxacin		30.2 μM		NE	G
) Results from interference testing with SARS	<del></del>				
Potential interfering substance	Concen	tration	Viral stra level <sup>m)</sup>	ain	Res- ults <sup>n)</sup>
Respirato	ory sample	es	•		
Mucin: bovine submaxillary gland, type I-S	100 μg/r	ml	SARS		POS
Blood (human) EDTA anticoagulated	5 % (v/v	·)	CoV-2	d	POS
Biotin	100 μg/ml		virus me	virus media	
Nasal spr	ays or dro	ps			
Neo-Synephrine (Phenylephrine)	10 % (v/	/v)	SARS		POS
Afrin Nasal Spray (Oxymetazoline)	10 % (v/v)		CoV-2 cultured virus media		POS
Saline Nasal Spray	10 % (v/v)				POS
Homeopathic alle	ergy relief	medicine			
Homeopathic Zicam Allergy Relief Nasal Gel	5 % (v/v	')	SARS CoV-2		POS
Sodium Cromoglycate	20 mg/n	nl	culture	d	POS
Olopatadine Hydrochloride	10 mg/n	nl	virus me	dia	POS
Anti-vi	ral drugs				
Zanamivir (Influenza)	5 mg/ml				POS
Oseltamivir (Influenza)	10 mg/n	nl	1		POS
Artemether-lumefantrine (Malaria)	50 μM		SARS	-	POS
Doxycycline hyclate (Malaria)	70 µM		CoV-2		POS
Quinine (Malaria)	150 µM		virus me		POS
Lamivudine (Retroviral medication)	1 mg/ml		0)		POS
Ribavirin (HCV) 1 mg/ml	1 mg/ml				POS
Daclatasvir (HCV)	1 mg/ml				POS
Anti-inflamma	<u> </u>				
Acetaminophen	200 μM		SARS CoV-2		POS
Acetylsalicylic acid	3.7 mM		culture	d	POS
Ibuprofen	2.5 mM		virus me	ula	POS

Lamivudine (Retroviral medication)	1 mg/ml		POS
Ribavirin (HCV) 1 mg/ml	1 mg/ml		POS
Daclatasvir (HCV)	1 mg/ml		POS
Anti-infla	mmatory medication	•	
Acetaminophen	200 μΜ	SARS- CoV-2	POS
Acetylsalicylic acid	3.7 mM	cultured	POS
Ibuprofen	2.5 mM	virus media	POS
	Antibiotic	'	
Mupirocin	10 mg/ml	SARS-	POS
Tobramycin	5 μg/ml	CoV-2	POS

Mupirocin	10 mg/ml	SARS- CoV-2 cultured virus media
Tobramycin	5 μg/ml	
Erythromycin	81.6 μM	
Ciprofloxacin	31 μM	

o) 1/800 dilution (1.25 X 1032 TCID50/ml)

m) in multiples of LoD

n) detected X/3

4. High-dose hook effect: SARS-CoV-2 cultured virus was spiked into sample. SARS-CoV-2 cultured virus did not show hook-effect at 1 X  $10^{6.2}$  TCID $_{50}$ /ml.

A point (period/stop) is always used in this Method Sheet as the decimal separator to mark the border between the integral and the fractional parts of a decimal numeral. Separators for thousands are not used.



Reference number



Batch code

Caution



in vitro diagnostic medical device



This product fulfills the requirements of the European Directive 98/79/EC



Consult instructions for use



Contains sufficient for <n> tests



Use-by date Temperature limit



Analyzers/Instruments on which reagents can be used



Do not re-use

Unique Device Identifier

Date of manufacturing

Global Trade Item Number



Do not use if package is damaged



Manufacture



Keep product dry

Keep away from sunlight

Additions, deletions or changes are indicated by a change bar in the margin.



## SD BIOSENSOR

Head office: C-4th&5th, 16, Deogyeong-daero 1556beon-gil, Yeongtong-gu, Suwon-si, Gyeonggi-do, 16690 REPUBLIC OF KOREA

Manufacturing site: 74, Osongsaengmyeong 4-ro, Osong-eup, Heungdeok-gu,
Cheongju-si, Chungcheongbuk-do, 28161 REPUBLIC OF KOREA

# www.sdbiosensor.com

POS POS Roche Diagnostics GmbH, Sandhofer Strasse 116, D-68305 Mannheim www.roche.com Roche order number: 09327592

EC REP Authorized Representative
MT Promedt Consulting GmbH Altenhofstrasse 80 66386 St. Ingbert Germany Phone: +49 6894 581020, Fax: +49 6894 581021