SAFEGUARD
MEDICAL™

Prometheus Omni Splint™

Lightweight Aluminium Foam Splint

Contents: 1Pcs

MD

Medical
Device

Does not
contain or
presence of
natural rubber
latex

Non-sterile

Do not re-use

Consult
instructions
for use

Manufacturer

Safeguard Medical Technologies Ltd.
The Old Rectory
Hope Under Dinmore
Herefordshire, HR6 0PW
United Kingdom

Importer

Authorized Representative
in the EC

Safeguard Technologies Ltd
Willow Grove, Delgany, A63 XY89
Co Wicklow, Ireland

Prometheus Omni Splint™

Lightweight Aluminium Foam Splint

EN

Intended Use

The Prometheus Omni Splint™ is intended to be used to aid support and immobilise a patient's limb(s).

User Profile

The Prometheus Omni Splint™ is indicated for use by laypersons and medical professionals.

Patient Profile

The Prometheus Omni Splint™ is indicated for use on patients with suspected limb fracture.

Clinical Benefits

The expected clinical benefit is the reduction of pain and damage caused by limb fracture for up to 24 hours, until definitive care is reached.

Product Claims

- Malleable
- Waterproof
- Lightweight and compact
- Natural Rubber Latex Free

Performance Characteristics

Secured in place with medical tape or bandage of choice.

Contraindications



N/A

Residual Risks/Warnings

- Do not re-use. Re-use risks cross-infection or loss of functionality.
- Dispose of immediately after use, as clinical waste as per local protocols. Improper disposal risks cross contamination.
- Use past shelf-life risks loss of functionality.
- Failure to monitor neurological function and distal circulation below the fracture risks loss of limb function or loss of the limb itself.
- Remove before diagnostic imaging. Use with X-ray or MRI risks excessive heating.
- Flammable, keep away from open flame. Exposure to open flame risks respiratory tract irritation.
- If irritation or allergy occurs, remove device immediately.
- Excessive length of immobilisation for over 24 hours risks deep vein thrombosis, joint stiffness or a decrease in bone mineral density.

Instructions for Use

1. Open and flatten the splint.
2. Create a 'C' shaped bend along with the length of the splint. This provides extra strength when immobilising the limb.
3. Apply the splint to the fractured limb so the fracture is between the two sides of the splint.
4. Adjust the shape of the splint to conform to the limb, if needed, and secure with medical tape or bandage of choice.
5. Monitor neurological function and distal circulation below the fracture before and after splintage.

Reporting

Any serious incident that has occurred in relation to this device should be reported to manufacturer vigilance@safeguardmedical.com and applicable national competent authority.

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