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

DE

10.7 CLOSING THE SLIDES (ONLY STAIR+)	21
10.8 PARKING BRAKES	21
10.9 CLOSING THE CHAIR	21
10.10 MOVING THE CHAIR ON THE GROUND	22
10.11 LIFTING THE CHAIR	22
10.12 MOVING THE CHAIR ON STAIRS	22
10.13 APPLICATION OF SEAT XL	22
10.14 APPLICATION OF HEADREST (ONLY STAIR+)	22
11. CLEANING AND MAINTENANCE	23
11.1 CLEANING	23
11.2 PRECAUTIONARY MAINTENANCE	23
11.3 PERIODIC MAINTENANCE	24
11.4 SPECIAL SERVICING	24
11.5 LIFE SPAN	24
12. TROUBLESHOOTING	24
13. ACCESSORIES	25
14. SPARE PARTS	25
15. DEMOLITION	25

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EN

1. MODELS	15
2. INTENDED USE	15
3. REFERENCE STANDARD	15
4. INTRODUCTION	15
4.1 USE OF THE MANUAL	15
4.2 LABELLING AND TRACKING CONTROL OF THE DEVICE	15
4.3 SYMBOLS	15
4.4 WARRANTY AND SUPPORT	15
5. WARNINGS	16
6. SPECIFIC WARNINGS	17
6.1 PHYSICAL REQUIREMENTS OF OPERATORS	18
7. RESIDUAL RISK	18
8. TECHNICAL DATA AND COMPONENTS	18
9. INSTALLATION AND START-UP	19
10. INSTRUCTIONS FOR USE	19
10.1 MEDICAL VEHICLE REQUIREMENTS	19
10.2 OPENING THE CHAIR	20
10.3 REAR HANDLES OPENING, ARMRESTS (WHERE PRESENT) AND FOOTREST	20
10.4 OPENING AND CLOSING OF TELESOPIC HANDLES	20
10.5 OPENING THE SLIDES (ONLY STAIR+)	20
10.6 EXTRACTABLE HANDLE OPENING (ONLY STAIR+)	21

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1. MODELS

The standard following models can undergo change, revision and implementation without any notice.

- 4BELL SILVER+
- 4BELL STAIR+
- 4BELL AVIO+

2. INTENDED USE

Intended use

Folding transport chairs are devices to be used to move and bring a patient into a sitting position in the ambulance, but not to be used for transporting the patient inside the ambulance. They must not be used for parking.

Intended users

The intended users are people trained in first aid procedures and the use of medical equipment in an EMS environment.

Among the possible users are also contemplated the fitters of emergency vehicles who can handle the product before putting it into service or during any maintenance of the vehicle on which the chair is installed.

Patients

The expected patients are those with temporary motor difficulties.

There are no known contraindications or side effects resulting from the use of the device, as long as it is used in accordance with the user manual

3. REFERENCE STANDARD

As a distributor or end user of products manufactured and/or marketed by Spencer Italia S.r.l. you are strictly required to have basic knowledge of any legal requirements applying to the devices contained in this supply that are in power in the goods final destination Country (including laws and norms regarding technical specifications and/or safety requirements) and therefore you are strictly required to have the necessary knowledge to guarantee all aspects regarding the total conformity of the products to the regulations in the relevant territory.

REFERENCE	DOCUMENT TITLE
UNI EN 1865-1 (not applicable to model Stair)	Patient handling equipment used in road ambulances - Part 1: Specification for general stretcher systems and patient handling equipment
UNI EN 1865-4 (Only mod. Stair)	Patient handling equipment used in road ambulances – Part 4: Foldable patient transfer chair

4. INTRODUCTION

4.1 USE OF THE MANUAL

This manual is intended to provide the health care operator with the all the necessary information for its safe and appropriate use as well as adequate maintenance of the device.

Note: The manual is an integral part of the device. It must be kept for the duration of the device and must accompany the device in case of change of ownership or destination. If the operating instructions received relate to products not received, you must immediately contact the manufacturer before use.

The Spencer product manuals can be downloaded from the website or can be requested by <http://support.spencer.it> or by contacting the manufacturer. Exceptions are items whose essentiality for reasonable and predictable use is such as to make it unnecessary to prepare instructions in addition to the following warnings and directions on the label.

Regardless of the level of experience gained in the past with similar devices, it is recommended that you carefully read this manual before installing, operating or using the product or any maintenance.

4.2 LABELLING AND TRACKING CONTROL OF THE DEVICE

Each device has got an identification label positioned on the device itself and/or on its box, which includes identification data about the manufacturer, the product, the CE mark, the serial number (SN) or lot number (LOT). It must never be removed or covered.

In case of damage or loss, request a duplicate from the manufacturer. Failure to do so will interrupt the validity of the guarantee as the device can no longer be traced.

The Regulation UE 2017/745 requires manufacturers and distributors of medical devices to keep track of the device location. If the device was in a different location to the address where it was sent or to where it had been sold, donated, lost, stolen, exported, or destroyed, permanently removed from use, or if the device had not been delivered directly from Spencer Italia S.r.l., register your device at <http://service.spencer.it> or inform the customer (see § 4.4).

4.3 SYMBOLS

Symbols	Meaning	Symbols	Meaning
	Class I Medical Device, compliant with the Regulation UE 2017/745		Caution: Federal Law restricts this device to sale by or on the order of a licensed practitioner
	Medical Device		General or specific warnings
	Manufacturer		See instructions for use
	Date of manufacture		Product code
			Serial Number
UDI	 (01)08057711230006(11)200626(10)1234567890	Production identification Alphanumeric code that identifies the production units of the device, composed of: (01)0805771123 company prefix 000 progressive GS1 6 control number (11)200626 date of production (YYMMDD) (10) 1234567890 lot/SN	

4.4 WARRANTY AND SUPPORT

Spencer Italia S.r.l. guarantees that products are without defects for a period of **one year from the date of purchase**.

For any information regarding the correct interpretation of the instruction manual, the use, maintenance, installation and restore of the product, contact the Spencer customer

IT

EN

DE

FR

ES

PT

PL

BG

RO

care service ph. +39.0521.541111, fax +39.0521.541222, e-mail service@spencer.it or visit <https://www.spencer.it/en/medical-emergency-company> to find the nearest service centre for assistance.

In order to facilitate the assistance service, please always indicate the lot number (LOT) or serial number (SN) shown on the label applied on the box or on the device.

Conditions for warranty and assistance can be viewed on <http://support.spencer.it>

Note: Record and store; these instructions, lot (LOT) or serial number (SN) if any, date and place of purchase, date of first use, date servicing, user names and comments.

5. WARNINGS

The warnings, notes and other important safety information are indicated in this section and are clearly visible throughout the entire manual.

 At least every 6 months, is important to verify if updated User Manuals are available or there is any change involving your product. This information is freely available on the website <http://support.spencer.it> in the section "User manuals" and "Product Updates".

User training

Note: Laboratory testing, post production tests, instruction manuals cannot always consider every possible scenario for use. This means that in some cases the performance of the product could be notable different from results to date obtained.

Instructions are continually being updated and are under tight surveillance of fully qualified staff with adequate technical formation.

- Regardless of the level of experience gained previously with similar devices, it is recommended that you carefully read this manual before installing, operating or using the product or carrying out any maintenance. If in doubt, contact Spencer Italia S.r.l. to obtain the necessary clarifications.
- The product must be used by trained personnel only, having attended specific training for this device and not for similar products.
- The suitability of the user to use the product may be attested by the records of training, where the names of those trained, of the trainers, dates and place are indicated. **This register which certifies the eligibility of the operators to use the Spencer device must be kept for a period of 10 years after the disposal of the device itself. This register will be made available to the competent Authorities and/or manufacturer if requested. In the absence of such documentation, sanctions will be applied.**
- Do not allow any untrained person to help during the use of the device, because they could cause damage to the patient or to themselves.
- The product should be operated only by personnel trained in the use of this product and not in others similar.

Note: Spencer Italia S.r.l. is always at your disposal to organise product training.

Installers training

The installation of the device must be carried out only by qualified and trained staff authorized by Spencer Italia S.r.l. Dates and procedures for participating in training courses will be arranged between the customer and our Commercial offices.

Product functionality

Use of the device in anyway other than described in this manual is forbidden.

- Before each use, the perfect operating state of the device must be thoroughly checked as specified in the instruction manual. If any damage or abnormalities which could in any way influence the correct functioning and the safety of the device are detected, the device must be immediately removed from service and the manufacturer must be contacted.
- If any failure or incorrect functioning of the device is detected, it must be immediately substituted with a similar item so that the rescue procedures are guaranteed without any interruption.
- The appliance must not in any way be tampered with (modification, adjustment, addition, replacement). In such cases all responsibility for any malfunctions or injuries caused by the appliance itself will be denied; moreover CE certification and product warranty will be considered void.
- Those who modify or have modified, prepare or have prepared medical appliances in such a way that they no longer serve the purpose for which they were intended, or no longer supply the intended service, must satisfy the valid conditions for the introduction onto the market.
- During use, position and adjust the device taking care not to cause any obstruction to rescuers and or any other rescue equipment.
- Ensure that all the necessary precautions are taken in order to avoid the hazards that can arise as the result of contact with blood or body fluids, when applicable.
- For devices intended for the transport of patients, always respect the maximum load capacity of the device, as indicated in this user's manual. Maximum load capacity means the total weight distributed according to the human anatomy. In determining the total loading weight of the product, the operator must consider the weight of the patient, the equipment and the accessories. Moreover, the operator must consider that the overall dimensions of the patient do not reduce the functionality of the device.
- For devices intended for the transport of patients, make sure, before lifting, that operators have appropriate physical condition as indicated in the manual.
- **The maximum weight supported by each rescuer must comply with the requirements prescribed by the law of the land, in the field of Health and Safety at Work.**
- **The warranty seals, where present, must not be removed;** in such case, the manufacturer will no longer recognize the product warranty and will accept no responsibility in case of incorrect operation or damage caused by the product itself.
- Avoid contact with sharp objects
- The installation of the device must be carried out only by qualified and trained staff authorized by Spencer Italia S.r.l. Timing and procedures for conducting training courses shall be agreed between the customer and our Commercial offices.
- Operating temperature: from -10°C to +50°C.

Storage

- The device should not be exposed to or come into contact with any source of combustion or inflammable agents. Store in a cool, dry, dark place and do not expose to direct sun.
- Do not store the device underneath any heavy objects which could cause structural damage.
- Store and transport the device in its original packaging. Failure to do so makes the warranty void.
- Storage temperature: from -20°C to +60°C.

Maintenance/cleaning

Spencer Italia S.r.l. disclaims any liability for any damage, direct or indirect, which is a result of improper use of the product and replacement parts and / or otherwise of any repairs made by an entity other than the authorized Spencer service centres; this will also invalidate the warranty.

- The operator must always wear adequate personal protection such as gloves and mask etc. during all checking, maintenance and cleaning procedures.
- Establish a maintenance program and periodic testing, identifying an employee responsible for overseeing. The person to whom the ordinary maintenance of the device is entrusted must ensure the basic requirements foreseen by the manufacturer in the user's manual.
- **The frequency of inspection is determined by factors such as legal requirements, the type of use, frequency of use, environmental conditions during use and storage.**
- Spencer Italia S.r.l. disclaims any liability for any damage, direct or indirect, that is the result of incorrect repairs or use of products made by Spencer Italia S.r.l. Repairs must necessarily be carried out by an authorized Spencer Italia service centre, which in using original spare parts will provide a quality repair service in strict accordance with the technical specifications given by the manufacturer. Spencer Italia S.r.l. disclaims any liability for any damage, direct or indirect, which is a result of improper use of spare parts and/or otherwise of any repairs made by an entity other than the Spencer service centres authorized to repair or make substitutions on this product and parts and/or otherwise of any repairs made by an entity other than the Spencer service centres authorized to do so; the warranty will also be invalidated.
- Use only original components, spare parts and or accessories, approved by Spencer Italia S.r.l., in order to carry out any operation without causing any alteration or modification to the device.
- For any operations that are not carried out directly by the manufacturer but by an authorised centre, we have to underline that a report regarding all operations carried out must be requested. This will permit both Spencer Italia S.r.l. and the end user to keep a log book regarding the operations carried out on the device.
- All maintenance and revision must be recorded and documented with the corresponding report for technical assistance; documentation shall be maintained for at least 10 years from the end of life and must be made available to the competent authorities and/or the manufacturer if requested.

- The cleaning schedule for reusable products must be performed in accordance with the directions provided by the manufacturer in the user manual, in order to avoid the risk of cross-infection due to the presence of secretions and/or residuals.
- The device and all its components, after washing, should be allowed to dry completely before storing.
- If required, lubrication must be carried out after cleaning and complete drying.

Regulatory requirements

As a distributor or end user of products manufactured and/or marketed by Spencer Italia S.r.l., you are strictly required to have a basic knowledge of any legal requirements applying to the devices contained in this supply that are in power in the final destination Country (including laws and norms regarding technical specifications and/or safety requirements) of the goods and therefore you are also strictly required to have the necessary knowledge to guarantee all aspects regarding the total conformity of the products to the regulations in the relevant territory.

- Promptly notify Spencer Italia S.r.l. (already during the first product enquiry) when requesting in details regarding any revisions to be made by manufacturer in order to guarantee the conformity of the products to the territory's legal specifications (including those resulting from rules and/or norms of other kind).
- Act, with all due care and diligence, and contribute to ensure conformity to general safety requirements of all devices marketed in the territory, by providing final users with all necessary information for carrying out periodical checks on their devices, exactly as specified in the relevant user's manual.
- **Actively contribute to safety checks on product**, sold, by communicating any relevant risk analysis information both to the manufacturer and to any competent authorities so that the necessary action can be promptly taken.
- The distributor or final user is aware that in the event of any failure to conform to the above mentioned requirements you will be deemed fully responsible for all damages that might occur. Therefore Spencer Italia S.r.l. expressly disclaims any responsibility and/or liability for your non-compliance with the present regulatory provisions.

General warnings for medical devices

The user must carefully read not only these general warnings, but also those listed below.

- It is not foreseen that the use of the device is prolonged beyond the time necessary for recovery operation described in the intended use and in this manual.
- When the device is being used, the assistance of qualified staff must be guaranteed and at least two operators must be present.
- Follow the procedures and protocols approved by the internal organization.
- The activities of disinfection and sterilization should be carried out in accordance with the parameters given in the validated cycle as specified in the technical standards.
- With reference to the Regulation UE 2017/745, we remind both public and private operators, , That in the exercise of their activity detect an accident involving a medical product are required to notify the Ministry of Health, under the terms and in the manner established by the relative ministerial decrees and also to the manufacturer. Health care providers whether public or private are required to communicate to the manufacturer, any other inconvenience that may allow for the adoption of measures to ensure the protection and health of patients and users.

6. SPECIFIC WARNINGS

For the use of the Transport Chair 4BELL+, , the user must have read, understood and follow carefully all the instructions described in this manual.

- Follow the procedures approved by the Emergency Medical Service for the immobilization and transport of patients.
- Follow the procedures approved by the Emergency Medical Service for the positioning and transport of patients.
- Do not use if the device or its parts are pierced, torn, frayed or excessively worn out.
- The device can be used with a patient only if the belts are intact and properly fastened
- Handling of the device on the wheels, is permitted only with all four wheels resting on the floor. Is not permitted to carry the chair tilted or resting on only two wheels.
- When the device is moved with the patient on it, is essential to have a firm grip on the handles according to the type of transport.
- The device can be lifted using front and rear handles. Before lifting, make sure that the front handles are fully extended and locked.
- Before use, ensure that all operators have a firm grip on the device.
- Avoid pulling the device on rough surfaces.
- Do not lift by crane or other mechanical lifters.
- Do not use drying machines.
- Do not use the device in combination with other equipment not expressly approved by the manufacturer.
- The device is intended to be used for patient transport as stated in the chapter Intended use and **cannot be used as a parking device or for the transport of patient inside the ambulance.**
- Practice with an empty chair and with a simulated load in order to become familiar with the manoeuvres.
- For patient loading techniques, for particularly heavy patients, for working on uneven ground or in special and unusual circumstances, the presence of more operators is recommended (not only 2 as required in standard conditions).
- Before each use, check the integrity of the belts and their hooks as specified in the user's manual. In case of malfunction or damage that may compromise the functioning and safety of the device, patient or operator, it is necessary to replace the belts.
- Make sure the belts are properly fastened to the device.
- Make sure that the tapes used to fix the seat fabric are properly secured.
- Always immobilize the patient, using the straps supplied by the manufacturer; lack of immobilization may cause serious damage. The manufacturer is not liable for any damage caused by improper or incorrectly applied belts.
- Do not operate in case the weight has not been distributed correctly.
- Before moving the patient with the device, inform him about the manoeuvres you are going to perform
- Use only the handles or knobs to lift the chair and not any other parts of the device.
- **Hold the chair firmly still when occupied by the patient.**
- **The Stationing brakes are aids for the operator, they do not in any way substitute his supervision.**
- Great attention must be paid to any possible obstacles (water, ice, debris, etc.) on the route of the stretcher/chair, because they could cause loss of balance of the operator and compromise the proper functioning of the device. If the path cannot be freed of obstacles, choose an alternative route.
- To avoid injury, before lifting the chair, always check that the lifting handles are securely locked and fully extended.
- It is preferable not to use the device in the event of suspected cervical or spinal trauma, or fractures.
- To avoid any risk to the safety of the patient and rescuer, during transport on the stairs, the presence of at least two operators is necessary.
- Condensation, water, ice and dust accumulation can affect the correct functioning of the device, making it unpredictable and causing a sudden change in weight that operators have to sustain.
- For obstacles higher than 10mm the device must be lifted using the handles. Impact against obstacles above this height, can compromise the integrity and security of the device.
- Do not alter or modify the device arbitrarily: the modification may cause unforeseeable functioning and damage to the patient or operators. In any case the warranty will be void and the manufacturer relieved from any liability.
- When closing the device, make sure nothing interferes with the handling systems. Carelessness during this operation could be the cause of crush injuries.
- Make sure that the slides, if present, are properly locked before moving the device on stairs.
- The 4BELL+ transport chair, is fulfils the requirements of the EN 1789 standard if used with a Spencer fastening system compliant to the standard. Therefore, the use of fastening systems is not permitted unless approved by the manufacturer. Fastening systems that have not been approved may alter the structural and functional characteristics of the device.
- **Do not use the device near MRI machines.**

IT

EN

DE

FR

ES

PT

PL

BG

RO

6.1 PHYSICAL REQUIREMENTS OF OPERATORS

The 4BELL+ transport chair is a device intended for professional use only. Each operator must be trained to transport patients safely and efficiently. Do not allow untrained persons to help operators during use of the product, as this may cause injury to themselves or to other people.

The operators that use the device must have the physical ability to use it and good muscle coordination, as well as presenting a strong back, arms and legs to raise and support and be able to grasp the device firmly with both hands.

Operators must be able to provide the necessary assistance to the patient.

Users should be able to lift and handle safely the weight of the chair and the patient as well as any other equipment used with the device.

During loading procedures of extremely heavy patients, operations on rough terrain or in particular circumstances, the presence of more operators may be needed (not only 2 as required in normal conditions). The ability of all operators must be considered before determining their role in the use of the chair.


7. RESIDUAL RISK

The residual risks listed below have been identified exclusively in reference to the intended use of the device.

- The use by untrained personnel may result in injury to the patient, rescuer or others.
- Inadequate disinfection procedures may involve cross-infection risks.
- Failure to immobilize the patient on the device or immobilization performed using incomplete belt systems or incorrectly fixation to the chair, can result in serious injury to the patient.
- Partial extraction of front handles could cause their sudden closing when lifting the device. Make sure the handles are fully extracted and locked before lifting.
- Patient intervention on the device and its closing and braking mechanisms can cause unlocking of the device from open position or unexpected movements with consequent risks for the patient. Make sure the patient does not intervene on any element of the device.
- Failure to close handles when their extraction is not required can cause obstruction during device handling. Make sure the handles are closed when moving the device on wheels.
- Partial opening of the chair can cause its unexpected closure if lifted with consequent risks for the patient. Always verify the proper insertion of the locking system.
- A slides locking failure, for models equipped with them, can cause their sudden closure during transport on stair with serious consequences for the patient. Always check that all locking mechanisms are properly engaged.
- The use of fasteners not approved by the manufacturer may result in instability of the device inside the ambulance, causing the risk of violent impact of the device against people and/or things.
- Failure to comply with warnings given to operators, can cause crushing risks related to the mechanisms of the device.
- Failure to respect the life span established by the manufacturer and the periodic revisions described in this manual, can result in serious damage for patient and operators linked to the use of a product that does not comply with the safety requirements established by the manufacturer.
- Failure to read and understand the instructions, can have consequences for the patient and operators.

8. TECHNICAL DATA AND COMPONENTS

Note: Spencer Italia S.r.l. reserves the right to make changes to specifications without prior notice.

N°	ELEMENT	DESCRIPTION AND MATERIALS
1	Backrest frame	Made of aluminum, it acts as a transport handle as well as a backrest for the patient
2	Armrests (optional)	Made of aluminum, they allow the patient to rest their arms offering greater comfort
3	Backrest belts	Made of polyester and equipped with a quick release buckle, restrains the patient's torso to the backrest frame
4	Seat belts	Made of polyester and equipped with a quick release buckle, restrains the patient's pelvis to the seat frame
5	Seat sheet	Made of PVC, it is the surface on which the patient sits
6	Unlocking lever	Made of aluminum, if activated allows to unlock the closure system of the chair
7	Telescopic handles	Made of aluminum with PVC knob, they can be extracted for lifting maneuvers
8	Footrest	Made of aluminum, it allows the patient's feet to rest
9	Front wheels	Made of rubber / PP / PA, they are pivoting and allow to easily maneuver the chair
10	Rear wheels	Made of PU, they allow movement of the chair
11	Rear brakes	Made of nylon, if activated they perform a braking action on the rear wheels
12	Lower rear handles (optional)	Made of aluminum, they allow the lifting of the chair

13	Slides (only STAIR+)	Made of steel, aluminum and PU, once open they allow the chair to slide downwards on stairs
14	Upper rear handles	Made of aluminum, they allow the lifting of the chair
15	Extractable handle	Made of aluminum, it can be extracted allowing operators to grab the chair while descending the stairs
16	Head immobilizer belt (Optional only on STAIR +)	Made of polyester, it allows to stabilize the patient's head
17	Seat XL	Made of aluminum, it offers a wider support surface

The technical specifications of the main models are as follow:

FEATURE	AVIO+	SILVER+	STAIR+
Length opened w/handles closed	760 ± 5 mm	760 ± 5 mm	760 ± 5 mm
Length opened w/handles opened	1120 ± 5 mm	1120 ± 5 mm	1120 ± 5 mm
Length w/slides and handles opened	/	/	1280 ± 5 mm
Width	465 ± 5 mm	490 ± 5 mm	490 ± 5 mm
Height opened	905 ± 5 mm	905 ± 5 mm	940 ± 5 mm
Height w/headrest risen	/	/	1475 ± 5 mm
Dimensions folded	290x465x905 ± 5 mm	290x490x905 ± 5 mm	290x490x950 ± 5 mm
Rear wheels	ø150 mm	ø150 mm	ø150 mm
Front wheels	ø100 mm	ø100 mm	ø100 mm
Weight	9,3 ± 0,2 kg	9,4 ± 0,2 kg	14,7 ± 0,2 kg
Loading capacity ¹	180 kg	180 kg	180 kg

¹ Maximum load capacity means the total weight distributed according to the human anatomy. In determining the total load of the weight on the product, the operator must consider the weight of the patient, the equipment and the accessories.

9. INSTALLATION AND START-UP

Before the first use verify that:

- The packaging is intact and has protected the device during transport.
- Check that are present all the components included in the accompanying list.
- Functionality of the device according to the user manual.
- The ambulance is equipped with a dedicated Spencer fastener and its installation will not hinder access to any other devices.
- Sheets and belts are properly fastened. It is important to make sure that the buckles under the seat are properly inserted and the strap bands fully adhere to each other.
- Upper belts are properly fixed to the backrest.

Do not modify for any reason any part or component of the chair as this may cause injury to the patient and/or rescuers.

⚠ If conditions mentioned above are not satisfied, the device cannot be considered secure, safe use is compromised and it is a possible cause of injury for patient and operators and can cause damage to the device.

Practice with a chair without patient before the regular start-up. For instructions on use after start-up, follow the operations described on paragraph 12.

If the above conditions are met, the device may be considered ready for use; otherwise you must immediately remove the device from service and contact the manufacturer.

Do not alter or modify in any way the appliance; any such interference could cause malfunctions and injury to the patient and/or rescuer. Moreover the manufacturer will no longer recognize the product warranty and will accept no responsibility.

10. INSTRUCTIONS FOR USE

Before transferring, lifting or transporting the patient, primary medical evaluations have to be performed. Once the diagnosis has been assured, it is preferable (if possible) to suggest to the patient that he actively collaborates during transfer onto the chair, and to make sure the patient is fully aware of all risks. Before loading the patient, place the device near the patient; always maintain the safety belts fastened during the manoeuvres.

10.1 MEDICAL VEHICLE REQUIREMENTS

When placed inside the ambulance, the chair 4BELL must be properly anchored in order to limit risks for occupants of the vehicle linked to violent movements resulting from rapid acceleration or deceleration of the vehicle. For this purpose, it is essential to use a dedicated fastener. For its installation, it is necessary that the emergency vehicle is equipped with suitable surfaces to hold safely the assembly consisting of fastener and chair.

The user manual of the fastener contains more details about requirements of support surfaces and installation.

⚠ Failure to use the dedicated fastener or its improper installation can seriously affect the patient's and operator's safety.

IT

EN

DE

FR

ES

PT

PL

BG

RO



10.2 OPENING THE CHAIR

The 4BELL transport chair is maintained in closed position thanks to the belts on the backrest. Before opening the chair, release the quick release buckle of the upper belts.

Firmly hold the transport handle, press on the frame of the seat as shown by the arrow, until all the wheels rest on the floor.

 The hand with which the seat frame opens must be placed on the upper part of the frame  pushing it until it is completely open. Do not place your hands in the area identified by the image  on the side, it is not needed to open the chair.

Firmly press down the seat frame until the locking system is activated.

10.3 REAR HANDLES OPENING, ARMRESTS (WHERE PRESENT) AND FOOTREST



Rear handles and armrests (if present), can be opened by turning them as shown by the arrows.

To close the upper rear handles, pull the release control towards you and move each handle to the closed position.

To close the rear lower handles and armrests, simply rotate them in the opposite direction to that of opening.

Similarly, open the footrest by rotating it outwards until it is completely opened.

After use, close the footrest by rotating in the opposite direction until it is completely closed.

10.4 OPENING AND CLOSING OF TELESCOPIC HANDLES

The telescopic handles have been made to facilitate lifting of the chair. They are kept in position by a mechanism that is inserted automatically when reaching the excursion limit both during opening and closing.


To unlock the mechanism allowing handle movement, press downwards and move simultaneously.

 Pull the handles until they are fully open. It is essential to verify that the locking system has been properly inserted.

To do this, try to slide repeatedly the handles without applying pressure downwards. There must be no movement.

DO NOT LIFT THE CHAIR IF THE HANDLES ARE NOT COMPLETELY EXTRACTED AND LOCKED.

Similarly it is possible to close the telescopic handles verifying that they are locked in closed position.


10.5 OPENING THE SLIDES (ONLY STAIR+)


Load the patient on the device and secure him using all the belts.

Grab the crossbar of the slide and push in the opening direction until the mobile frame that connects the slide to the frame of the chair reaches its end of stroke.


 Make sure that the mobile frame that constrains the movement of the slide has reached the end of its stroke. If this does not happen, it could be the source of a sudden closure of the slide during movement on stairs.

10.6 EXTRACTABLE HANDLE OPENING (ONLY STAIR+)

To open the pull-out handle, grasp it at the end and lift it until the catches come out allowing the handle to be blocked. Make sure that the red latches have completely come out of their seat and that the handle is locked.


To close the handle, push the sliding rings down to ensure that the locks return to their seat. After the latches have been closed, the pull-out handle can be pushed down until completely closed.

10.7 CLOSING THE SLIDES (ONLY STAIR+)

To close the slides, press with a foot on the crosspiece of the mobile frame, then return the slide to the closed position.

It is recommended to always close the slides when no longer needed.

10.8 PARKING BRAKES

4Bell+ chairs are equipped with parking brakes that block the rear wheels. To lock the rear wheels, press the brakes on the operator's side until the wheel stops.


To unlock the wheels, lift the brake pedal in the opposite direction to that of activation.

Always insert both brakes.

During any handling or working operation with the patient on the chair, the operator must have a firm grip on the device, regardless of the brakes having been activated or not.

10.9 CLOSING THE CHAIR

Place all the handles, headrest, armrests and slides in closed position.



Using one hand, raise in the direction shown by the arrow, until the part A of the frame is almost in contact with the part B of the seat frame. This movement will disengage the locking system of the chair.

Then lift the seat until it touches the backrest. The front wheels will be positioned automatically.

As for opening, even during closing do not place your hands between the tubulars "A" and "B", but move upwards by pressing only on the tubular "B".

Before placing the chair on the fixing system, secure it in closed position using the integrated belts anchored to the chair backrest.

The belts must pass through free space between the handles and the frame of the backrest and must be attached at the level of the grips of the telescopic handles.

IT

EN

DE

FR

ES

PT

PL

BG

RO

10.10 MOVING THE CHAIR ON THE GROUND

All the chairs 4Bell+, can be easily moved on the ground thanks to the wheels fitted in each model.

Use the transport handle to move the chair and, once on the rescue place, with or without the patient on the device, always activate the parking brakes.

⚠ For ground movements of all models, it is necessary to close rear and front telescopic handles as they would hinder device's handling, being unused items during this type of transport. **During patient transport, it is essential to fasten all the belts.**

For models equipped with slides, keep them closed until reaching any stairs, as they would unnecessary interfere with operators movements. If you are not carrying a patient on the device, it is suggested to carry the chair in closed position to limit overall dimensions and for a more practical handling.

10.11 LIFTING THE CHAIR

To carry the chair on stairs or on surfaces where the use by wheels could be difficult, it is possible to lift the device using the telescopic handles and rear handles. For this method of transport the presence of more operators is required with at least one at the front side and at least one on the back side of the chair.

⚠ Before lifting it is fundamental to verify that telescopic handles are fully extended and locked as described in par. 11.4. Failure to control could cause serious injury to patient and operators due to the sudden closure of the telescopic handles.

Once on a surface allowing transport on wheels, is recommended to place the chair on the ground in order to avoid unnecessary load on operators.

10.12 MOVING THE CHAIR ON STAIRS

The chair can be transported upstairs or downstairs by lifting the chair as described in the previous paragraph. The descent of the stairs can be easier by using the slides available on Stair + model.

To move the chair downstairs by means of the slides, after the slides have been opened as described before, proceed as follows:

⚠ To move the chair downstairs, approach the first step of the stairs, tilt the chair by holding the headrest, keeping it balanced accompanying it during the descent until the slides are lying on at least two steps. The operator on the rear side of the chair, must always maintain a firm grip on the extractable handle, accompanying the device regulating the descent speed.

The second operator should never be in front of the device, but should be placed at the side of the device and a few steps lower, and has to guide the descent optimally and must be ready to intervene in the event of need, without compromising its safety. This operator must act on the front handles to direct the chair by making it follow the safest and most linear path during the descent of the stairs.

Once the flat floor is reached, bring the chair back to the vertical position, so that the front wheels adhere to the ground, continuing the transport with standard ground handling procedures.

To carry the chair upstairs, do not use the slides, but lift the chair like described before.

When lifting the chair and going upstairs, the first operator should be placed behind the device grabbing the rear handles and facing the back to the patient. The second operator should be placed in front of the chair facing the patient and grabbing the telescopic handles of the chair making sure they are properly locked in open position.

Using the proper lifting technique, the operators should lift simultaneously the device and start the transport.

⚠ For this type of transport is recommended the presence of a supporting operator which will guide the maneuvers of other operators optimally.

10.13 APPLICATION OF SEAT XL

The XL seat is equipped with 4 plastic elements that must be fitted onto the free portions of the seat frame.

Once the correct alignment has been identified, simply press lightly until it clicks into place.

10.14 APPLICATION OF HEADREST (ONLY STAIR+)

The 4Bell Stair + chair can be equipped with an optional head immobilizer. To apply it, proceed as follows:


Position yourself on the back of the chair and open the extractable handle. Position the head immobilizer band on the patient side of the removable handle.

Place one of the loops near the tube of the extractable handle and start wrapping the free end of the tape around the same tube.


Wrap the tape completely around the tube and pass it through the loop from the bottom side.


Pull the tape so that the tape is tightened around the tube.

Make the tape pass through the other loop and pull the tape.

Once you have carried out the same operations on the other side and have checked that the tape is sufficiently tensioned to remain in position on the extractable handle, make the tear-off portions of the tapes adhere. You may need to make some adjustments. These can be done by sliding the tape inside the loops, thus compensating for excess or shortage of length in the desired portions.

After positioning the patient, the head can be wrapped with the band.

IT

EN

DE

FR

ES

PT

PL

BG

RO

11. CLEANING AND MAINTENANCE

11.1 CLEANING

Failure to carry out the correct cleaning routine could increase the risk of cross infection, due to presence of body fluids and/or residuals. The operator must always wear adequate personal protection such as gloves and mask, etc. during all checking and cleaning procedures.

The exposed metal parts are usually treated and/or painted in order to increase their resistance to external agents. Clean the exposed parts with water and a delicate soap; **never use solvents or stain removers.**

Rinse thoroughly with warm water making sure that you have removed all traces of detergent, which could ruin or compromise the integrity and durability of the device. **The use of high pressure water should be avoided,** because water penetrates in the joints removing the lubricant and creating the risk of corrosion of components. Allow the device to dry thoroughly before storing. Drying after washing or after use in wet environments must be natural and not forced; do not use flames or other sources of direct heat. If the device needs to be **disinfected**, use products that do not have corrosive or solvent action on the materials of which the device is made. Be sure to take every precaution to ensure that there is no risk of cross-infection or contamination for patients and / or operators.

To maintain the polished appearance of the frame parts we recommend the use of **Spencer STX 99 polish cleaner** or in alternative creme or wax normally used for polishing car bodywork.

11.2 PRECAUTIONARY MAINTENANCE

Establish a maintenance programme and periodic testing routine and identify an employee responsible for this. The person to whom the ordinary maintenance of the device is entrusted must ensure that the basic requirements foreseen by the manufacturer in following paragraphs are inspected. All maintenance and periodic servicing activities must be registered and kept together with the servicing reports. These documents have to be kept for a period of 10 years after the disposal of the device itself. This register will be made available to the competent authorities and/or manufacturer if requested.

Routine maintenance of the device must be carried out by operators in possession of specific qualifications, trained and experienced in the use and maintenance of the device. Such training must be recorded and maintained for at least 10 years from the end of life of the product and shall be made available to the competent authorities and / or the Manufacturer upon request.

The operator must always wear adequate personal protection such as gloves and mask, etc. during all checking and cleaning procedures.

Checks to be carried out before and after each use and at deadline indicated above, are as follows:

- General functionality of the device
- Cleanliness of the device (remember that the failure to clean could be the cause of cross infections)
- Correct fixation of all nuts, bolts and screws

IT

EN

DE

FR

ES

PT

PL

BG

RO

- Absence of cuts, holes, tears on the structure, including the straps
- None of the tubes or metal sheets present bends or cracks
- The seat or the backrest has no structural damage or cracks
- The seat belts, sheets, moving parts, wheels and handles are intact and functioning properly.
- Lubrication of moving parts
- Wearing on wheels and breaking system
- Wear of slides and their proper locking in the open position. Verify the correct functionality of the locking spring.
- The belts of the slides have a good sliding action
- The wheels are correctly fixed, they are stable and turn properly
- The wheels are free from debris.
- The device opens and locks properly
- The device unlocks and closes properly
- Functioning of springs
- The patient's immobilization belts are present and they are intact and functioning.
- Telescopic handles can be positioned, folded and locked properly in all positions as described in this manual.
- Headrest (if present) works and locks properly
- Verify the functionality of the system which locks the chair in open position. If you try to close the chair without actioning the unlocking system, the chair must not close.
- The emergency vehicle is equipped with a Spencer fastening system intended for the chair
- The connection of the two devices is always correctly inserted and safe.

The inspection frequency is determined by factors such as local legal requirements, the type of use, frequency of use, environmental conditions during use and storage.

Before lubricating any part, if needed, be sure to remove any dust or dirt in the area.

Please note that you must do the cleaning as described in this manual and verify functionality before and after each use. Spencer Italia S.r.l. declines any responsibility for the improper functioning or injury caused to the patient or user by the use of devices that have not been subjected to a routine maintenance programme which will void the warranty and the compliance to the Regulation UE 2017/745.

 **Every 2 years** the seat and backrest PVC sheets and patient belts must be substituted.

Use only accessories/original spare parts approved by Spencer Italia S.r.l. Failing to do so, we will accept no responsibility for the incorrect functioning and/or damage caused by the use of any device which has not been repaired, or certified on expiry date by the manufacturer or by one of the manufacturer's authorised service centres. Warranty will be considered void in compliance with the Regulation UE 2017/745.

11.3 PERIODIC MAINTENANCE

The device must be serviced by the manufacturer or by an authorised centre, every year.

If the correct revision is not carried out, the device MUST BE PUT OUT OF SERVICE, because the CE branding will no longer be considered valid and consequently it is no longer compliant with the safety standards declared by the manufacturer at time of purchase.

Spencer Italia S.r.l. will take no responsibility for the incorrect functioning or any damage caused by a device that has not undergone regular revision.

The belts, seat and backrest sheets must be necessarily replaced **every 2 years**.

For any operations that are not carried out directly by the manufacturer but by an authorised centre, we have to underline that a report regarding all operations carried out must be requested. This will permit both Spencer Italia S.r.l. and the end user to keep a log book regarding the operations carried out on the device.

11.4 SPECIAL SERVICING

Only the manufacturer or centres with written authorisation are authorised to complete any special servicing operations.

For any operations that are not carried out directly by the manufacturer but by an authorised centre, we have to underline that a report regarding all operations carried out must be requested. This will permit both Spencer Italia S.r.l. and the end user to keep a log book regarding the operations carried out on the device.

The end user is authorised to replace only the spare parts indicated in the paragraph 15.

11.5 LIFE SPAN

The device, if used as indicated in the following instruction manual, has a life span of 5 years starting from the purchase date. The life span can be extended following the annual revision. PVC sheets, seat belts and backrest have to be replaced every 2 years from the purchase date.

General revisions must be carried out by the manufacturer or by a centre authorised by the manufacture. **If such annual revisions are not carried out, the device MUST BE DISPOSED ACCORDING TO THE PROCEDURES SPECIFIED IN PARAGRAPH 16 AND THIS EVENT MUST BE NOTIFIED TO THE MANUFACTURER.**

The life span can be extended, by unquestionable judgment of Spencer or authorized service center, if the safety requirements of the device are still guaranteed. Spencer Italia S.r.l. will accept no responsibility for the incorrect functioning and/or damage caused by the use of any device which has not been serviced by the manufacturer or authorised centre, or of any device for which the life span is expired.

12. TROUBLESHOOTING

PROBLEM	CAUSE	SOLUTION
The seat sheet doesn't support the patient's weight	The quick release buckles are not properly engaged	Verify the proper closure of the buckles placed on the lower side of the sheet
	Worn cloth	Replace the sheet as described in the manual
The telescopic handles don't lock in closed or opened position	The telescopic handles have not been completely extracted or inserted.	Extract or insert the telescopic handles until the end position is reached
	The locking system is damaged.	Put the chair immediately out of service and contact the service centre
The headrest doesn't lock in open position	The locking system has not been correctly inserted	Make sure that the locking handle is properly turned. Check also that no element of the locking lever is damaged. If anomalies are present, put the device immediately out of service and contact the service centre
The brakes do not work	Worn wheels or damaged mechanism	Put the chair immediately out of service and contact the service centre
The rear handles are unstable or deform under loading	Handle assembly damaged	Put the chair immediately out of service and contact the service centre
The device is difficult to move on wheels or is not stable	One or more wheels are damaged	Put the chair immediately out of service and contact the service centre

When placed on the fastener, the chair is not stable or doesn't remain closed	The coupling between the chair and the fastener was not done properly.	Consult the users manual of the fastener
	The fixing system in use is inappropriate for the chair.	Put the device out of service and replace the fixing system with an appropriate one
The spare sheets haven't enough tension or can't be installed because are too tight	The purchased sheets are for the narrower or the wider models	Contact your dealer to obtain the correct item
Slides do not lock properly	The slides are not completely opened	Verify that the slides are completely opened like described in this manual
	The slides or their locking systems are damaged	Put the chair immediately out of service and contact the service centre
The belt of the slides does not run smooth	Accumulated dirt, wear or failure of the adjustment mechanisms	Attempt to clean with compressed air, if the problem persists, immediately put the device out of service and contact assistance

13. ACCESSORIES

ST10453B	4BELL+ MAX FASTENER (COMPATIBLE WITH 4BELL+)	ST10451	REAR AUXILIARY HANDLES 4BELL+
ST10455B	4BELL STAIR+ MAX FASTENER (COMPATIBLE WITH 4BELL STAIR+)	ST10452	SEAT XL
ST10450	ARMRESTS	ST10454	HEADREST 4BELL STAIR+

14. SPARE PARTS

RISTO140	SPARE SEAT COVER 4BELL+
-----------------	-------------------------

15. DEMOLITION

When the device is no longer suitable for use, if they haven't been contaminated by any particular agents, they can be disposed of as normal solid waste. Otherwise follow the current regulations for demolition.

IT

EN

DE

FR

ES

PT

PL

BG

RO

Warning

The information contained in this manual is subject to change without notice.
The Diagrams are inserted only for reference and may vary slightly from the actual device. Spencer Italia S.r.l. assumes no responsibility for any errors contained herein or for damage, accidents or consequences connected with the supply, performance or use of this manual.

The information contained in this document could be modified without any warning and is not to be intended as a commitment on behalf of Spencer Italia S.r.l.
Spencer products are exported to many countries and the same identical regulations are not always valid. For this reason there could be differences between the description here described and the product actually delivered. Spencer continually strives to reach the perfection of all items sold. We therefore hope you will understand if we reserve the right, at any time, to modify the shape, equipment, lay-out or technical aspects that are herein described.

© Copyright Spencer Italia S.r.l.

All rights reserved. No part of this document can be photocopied, reproduced or translated into another language without the written approval of Spencer Italia S.r.l.

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