**HEINE** QUALITY  
MADE IN GERMANY

## HEINE® EN200

CE

**HEINE Optotechnik GmbH & Co. KG**  
Kientalstr. 7 · 82211 Herrsching · Germany  
Tel. +49(0)81 52/38-0  
Fax +49(0)81 52/38-202  
E-Mail: [info@heine.com](mailto:info@heine.com) · [www.heine.com](http://www.heine.com)  
**MED 113066**      **2016-11-25**

# HEINE® EN200

DEUTSCH	3
ENGLISH	6
FRANÇAIS	9
ESPAÑOL	12
ITALIANO	15
SVENSKA	18
NEDERLANDS	21
DANSK	24
NORSK	27
SUOMI	30
PORTUGUÊS	33

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 Please read and follow these instructions for use and keep them for future reference.

### Intended Use

The HEINE® EN200 Wall Transformer is a two-channel, dimmable 3.5 V power supply with a 5 V output for medical HEINE Diagnostic Instruments. It should only be operated by trained medical staff within a professional health facility.

 **For U.S. only:**  
Federal law restricts this device to sale by or on the order of a Physician or Practitioner!

### Warnings and Safety Information

 **CAUTION!** Indicates potential hazardous situations. Ignoring the corresponding instructions may lead to dangerous situations of mild to moderate extent. (Background color yellow; foreground color black).

 **NOTE!** Indicates valuable advice in terms of installation, operation, maintenance or repair. Notes are important, but not related to hazardous situations.

### Assembly and disassembly

Assembly accessories

Allen wrench 1.3

 To ensure all-pole disconnection from the mains at all times, the ME device must be installed so that the switch mode power supply is accessible and disconnectable.

Before mounting the device, check whether you need to use a special wall plug for the subsurface of your wall and whether the subsurface exhibits sufficient load-bearing capacity for the mechanical stress.

 The wall plugs supplied with the product are universal wall plugs and suitable for most building materials (e.g. concrete, solid brick, brick).

1. Hold the drill template horizontally in place on the desired position and mark the drill holes. Drill a hole of at least 40 mm deep using a  $\varnothing$  6 mm drill.
2. Fit wall plugs.
3. Use a suitable screwdriver to tighten the two upper screws with a distance of approx. 3 mm to the wall.
4. Hang the device – disconnected from the mains – on the two screw heads, press firmly and push downwards. Check to make sure that the device is correctly positioned on both screw heads.
5. Screw in the two lower screws.

To disassemble, disconnect the device from the mains first, loosen the two lower screws, push upwards and take off the wall.

### USB connector

The device has an USB interface on its rear side as power supply for additional medical diagnostic instruments from HEINE.

Connect the HEINE USB cord and run the cable via the housing canal to the outside on the left or right hand side. To do so, press the cable into the clamping hubs, which act as strain relief.

### Setting up

Primary plug

EU

USA / Japan

UK

Australia

 The device may only be operated with the mains voltage specified on the nameplate of the power supply.

- Attach the primary plug to the plug-in transformer until it engages.

- Connect the device to the mains with the plug-in transformer. The green indicator light on the EN200 illuminates.

### Connecting instruments

The HEINE® EN200 cable handles are fitted with HEINE automatic connectors for use with all similarly equipped HEINE instruments.

- Hold the handle in one hand and insert the instrument connector into the guides in the handle head.
  - Push the instrument into the handle until it locks in place with an audible "click".
  - Instruments can be locked to prevent unauthorized removal by using the enclosed Allen wrench according to fig. 1.
  - To unlock, the screw may only be unscrewed so far until the connector can be turned.
  - The instrument is ready for operation when taken off the unit.
  - To adjust the brightness of the instrument activate the adjustment ring on the handle.
  - Both handles can be simultaneously used independent of each other.
  - To remove the instrument turn the upper conical ring to the right and pull the instrument off.
-  The device enters power saving mode, if the handles are taken off the unit and are not used for more than 15 minutes. The handles are activated by adjusting the brightness or by placing the handle back and taking it off the unit again.


Fig. 1

The setting up and operation of the HEINE instrument heads are described in a separate instruction of use.

### Hygienic Reprocessing

Instructions on hygienic reprocessing must be adhered to, based on national standards, laws and guidelines.

Classification according to KRINKO: non-critical

Spaulding Classification USA: non-critical

 Before cleaning, disconnect the device from the power source.

In the event of suspected contamination, carry out hygienic preparation of the instrument.

The described cleaning and disinfection measures do not replace the specific rules applicable for the establishment.

HEINE Optotechnik only approves the resources and procedures named in these instructions for use.

Cleaning and disinfection may only be carried out by personnel with sufficient hygienic knowledge.

Observe the instructions of the manufacturer of the reprocessing media.

Do not use any spray disinfectant, or any dripping-wet or heavily-foaming cloths. Do not prepare mechanically.

### Procedure

#### Cleaning and disinfection

Clean and disinfect the device manually (wipe clean and wipe disinfect).

When wiping the handles, hold them down to prevent the ingress of liquid.

Ensure that the instrument is completely dry after reprocessing before you use it again.

#### Recommended agents

Cleaning agent: enzymatic (e.g. neodisher® MediClean)

Disinfectant agent: alcoholic (e.g. Incides® N), agent ortho-phthalaldehyde (e.g. Cidex® OPA) or hydrogen peroxide (e.g. Accel® TB).

## Maintenance

The EN200 and the power supply are maintenance free. We as manufacturer do not stipulate a safety inspection in accordance with § 6 Medical Device Operator Ordinance (MPBetreibV), ref. attachment 1, for the HEINE product EN200.

The expected operational life of the HEINE® EN200 is 5 years.

## Service

The device does not require service.

## General Warnings

 Check the correct operation of the device before use! Do not use the device if there are visible signs of damage.

The power supply must not be opened (risk of electrical hazard).

Use power supply in dry environment only. Avoid splashing water. Do not immerse in liquids.

Do not use the device in fire- or explosive risk area (e.g. oxygen saturated or anaesthetic environments).

The device must not be placed near strong Magnetic fields, e.g. MRI units.

Do not modify the device.

Use only original HEINE parts, spare parts, accessories and power sources.

Repairs shall only be carried out by qualified persons.

## General Notes

 The warranty for the entire product is invalidated if non-genuine HEINE products or non-original parts are used and if repairs or modifications are made to the device by persons not authorized by HEINE. For more information, please visit [www.heine.com](http://www.heine.com).

## Disposal

 The product must be recycled as separated electrical and electronic devices. Please observe the relevant state-specific disposal regulations.

## Electromagnetic Compatibility

Medical electric devices are subject to special precautionary measures with regard to electromagnetic compatibility (EMC). Portable and mobile high frequency communication equipment can affect medical electric devices.

 This is a device in the domestic environment, this device may cause radio interference, so that it may be necessary in this case, to take appropriate remedial measures, as e.g. orientation, new arrangement or shielding of the device or restrict the connection to the site.

The use of accessories, converters or cables other than the ones specified by HEINE might lead to increased emission reduced electrical immunity of the medical equipment.

The device may not be stacked directly near or used directly beside other devices. If the device is to be operated in a stack or with other devices, the device should be watched to ensure it operates properly in this location.

The appendix contains following tables

- Guidance and Manufacturer's declaration – Electromagnetic immunity
- Technical specification
- Explanation of the used symbols

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**Guidance and manufacturer's declaration – electromagnetic emissions**

The device is intended for use in the electromagnetic environment specified below.  
The customer or the user of the device should assure that it is used in such environment.

<b>Emission test</b>	<b>Compliance</b>	<b>Electromagnetic environment – Guidelines</b>
RF emissions CISPR11	Group 1	The device uses RF energy only for its internal function. Therefore, RF-emission is very low and it is unlikely that any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The device is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic Emissions IEC 61000-3-2	Class A	
Voltage Fluctuations/ Flicker Emissions IEC 61000-3-3	Passed	

**Guidance and manufacturer declaration - Electromagnetic immunity**

The device is intended for use in the electromagnetic environment specified below.  
The customer or the user of the device should assure that it is used in such environment.

<b>Immunity test</b>	<b>IEC 60601 test level</b>	<b>Compliance level</b>	<b>Electromagnetic environment – Guideline</b>
Electrostatic discharge (ESD) acc. to IEC 61000-4-2	± 6 kV contact discharge ± 8 kV air discharge	± 6 kV contact discharge ± 8 kV air discharge	Floors should be wood, concrete or covered with ceramic tiles. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transient/ burst IEC 61000-4-4	± 2 kV for mains cables ± 1 kV for input and output lines	± 2 kV for mains cables ± 1 kV for input and output lines	The supply voltage quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	± 1 kV voltage phase – phase, ± 2 kV voltage phase – earth	± 1 kV voltage phase – phase ± 2 kV voltage phase – earth	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	< 5% U <sub>T</sub> , (>95% dip in U <sub>T</sub> ) for 1/2 period 40% U <sub>T</sub> , (60% dip in U <sub>T</sub> ) for 5 periods 70% U <sub>T</sub> , (30% dip in U <sub>T</sub> ) for 25 periods <5% U <sub>T</sub> , (>95% dip in U <sub>T</sub> ) for 5 seconds	< 5% U <sub>T</sub> , (>95% dip in U <sub>T</sub> ) for 1/2 period 40% U <sub>T</sub> , (60% dip in U <sub>T</sub> ) for 5 periods 70% U <sub>T</sub> , (30% dip in U <sub>T</sub> ) for 25 periods <5% U <sub>T</sub> , (>95% dip in U <sub>T</sub> ) for 5 seconds	Mains power quality should be that of a typical commercial or hospital environment. If the user of the device requires continued operation during power mains interruptions, it is recommended that the device be powered by a UPS (uninterruptible power supply) or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

Comment: U<sub>T</sub> is the a.c. supply voltage prior to application of the test level.

### Guidance and manufacturer's declaration – electromagnetic immunity

The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment Guidelines
Conducted RF IEC 61000-4-6	3 Veff 150 kHz to 80 MHz	3 V eff	Portable and mobile RF communication equipment should be used no closer to any part of the device, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.
Radiated HF IEC 61000-4-3	3 V/m 80MHz to 2,5GHz	3 V/m	<p><b>Recommended separation distance:</b>  <math>d = 3,5/3 \cdot \text{SQRT}(P/W)</math>  <math>d = 3,5/3 \cdot \text{SQRT}(P/W)</math> 80 MHz to 800 MHz  <math>d = 7/3 \cdot \text{SQRT}(P/W)</math> 800 MHz to 2,5 GHz</p> <p>where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m).</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey<sup>a</sup>, should be less than the compliance level in each frequency range.<sup>b</sup></p> <p>Interference may occur in the vicinity of equipment marked with the following symbol:</p> 

Note 1: At 80MHz and 800MHz, the higher frequency range applies.

Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

**a** Field strengths from fixed transmitters, such as base stations for radio (cellular / cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy.

To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the device is used exceeds the applicable RF compliance level above, the device should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the device.

**b** Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V / m.

#### Recommended separation distances for portable and mobile RF communication equipment and the device

The device is intended for use in the electromagnetic environment in which radiated RF disturbances are controlled. The customer or user of the device can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communication equipment (transmitters) and the device as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter (W)	Separation distance according to frequency of transmitter (m)		
	150 kHz to 80 MHz $d = 3,5/3 \cdot \text{SQRT}(P)$	80 MHz to 800 MHz $d = 3,5/3 \cdot \text{SQRT}(P)$	800 MHz to 2,5 GHz $d = 7/3 \cdot \text{SQRT}(P)$
0.01	0.1	0.1	0.2
0.1	0.4	0.4	0.7
1	1.2	1.2	2.3
10	3.7	3.7	7.4
100	11.7	11.7	23.3

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in metres (m) can be determined using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

Note 1: At 80 MHz and 800 MHz, the higher frequency range applies.

Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

## Technical specification

Environmental conditions for operation	+10 °C to +35 °C 10 % to 75 % rel. humidity 700 hPa to 1060 hPa
Environmental conditions for storage	+5 °C to +45 °C 45 % to 80 % rel. humidity 500 hPa to 1060 hPa
Environmental conditions for transport	-20 °C to +50 °C 10 % to 90 % rel. humidity 500 hPa to 1060 hPa
Mains power	100 – 240 V AC / 50 – 60 Hz
Current consumption	300 – 150 mA
Output (USB)	typ. 5 V
Output (AV)	typ. 1.6 V – 3.6 V
Class	II
Degrees of protection	IP 40
Dimensions	180 x 234 x 99 mm
Weight	940 g

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### Explanation of utilized symbols

The following symbols are used on the device or on the packaging:

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