CareFusion

MicroLab

Operating Manual – English

Operating Manual – USA

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
Contents

I.	Introduction - English.....	7
	Contraindications.....	8
	Warnings and Cautions.....	8
	Overview	10
	Indication for Use	11
	Getting Started	12
	Calibration Check	22
	Customisation	24
	Administration Mode	26
	Paper Loading	27
	External Printing.....	28
	Switching Off	28
	Charging Procedure	28
	PC connection using SPCS	30
	Looking after your Spirometer.....	31
	Cleaning Instructions	31
	External Surfaces of the Spirometer.....	32
	Cleaning accessories	32
	Cleaning the Transducer.....	33
	Servicing	34
	Trouble Shooting Information.....	34
	Safety Designation per IEC 60601-1	36
	Electromagnetic Compatibility (EMC) to EN60601-1:2007.....	38
	Specification of the MicroLab.....	44
	Spirometry Measurements:.....	45
	Consumables/Accessories	48
	Customer contact information	50
II.	Introduction – English - USA	54
	Contraindications.....	55
	Warnings and Cautions.....	55
	Overview	57

Indication for Use	58
Getting Started	59
Calibration Check	69
Customization	71
Administration Mode	73
Paper Loading	74
External Printing.....	75
Switching Off	75
Charging Procedure	75
PC connection using SPCS	77
Looking after your Spirometer.....	78
Cleaning Instructions	78
External Surfaces of the Spirometer.....	79
Cleaning accessories	79
Cleaning the Transducer.....	80
Servicing	81
Trouble Shooting Information.....	81
Safety Designation per IEC 60601-1	83
Electromagnetic Compatibility (EMC) to EN60601-1:2007.....	85
Specification of the MicroLab.....	91
Spirometry Measurements	92
Consumables / Supporting Products.....	95
Customer contact information	96
 III. Introduction - Français	 99
Contre-indications	100
Avertissements et mises en garde	100
Présentation	102
Indications d'utilisation	103
Démarrage	104
Vérification de l'étalonnage.....	114
Personnalisation	116
Mode d'Administration.....	118
Chargement du Papier.....	119

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Operating manual

I. Introduction - English

The MicroLab is mains/battery operated portable spirometer with the unique combination of ease of use and sophistication. Ease of use is assured through the use of context sensitive help screens, accessed at a touch of a button, that explain every MicroLab feature.

The MicroLab uses the CareFusion Digital Volume Transducer, an extremely stable form of volume transducer, which measures expired air directly at B.T.P.S (Body Temperature and Pressure with Saturated water vapour) thus avoiding the inaccuracies of temperature corrections. The transducer is insensitive to the effects of condensation and temperature and avoids the need for individual calibration prior to performing a test.

Test results may be uploaded to a PC using the optional Spirometry PC Software (SPCS) and patient details may be downloaded to the MicroLab.

Stored data may be printed to the integral thermal printer or to an external HP printer using the USB cable supplied.

Contraindications

- Acute disorders affecting test performance (e.g Vomiting , nausea, vertigo)
- Recent eye surgery (increases in intraocular pressure during spirometry)
- Oral or facial pain exacerbated by a mouthpiece
- recent myocardial infarction
- post- operative thoracic surgery patients
- hyperventilation syndrome

Note: Extensive exhalation might lead to syncope

Warnings and Cautions

The following terms are used as follows in this manual

CAUTION: Possibility of injury or serious damage

WARNING: conditions or practices that could result in personal injury

Note: Important information for avoiding damage to the instrument or facilitating operation of the instrument.

Note: Patients below 4 years of age may struggle to perform spirometry correctly and reproducibly

Note: The device should only be used by trained and qualified personnel.

	CAUTION: Read the manual before use
WARNING: The instrument is not suitable for use in the presence of explosive or flammable gases, flammable anaesthetic mixtures or in oxygen rich environments.	
CAUTION: Mouthpieces are single patient use. If used on more than one patient there is a risk of cross-infection. Repeat use may degrade materials and lead to an incorrect measurement.	
CAUTION: Pulmonary filters are single patient use. If used on more than one patient there is a risk of cross-infection. Repeat use may increase air resistance and lead to an incorrect measurement.	
	PLEASE NOTE: The product you have purchased should not be disposed of as unsorted waste. Please utilise your local WEEE collection facilities for the disposal of this product.
PLEASE NOTE: Degree of protection against Ingress of Water is IPX0.	
WARNING: To avoid risk of electric shock, this equipment must only be connected to a supply mains with protective earth	
WARNING: Do not connect devices that are not specified as part of the system.	

Overview

When the unit is turned on for the first time the system customise menu is displayed. Touch the name of the required region with the stylus stored on the left hand side of the unit and then touch 'Finish'. This will set the regional defaults including language, if available.

If this screen is not displayed then touch the 'Customise' icon on the main menu followed by the 'System' icon.


The MicroLab uses a touch screen with icons representing each function available. A stylus, housed in the left hand side of the unit, is provided for icon screen activation.

Touch the displayed time to adjust time and date. Touch the toolbox icon to adjust volume and brightness.

Unused icons may be disabled by touching the blue background and selecting from the list displayed.

Four levels of battery charge are indicated by the segmented battery icon. When this icon turns red the battery is nearly exhausted and the batteries must be charged – see Charging Procedure.

The complete functionality is described on the help screen.

This is obtained by pressing the help button (?).

Help text exists for every screen viewed during the operation of the MicroLab.

You are recommended to make full use of the extensive Help screens provided.

Indication for Use

The MicroLab spirometer is intended, for prescription use only, to measure the maximal volume and flow of air that can be moved in and out of a patient's lungs and for pulse oximetry measurements. The system is intended for use with paediatric (4 to 17 years) and adult (18 to 99 years) patients in hospitals, physician offices, laboratories and occupational health testing environments.

Note: Pulse Oximetry is an optional measurement mode. Therefore the sensor for pulse oximetry measurements is not part of the standard MicroLab package and needs to be purchased separately.

Getting Started

When performing a spirometry test the recommended workflow is to enter the patient's details, or retrieve them from memory, perform the required test and then print and save the results.

Please ensure that the turbine transducer, is plugged in to either of the first two sockets on the right hand side of the instrument.

Patient Select

Search	Database Usage
ID:	1%
Name:	

ID	Name
123	Smith, Adam
234HI889	Lawson, Christopher
54tjn666889	Denton, James

Cancel Add Today

10:51

Touch the 'Patients' icon to enter the patient database.

The required patient may be selected from the stored patient list.

If the patient details have not been previously stored then touch **Add** to enter the new patient's details. The patient details may also be downloaded from the optional SPCS software.

Once selected, the patient's name will appear at the bottom of the screen.

Use the help button to obtain further information.

New Patient

Patient Details

ID

Last Name

First Name

Sex Male

Origin Caucasian

Height (cm)

1	2	3	4	5	6	7	8	9	0
q	w	e	r	t	y	u	i	o	p
↑	a	s	d	f	g	h	j	k	l
	z	x	c	v	b	n	m		
⬅	⬅	⬅	⬅	⬅	⬅	⬅	⬅	⬅	⬅

Cancel Finish

10:16

To add a patient to the database use the on screen keyboard to type a unique patient ID and then touch the enter key.


You will then be prompted for Last Name, First Name, Sex, Ethnic Origin, Height, Weight, date of Birth and Factor.

A factor can be applied when testing individuals of other ethnic origins who would not normally be tested against the countries set of predicted values eg ECCS predicted values are used but Non Caucasian individuals tested. The factor alters the predicted value set on volume indices by the percentage applied.

The following factors are recommended when using ECCS normal values:


Hong Kong Chinese	100%
Japanese American	89%
Polynesians	90%
North Indians and Pakistanis	90%
South Indians and those of African descent	87%

Ref: Lung Volumes and Forced Ventilatory Flows. P.H.Quanjer et al.
Eur Respir J, 1993, 6, Suppl. 16p5-40

Once all the patient details are added the patient is added to the database and the main menu is displayed with the patient name displayed at the bottom of the screen.

From the main menu select the required test, by touching the icon with the stylus. If the displayed patient is not required for testing then touch the patients name and options to change or remove the current patient will become available

If Relaxed Spirometry is selected then a volume/time graph will be displayed. Note that the unit may be customised to perform a relaxed Vital Capacity with tidal breathing or from a single expiration or single inspiration.

When a manoeuvre has been obtained touch 'Results' to view the indices, 'Again' to repeat the manoeuvre, 'Reject' to delete the manoeuvre or 'Done' to end the test.

Relaxed Result

Base
1*


Indice	Value	%Pred	[Min	Pred	
EVC	4.70				
IC	3.61				
TV	0.36				
ERV	1.09				
IRV	3.25				
FR	27.0				
Ti	1.16				
Te	0.97				
Ti/Ttot	54				

Graph
Done

Christopher Lawson
13:28

All the active indices are displayed for any of the manoeuvres selected together with an option to review the volume/ time curves. The active indices listed can be changed by using the customisation option.


Select 'Done' to proceed to the Spirometry Main Menu.

From this menu the results of the test may be viewed, saved, or printed and notes may be added.

It is also possible to proceed to a forced baseline spirometry test, or a post medication relaxed spirometry test.

Select 'Exit' when all the required functions have been used.


If forced spirometry is selected the default graph will be displayed. This may be changed by touching the arrows at the top of the screen. Flow/Volume, Volume/time or child incentive default displays may be selected using the customise option from the main menu.

When the spirometry manoeuvre has been completed options to repeat the test, reject the test, and view results will be available.


At the end of the test options to view results, save results, print results, and to add notes will be available from the spirometry main menu.

'Touch' the MVV icon to select this mode of testing and the display will instruct the patient to start breathing hard to commence the test.


It is recommended that the patient perform 3 tidal breathing manoeuvres prior to performing hard and fast rapid breathing (required for the MVV manoeuvre).

The patient should be instructed to tidally breathe. The tidal breaths are automatically detected prior to commencing the MVV manoeuvre.

Once tidal breathing is complete the display will change and an audible beep heard to instruct the patient to start rapid fast breathing. The start button should be touched using the stylus to start registering the MVV manoeuvre.

The current manoeuvre will be displayed in black, during the manoeuvre the breath rate (BR) will be displayed in green if the breath rate is acceptable (> 65 breaths per min) if the breath rate falls below this level it will be displayed in red to show the operator that the patient needs to be instructed to breath harder and faster during the manoeuvre. After 12 seconds of hard, fast and rapid breathing the display will show a green line indicating 12 seconds of the manoeuvre have elapsed – the patient should be encouraged to continue until the display changes to signify the end of the test. The MVV rate, the % variation between manoeuvres the breath rate and an ATS quality warning for the manoeuvre will be displayed.

Note: Patient efforts is acceptable when:

Patient made a maximum effort. Indicated to the user by the breath rate being displayed in green (>65 breaths per minute).


Lasts a full 12 seconds manoeuvre. Indicated to the user by a green line being displayed – the patient should ideally continue until the test is automatically terminated at 15s.

Interruption ie. Did not cough.

Once the test has finished the display will show current test (shown in black – if more than one manoeuvre has been performed then the best manoeuvre will also be displayed in blue) the MVV rate, the % variation between manoeuvres the breath rate and an ATS quality warning for the test session.

'Touch' Again to repeat the manoeuvre, Reject to reject the current manoeuvre, Results to display a list of indices, the values obtained, % predicted where applicable and also a quality statement concerning the test session.

To meet the ATS quality criteria for a Good blow the manoeuvre should last 15 seconds with a breath rate greater than 65 breaths per minute. The ATS reproducibility criterion is two manoeuvres with a Good Blow and the MVV variability between manoeuvres should not exceed 20%.

Note: The MVV test is an exhausting test. It should not be repeated without a rest period. Some elderly or ill people cannot repeat this test even after the rest period.

MVV Results

Base

1

2

3*

Indice	value	%pred.	[min	pred.	n
MVV	225	139	114	162	20
BR	75				
B.T.	15				
VT	0.00				
Te	0.00				
Ti	0.00				
TTOT	0.00				
Ti/Te	0.00				
Te/TTOT	0.00				
Ti/TTOT	0.00				

Good blow. ATS Criterion Met

Back

||||

Quick Patient

🕒 14:45

'Touch' Back to return to testing and the current manoeuvre.

Note: If the breathing rate is insufficient (less than 65 breaths per minute) then the BR value will be displayed in red – an MVV value will be calculated and a message displayed that the MVV result was extrapolated from a manoeuvre with a poor breath rate.

Once the number of manoeuvres has been completed and the test session has finished select Done and the results with selected indices will be displayed. Each manoeuvre will be numbered and the best manoeuvre highlighted with an asterix (*). Touch Graphs to view the graphs of the currently selected manoeuvre and best manoeuvre. Select Set best to manually select the best manoeuvre. Select Done to return to the main MVV menu.


Once testing is complete the MVV main menu will be displayed 'Touch' the appropriate icon to allow a Post 1 MVV test to be performed, View Results, Print Results, Add notes for the patients examination, Save the tests or Exit to return to the main spirometer menu.

Calibration Check

The spirometer is calibrated to read in litres at body temperature, Barometric Pressure Saturated with water vapour (BTPS).

The calibration should remain stable indefinitely, unless the transducer is physically damaged, and the unit should not require re-calibration. However, to ensure the correct functioning of the unit it is recommended that a calibration check is performed after the transducer has been removed for cleaning.


Connect a 3L syringe to the transducer with the minimum of adapters and empty it by pushing the handle fully in.

Note: It is recommended that the transducer is disinfected prior to a calibration check or a MicroGard filter is used during the procedure.

Select 'Calibration Check' from the main menu and then select 'check Calibration'.

Fill the syringe by pulling the handle at a constant rate until the end stop is reached and then immediately empty the syringe completely. Try to maintain a flow rate that keeps the trace within the grey bands on the display.


Press 'Reject' to retry the calibration check at the required flow rate.

Press 'Again' to repeat the calibration check at a low flow rate.

Press 'Again' to repeat the calibration check at a high flow rate.

When a calibration check at all three flow rates has been completed press 'Done' to view the calibration check report screen.

The calibration error for expiration and inspiration at each flow rate are displayed. The calibration error should be less than 3%. If a greater error is shown, repeat the procedure ensuring that the syringe is emptied and filled in a smooth manner without jerking the handle. If an error greater than 3% is still shown, inspect the turbine transducer and clean if necessary.

Customisation

The 'Customise' option from the main menu may be used to configure many of the features of your MicroLab and are divided into system, spirometry and MVV options.

System options allow you to configure the following:

- Language.
- Height and weight units.
- Date format.
- Date separator.
- Colour or monochrome printing (on external printer).
- Personalised printout heading.

Spirometry options allow you to configure the following:


- Relaxed spirometry mode (with or without tidal breathing).
- Predicted value sets.
- Predicted area or line display.
- Display default.
- Incentive display type.
- Printed graphs.
- Best test criteria
- Interpretation and Lung Age indication.
- Dyspnoea score and smoking status.
- Daily calibration reminder.
- Manual temperature adjustment.
- Indices selection.

MVV options allow you to configure the following:


- Choice of predicted values
- Display ambient temperature during MVV test
- Include graph of MVV manoeuvre in the final printout

Note: when the language is selected, the height and weight units, date format, and date separator will be automatically changed. However this automatic selection may be overridden manually.

Administration Mode

Administration mode allows the administrator to restrict the availability of functions to the user by disabling icons on the main menu. For example, after the unit has been configured to the administrator's requirements, disabling of the 'Customise' icon will prevent any further adjustment by the user. Similarly, disabling of the 'Database Management' icon will prevent the user from deleting any patient details or test results.


To enter administration mode turn the unit on whilst holding down the help key. The default access code is 0000. Type this number in using the on-screen keyboard. A number of functions are now available.

Touch the 'Change Access Code' icon in order to enter your personal access code and thereby restrict the use of administration to authorised personnel.

Press the help button to obtain a full description of the functions.

Paper Loading

To load a new roll of thermal paper lift the paper cover using the side levers, place the paper into the compartment as shown and close the cover firmly. It is recommended that only CareFusion thermal printer paper (Cat No. 36-PSA1600) is used with the MicroLab to avoid damage to the thermal printer head.

To tear off the paper pull the paper towards you and to the right as shown below:

External Printing

When printing, the MicroLab will default to an external USB printer, if connected.

Using the USB cable provided connect the mini A/B socket on the right hand side of the MicroLab to the input socket on the printer. For a list of compatible printers refer to the web site www.carefusion.com/micromedical or contact your local dealer.

It is recommended that during printing the batteries are on charge with the unit connected to the mains adapter.

Note: Keep the printer out of reach of the patient at all times.

Note: Disconnect the printer during live measurements.

CAUTION: Connect the MicroLab only to printers and computers that comply with EN 60601-1 – ‘Medical Electrical Equipment’ Part 1: General requirements for basic safety and essential performance.

Switching Off

The unit is switched off by pressing the On/Off button.

The unit can be disconnected from the mains by unplugging the charger from the mains socket or unplugging the USB cable if no mains supply is connected.

Charging Procedure

The MicroLab’s internal batteries are discharged when shipped from the factory and should be fully charged on first use. Plug the AC adapter into the mains supply and plug the adapter output plug into the power input socket on the right hand side of the instrument. The orange charging light next to the power input socket will flash to indicate charging and will turn on constantly to indicate full charge.

The batteries will take approximately 4 hours to become fully charged.

Note: Use only the AC adapter supplied. Use of any other type may cause permanent damage to the MicroLab and cause a fire or electrical hazard. Do not plug in and remove the power lead from the AC adapter repeatedly.

Please Note: Dispose of the waste battery in accordance with EU Waste Battery Regulations.

PC connection using SPCS

SPCS is an easy to use PC based windows application that interfaces to the MicroLab via a USB or serial port. It incorporates a database into which patient details can be entered and downloaded to the MicroLab or test results may be uploaded from the MicroLab to the PC.

Using SPCS and the MicroLab, live blows can be performed with the PC directly controlling the operation of the MicroLab.

The results and graphs produced are displayed directly on the PC screen.

The spirometer is connected from the serial or USB port on the PC, to the port on the right hand side of the instrument using the USB or serial cable provided with SPCS.

Note: Keep the PC and Monitor out of reach of the patient at all times.

It is recommended that whilst the unit is connected to a computer the mains adapter is used.

Looking after your Spirometer

Please observe the following precautions:

- Do not touch the screen with fingers. Use only the stylus provided.
- Use only a damp, lint free, cloth to clean the screen.
- Do not keep the spirometer in a damp place or expose it to extremes of temperature.
- Do not direct the transducer holder towards a strong light source whilst operating the spirometer.
- Check the AC charger for compatibility with local power rating.

Cleaning Instructions

Disinfection of contaminated parts is only effective after having them carefully pre-cleaned. CareFusion recommend the tested solution of PeraSafe sterilizing powder (36-SSC5000A) for pre-cleaning and disinfection. If a different solution is used, please follow the given manufacturer's instructions.

The device must not be wiped with any aqueous solutions and must not be exposed to solvents e.g. alcohol, chloride solutions as there are electronic components inside that will be permanently damaged.

CAUTION: Switch off the device and always unplug your MicroLab before cleaning.

External Surfaces of the Spirometer

CAUTION: Do not attempt to wash or immerse the MicroLab transducer housing in water or cleaning fluid, as there are electronic components inside that will be permanently damaged

The external housing of the spirometer may be wiped with sterile wipes or a damp cloth that has been immersed in a cold sterilising solution, when required.

CAUTION: Do not wipe the touch screen.

Cleaning accessories

With the use of a filter MicroGard® IIB (V-892381), MicroGard® IIC (V-892384) or a One-way Safety mouthpiece (36-SST1250) for each patient, cleaning of the transducer is recommended once a month.

When using the disposable cardboard mouthpiece (adult: 36-PSA1000, pediatric: 36-PSA1200) without a filter and under the prerequisite that the patient was instructed only to exhale into the transducer, the following parts have to be cleaned once a day: adult/ pediatric adapter, transducer.

With any other use as described above all contaminated parts must be disinfected between patients.


Important note: Used single patient nose clip pads, mouthpieces, MicroGard® IIB filter and MicroGard® IIC filter must be disposed off immediately after the use.

If there are changes on the material surfaces (cracks, brittleness) the respective parts must be disposed off.

Cleaning the Transducer

The transducer requires no routine maintenance or servicing. However, if you wish to disinfect or clean the transducer it may be removed by means of the following procedure:

1. Rotating the turbine transducer anti-clockwise until the locating pip lines up with the small rectangular cut-out in the housing as shown below.
2. Gently pull the transducer away from the housing.
3. The transducer may now be immersed in warm soapy water for routine cleaning or immersed in cold disinfecting solutions e.g. Perasafe (Cat No. 36-SSC5000A) for a maximum of 10 minutes (Alcohol and chloride solutions should be avoided).
4. After cleaning/disinfecting, the transducer should be immersed briefly in distilled water and left to dry.
5. Re-assemble the mouthpiece holder

CAUTION: Do not attempt to wash or immerse the MicroLab transducer housing in water or cleaning fluid, as there are electronic components inside that will be permanently damaged.

Servicing

There is no routine maintenance required for the MicroLab and there are no user serviceable parts in this spirometer. Please return the unit to CareFusion or an authorised agent, if servicing is required.

If your spirometer requires service or repair please see page 50 for contact details.

Trouble Shooting Information

Should you encounter problems operating your MicroLab meter, please consult the table below:

Problem	Possible Cause	Solution
display 'freezes' and the unit does not respond to any key presses	Multiple icons have been selected or accidentally pressed	hold the on/off button down for 10 seconds until the unit switches off and then turn on again
Does not register a blow	Head assembly or cable broken	Replacement of head assembly or return the unit for servicing
No display present	Charger not connected or battery is exhausted	Connect charger to the mains and leave unit to fully charge, or return the unit for servicing
Blows are inverted on the display	Head assembly or cable broken	Replacement of head assembly or return the unit for servicing
Blows tracing ends abruptly although patient still exhaling	Turbine sticking	Clean turbine in warm soapy water or sterilising solution, if problem persists a replacement turbine may be required
Battery does not hold charge	Exhausted battery	Return the unit for servicing
	Mains charger faulty	Replace the mains charger

Paper not printing	Check paper is housed in printer compartment correctly	Refer to section “paper loading” of this manual, and follow the instructions.
	Incorrect thermal paper being used	Ensure you are using CareFusion recommended thermal paper (See Consumables /Accessories section)
Stylus does not register icons on the display	Touch screen display requires calibration	Select the calibration check icon and choose touch screen, and follow the instructions.
Icons missing from the display	Icon has been de-selected	Hold stylus on to the blue area of the display, a list will appear, ensure required icon is selected.
Calibration procedure failed or cannot be completed	Turbine may be faulty	Repeat calibration procedure, if problem still persists, replace turbine or return unit for servicing
	Turbine not fitted tightly to calibration syringe	Ensure the syringe is fitted to the turbine using the adapter supplied part # 861427. (supplied with the syringe)
	Calibration syringe does not have an inspiratory seal or seal is leaking.	Ensure you are using a CareFusion syringe.
	Shaft of the syringe being pushed down	The syringe should be emptied and filled with one smooth stroke, avoid pushing down on the shaft or banging at the end of each manoeuver.

Safety Designation per IEC 60601-1

Type of protection against electrical shock	Internally powered Equipment and Class I
Degree of protection against electrical shock	Type B applied part
Power Equipment	Battery type: NiMH battery pack, 8.4V, 1100mAh
Degree of Electrical connection between equipment and Patient	Equipment designed as non-electrical connection to the patient.
Degree of mobility	Transportable
Mode of operation	Continuous

Classifications according to IEC 60601-1

MicroLab	Applied part, type B
Volume Transducer	Applied part, type B
Option SPO2 Finger Probe (Nonin IPOD® Integrated Oximetry Device)	Applied part, type BF

WARNING: No modification of this equipment is allowed.

Note: When you connect other equipment to the unit, always make sure that the whole combination complies with the international safety standard IEC 60601-1 for medical electrical systems. During measurements, connect the MicroLab only to printers and computers that comply with IEC/EN 60601-1 / UL 60601-1.

WARNING: The user must not touch any voltage-carrying parts and the patient at the same time.

During database upload the MicroLab may be connected to a computer that complies with EN 60950 – ‘Information technology equipment - Safety - Part 1: General requirements’.

IMPORTANT: Only use the mains adapter supplied (36-PSU1012 12V DC 2.5A) The adapter contains a transformer. Do not cut off the adapter to replace it with another plug as this causes a hazardous situation.

- The adapter transforms the mains voltage (100-240 Volts) to a safe voltage (12V DC)
- Make sure the adapter does not get wet
- Do not use a damaged adapter
- Always unplug your MicroLab before cleaning

WARNING: Do not connect devices that are not specified as part of the system.

Note: If an MPSO (Multiple Portable Socket Outlet) is used with the system, the maximum permitted load should not be exceeded. Do not connect electrical equipment that has not been supplied as part of the system.

Electromagnetic Compatibility (EMC) to EN60601-1:2007.

WARNING: use of portable phones or other radio frequency (RF) emitting equipment near the system may cause unexpected or adverse operation

The MicroLab has been tested to EN60601-1-2:2007, regarding its ability to operate in an environment containing other electrical/electronic equipment (including other medical devices).

The purpose of this testing is to ensure that the MicroLab is not likely to adversely affect the normal operation of other such equipment and that other such equipment is not likely to adversely affect the normal operation of the MicroLab.

Despite the testing of the MicroLab that has been undertaken, normal operation of the MicroLab can be affected by other electrical/electronic equipment and portable and mobile RF communications equipment.

As the MicroLab is medical equipment, special precautions are needed regarding EMC (electromagnetic compatibility).

It is important that the MicroLab is configured and installed/put into service, in accordance with the instructions/guidance provided herein and is used only in the configuration as supplied.

Changes or modifications to the MicroLab may result in increased emissions or decreased immunity of the MicroLab in relation to EMC performance.

The MicroLab should be used only with the accessories (USB cables, mains adapter and turbine transducer) supplied (which are referenced in the accessories section of this manual). None of the MicroLab cables should be extended in length by the user.

If any cables are extended by the user or non approved accessories are used, this may result in an increased level of emissions or decreased level of immunity, in relation to the MicroLab's EMC. None of the MicroLab accessories should be used with other devices, as this may result in an increased level of emissions or decreased level of immunity, in relation to the other devices's EMC.

The MicroLab has an essential performance – the product should continue to operate correctly. In the unlikely event of a Fast Transient / ESD event occurring, the device should be reset and located away from the source of interference.

WARNING: The MicroLab should not be used adjacent to or stacked with other equipment. If adjacent or stacked use with other equipment is necessary, the MicroLab and the other equipment should be observed / monitored, to verify normal operation in the configuration in which it will be used.

Guidance and Manufacturer's Declaration – Electromagnetic Emissions

The MicroLab is intended for use in the electromagnetic environment specified below. The customer or the user of the MicroLab should assure that it is used in such an environment.

Emission Test	Compliance	Electromagnetic Environment - Guidance
RF emissions CISPR 11	Group 1	The MicroLab uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Group B	The MicroLab is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes
Harmonic emissions IEC61000-3-2	Class A	
Voltage fluctuations / flicker emissions IEC61000-3-3	Complies	

Guidance and Manufacturer's Declaration – Electromagnetic Immunity

The MicroLab is intended for use in the electromagnetic environment specified below. The customer or the user of the MicroLab should assure that it is used in such an environment.

Immunity Test	IEC 60601 test level	Compliance Level	Electromagnetic Environment -Guidance
Electrostatic discharge (ESD) IEC61000-4-2	± 6 kV contact ± 8 kV air	± 6 kV contact ± 8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient / burst IEC61000-4-4	± 2 kV for power supply lines ± 1 kV for input / output lines	± 2 kV for power supply lines Input/output line tests not applicable (<3m)	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC61000-4-5	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	Mains power quality should be that of a typical commercial or hospital environment.

Voltage dips, short interruptions and voltage variations on power supply input lines IEC61000-4-11	<5% U_T (>95 % dip in U_T) For 0.5 cycle 40% U_T (60 % dip in U_T) for 5 cycles 70 % U_T (30 % dip in U_T) for 25 cycles <5% U_T (>95 % dip in U_T) For 5 s	<5% U_T (>95 % dip in U_T) For 0.5 cycle 40% U_T (60 % dip in U_T) for 5 cycles 70 % U_T (30 % dip in U_T) for 25 cycles <5% U_T (>95 % dip in U_T) For 5 s	Mains power quality should be that of a typical commercial or hospital environment. If the user of the MicroLab requires continued operation during power mains interruptions, it is recommended that the MicroLab be powered from an uninterruptable power supply or a battery
Power frequency (50/60Hz) Magnetic field IEC61000-4-8	3 A/m	3 A/m	If incorrect operation occurs, it may be necessary to position the MicroLab further from sources of power frequency magnetic fields or to install magnetic shielding. The power frequency magnetic field should be measured in the intended installation location to assure that it is sufficiently low.
NOTE U_T is the a.c. mains voltage prior to application of the test level.			

Guidance and Manufacturer's Declaration – Electromagnetic Immunity

The MicroLab is intended for use in the electromagnetic environment specified below. The customer or the user of the MicroLab should assure that it is used in such an environment.

Immunity Test	IEC 60601 test level	Compliance Level	Electromagnetic Environment - Guidance
Conducted RF IEC61000-4-6	3 Vrms 150 kHz to 80 MHz	3 Vrms	Portable and mobile RF communications equipment should be used no closer to any part of the MicroLab, including any cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance (d) $d = 1.2\sqrt{P}$

Radiated RF IEC61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m	$d = 1.2\sqrt{P}$ 80 MHz to 800 MHz $d = 2.3\sqrt{P}$ 800 MHz to 2.5 GHz Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m). Fields strengths from fixed RF transmitters, as determined by an electromagnetic site survey, ^a should be less than the compliance level in each frequency range. ^b Interference may occur in the vicinity of equipment marked with the following symbol: 
NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies. NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			
^a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the MicroLab is used exceeds the applicable RF compliance level above, the MicroLab should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orientating or relocating the MicroLab.			
b. Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.			

Recommended separation distances between portable and mobile RF communications equipment and the MicroLab

The MicroLab is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the MicroLab can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the MicroLab as recommended below, according to the maximum output power of the communications equipment.

Rated Maximum Output Power of Transmitter in Watts (W)	Separation Distance in Meters (m) according to Frequency of Transmitter		
	150 KHz to 80 MHz $d = 1.2\sqrt{P}$	80 MHz to 800 MHz $d = 1.2\sqrt{P}$	800 MHz to 2.5 GHz $d = 2.3\sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23.3

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures objects and people.

Symbols

Type B device

In accordance with Directive 93/42/EEC

Disposal in compliance with WEEE

Consult the instructions for use

Caution: Consult the accompanying documents

Manufacturer

Date of Manufacture

Serial number

Direct current

Single patient use

Federal U.S. law restricts this device to sale by or on the order of a physician. (Rx only)

Nationally Recognised Test Laboratory (NRTL) officially recognized by the American Occupational Safety and Health Administration (OSHA) for OSHA Electrical safety and compliance.

Specification of the MicroLab

General:

Storage:	>2000 tests including Flow/Volume loops and Volume/Time curves
Printer Output:	Compatible Hewlett Packard USB printers (Specified models only).
Display:	Colour 1/4VGA LCD.
Power supply:	Input 100 to 240V, 47 to 63Hz. Output 12V 2.5A (Class 1) Type: MENB1030A1200F03
Battery Pack:	Rechargeable NiMH 8.4V 1A-hours
Battery life:	Approximately 30 hours with a fully charged new battery
Dimensions:	35x255x120 mm. Transducer 50x60x90mm.
Weight:	630g.
Operating Temperature:	0 to +40 °C
Operating Humidity:	30% to 90% RH

Transport and Storage -20 to + 70 °C

Temperature:

Transport and Storage 10% to 90% RH

Humidity:

Spirometry Measurements:

Relaxed Expiratory Vital Capacity (VC)

Forced Expired Volume in 0.75 seconds (FEV_{0.75})

Forced Expired Volume in 1 second (FEV₁)

Forced Expired Volume in 3 second (FEV₃)

Forced Expired Volume in 6 seconds (FEV₆)

Forced Vital Capacity (FVC)

Peak Expiratory Flow Rate (PEF)

FEV_{0.75} as a percentage of VC (FEV_{0.75}/VC)

FEV_{0.75} as a percentage of FVC (FEV_{0.75}/FVC)

FEV₁ as a percentage of VC (FEV₁/VC)

FEV₁ as a percentage of FVC (FEV₁/FVC)

FEV₃ as a percentage of VC (FEV₃/VC)

FEV₃ as a percentage of FVC (FEV₃/FVC)

FEV_{0.75} as a percentage of FEV₆ (FEV_{0.75}/FEV₆)

FEV₁ as a percentage of FEV₆ (FEV₁/FEV₆)

Maximum Expired Flow at 75% of FVC remaining (MEF₇₅)

Maximum Expired Flow at 50% of FVC remaining (MEF₅₀)

Maximum Expired Flow at 25% of FVC remaining (MEF₂₅)

Mean Mid-Expiratory Flow Rate (MMEF)

Forced expiratory flow at 50% of volume as a percentage of VC (FEF₅₀/VC)

Forced expiratory flow at 50% of volume as a percentage of FVC (FEF₅₀/FVC)

Maximal voluntary ventilation indicated (MVV (ind))

Forced inspired volume in 1 second (FIV₁)

Forced inspiratory Vital Capacity (FIVC)

Peak Inspiratory Flow Rate (PIF)

FIV₁ as a percentage of FIVC (FIV₁/FIVC)

Forced inspiratory flow at 25% of inhaled volume (FIF25)	
Forced inspiratory flow at 50% of inhaled volume (FIF50)	
Forced inspiratory flow at 75% of inhaled volume (FIF75)	
Forced expiratory flow at 50% of volume as a percentage of FIF50 (FEF50/FIF50)	
The time taken between 25% and 75% of the forced expired volume (MET2575)	
Forced Expiratory Time (FET)	
Tidal Volume (TV)	
Expiratory reserve volume (ERV)	
Inspiratory reserve volume (IRV)	
Inspiratory capacity (IC)	
Expiratory Relaxed vital capacity (EVC)	
Inspiratory vital capacity (IVC)	
Breathing frequency rate (FR)	
Inspiratory time (Ti)	
Expiratory time (Te)	
Ti as a % of total breath time (Ti/Ttot)	
Tidal volume as a % of Ti (TV/Ti)	
Breath Rate	BR
Breathing Time	B.T
Volume Tidal	VT
Expiratory Time - average time of expiration per breath in seconds	Te
Inspiratory Time - average time of inspiration per breath in seconds	Ti
Total Tidal Breath Time in Seconds	TTOT = Ti + Te
Ratio of Average Expiratory and Inspiratory Breaths	Ti/Te
Average Time of Expiration per Breath as a ratio to the Total Tidal Breath Time	Ti/TTOT

Tests per subject:

5 VC manoeuvre
8 FVC manoeuvres

Predicted Values:

Various – depends upon national preference

Transducer:

CareFusion Bi-Directional Digital Volume.

Resolution:

10ml volume 0.03l/s flow

Accuracy:

+/-3%. To ATS recommendations –
Standardisation of spirometry 1994
update for flows and volumes.

Consumables/Accessories

Cat. No.	Description
36-PSA1600	Thermal Printer Paper (pack of 5 rolls)
36-PSA1000	Adult Disposable Mouthpieces (500 per box)
36-SST1000	One-way Safety Mouthpieces (500 per box)
36-SST1250	One-way Safety Mouthpieces (250 per box)
36-PSA1200	Paediatric Disposable Mouthpieces (250 per box)
36-PSA1100	Paediatric Adaptor
36-MFA1010	MicroGard Filter Adapter (Pack of 10)
36-MFA1050	MicroGard Filter Adapter (Pack of 50)
V-892381	MicroGard® IIB with integrated mouthpiece (50 per box)
V-892382	MicroGard® IIB with integrated mouthpiece (100 per box)
V-892384	MicroGard® IIC (50 per box)
V-892385	MicroGard® IIC (100 per box)
36-SSC5000A	PeraSafe Sterilising Powder 81g (to make up 5 litres of solution)
36-VOL2104	Nose clip (pack of 5)
892121	Nose clip pads "foam material", 100 pcs.
36-SM2125	3 Litre Calibration Syringe
V-861449	Silicone adapter "oval" for MicroGard® II B
861427	Calibration syringe adapter
36-SPC1000	Spirometry PC Software
36-CAB1048	USB Lead
36-PSU1012	Mains Adapter
36-TDX1048	TurbineTransducer

For further information or to place an order for Disposables/
Supporting products please contact CAREFUSION, your local
distributor or view our website www.carefusion.com/micromedical

NOTE: USE ONLY CAREFUSION ACCESSORIES

Note: Information in this manual is subject to change without notice and does not represent commitment on the part of CareFusion 232 UK Ltd. The software may be used or copied only in accordance with the terms of that agreement. No part of this manual may be reproduced or transmitted in any form or by any means, electronic or mechanical, including photocopying and recording for any purpose without the written permission of CareFusion 232 UK Ltd.

Customer contact information

UK Customers only

For all Sales Order processing for products, training and spare parts, Service and Technical Support enquiries, please contact the following:

CareFusion UK 232 Ltd
UK Customer Service & Support
The Crescent
Jays Close
Basingstoke
RG22 4BS

Customer Service Sales Enquiries:

Telephone: 01256 388550

Email: micro.uksales@carefusion.com

Factory Repair and Administration Enquiries:

Telephone: 01256 388552

Email: micro.ukservice@carefusion.com

Technical Support Enquiries:

Telephone: 01256 388551

Email: support.rt.eu@carefusion.com

International customers only

For all Sales Order processing for products and Spare parts, Service and Technical Support enquiries, please contact the following:

CareFusion Germany 234 GmbH
Customer Service & Support International
Leibnizstrasse 7
97204 Hoechberg
Germany

Customer Service Sales Enquiries:

Telephone: 0049 931 4972 670

Email: micro.international-sales@carefusion.com

Factory Repair and Administration Enquiries:

Telephone: 0049 931 4972 867

Email: support.admin.eu@carefusion.com

Technical Support Enquiries:

Telephone: 0049 931 4972 127

Email: support.rt.eu@carefusion.com

MicroLab

Operating manual

Federal (USA) law restricts this device to sale by or on the order of a physician or licensed practitioner.

CareFusion,
22745 Savi Ranch Parkway,
Yorba Linda
CA 92887-4668
USA

MicroLab

Operating manual

II. Introduction – English - USA

The MicroLab is mains/battery operated portable spirometer with the unique combination of ease of use and sophistication. Ease of use is assured through the use of context sensitive help screens, accessed at a touch of a button, that explain every MicroLab feature.

The MicroLab uses the CareFusion Digital Volume Transducer, an extremely stable form of volume transducer, which measures expired air directly at B.T.P.S (Body Temperature and Pressure with Saturated water vapor) thus avoiding the inaccuracies of temperature corrections. The transducer is insensitive to the effects of condensation and temperature and avoids the need for individual calibration prior to performing a test.

Test results may be uploaded to a PC using the optional Spirometry PC Software (SPCS) and patient details may be downloaded to the MicroLab.

Stored data may be printed to the integral thermal printer or to an external HP printer using the USB cable supplied.

The MicroLab utilizes a **single patient use** disposable mouthpiece that must be disposed of after use.

The MicroLab provides a suggested interpretation **that must be supported by clinical judgment.**

Contraindications

- Acute disorders affecting test performance (e.g. Vomiting, nausea, vertigo)
- Recent eye surgery (increases in intraocular pressure during spirometry)
- Oral or facial pain exacerbated by a mouthpiece
- recent myocardial infarction
- post- operative thoracic surgery patients
- hyperventilation syndrome

Note: Extensive exhalation might lead to syncope.

Warnings and Cautions

The following terms are used as follows in this manual

CAUTION: Possibility of injury or serious damage.

WARNING: Conditions or practices that could result in personal injury.

Note: Important information for avoiding damage to the instrument or facilitating operation of the instrument.

Note: Patients below 4 years of age may struggle to perform spirometry correctly and reproducibly.

Note: The device should only be used by trained and qualified personnel.

CAUTION: Read the manual before use.

WARNING: The instrument is not suitable for use in the presence of explosive or flammable gases, flammable anesthetic mixtures or in oxygen rich environments.

CAUTION: Mouthpieces are single patient use. If used on more than one patient there is a risk of cross-infection. Repeat use may degrade materials and lead to an incorrect measurement.

CAUTION: Pulmonary filters are single patient use. If used on more than one patient there is a risk of cross-infection. Repeat use may increase air resistance and lead to an incorrect measurement.

PLEASE NOTE: The product you have purchased should not be disposed of as unsorted waste. Please utilize your local  WEEE collection facilities for the disposal of this product.

PLEASE NOTE: Degree of protection against Ingress of Water is IPX0.

WARNING: To avoid risk of electric shock, this equipment must only be connected to a supply mains with protective earth.


WARNING: Do not connect devices that are not specified as part of the system.

Overview

When the unit is turned on for the first time the system customize menu is displayed. Touch the name of the required region with the stylus stored on the left hand side of the unit and then touch 'Finish'. This will set the regional defaults including language, if available.

If this screen is not displayed then touch the 'Customize' icon on the main menu followed by the 'System' icon.


The MicroLab uses a touch screen with icons representing each function available. A stylus, housed in the left hand side of the unit, is provided for icon screen activation.

Touch the displayed time to adjust time and date. Touch the toolbox icon to adjust volume and brightness.

Unused icons may be disabled by touching the blue background and selecting from the list displayed.

Four levels of battery charge are indicated by the segmented battery icon. When this icon turns red the battery is nearly exhausted and the batteries must be charged – see Charging Procedure.

The complete functionality is described on the help screen.

This is obtained by pressing the help button (?).

Help text exists for every screen viewed during the operation of the MicroLab.


You are recommended to make full use of the extensive Help screens provided.

Indication for Use

The MicroLab spirometer is intended, for prescription use only, to measure the maximal volume and flow of air that can be moved in and out of a patient's lungs and for pulse oximetry measurements. The system is intended for use with paediatric (4 to 17 years) and adult (18 to 99 years) patients in hospitals, physician offices, laboratories and occupational health testing environments.

Note: Pulse Oximetry is an optional measurement mode. Therefore the sensor for pulse oximetry measurements is not part of the standard MicroLab package and needs to be purchased separately.

Getting Started

When performing a spirometry test the recommended workflow is to enter the patient's details, or retrieve them from memory, perform the required test and then print and save the results.

Please ensure that the turbine transducer is plugged in to either of the first two sockets on the right hand side of the instrument.

Patient Select	
Search	Database Usage
ID:	1%
Name:	
ID	Name
123	Smith, Adam
234HI889	Lawson, Christopher
54tjn666889	Denton, James
Cancel	Add Today

10:51

Touch the 'Patients' icon to enter the patient database.

The required patient may be selected from the stored patient list.

If the patient details have not been previously stored then touch **Add** to enter the new patient's details. The patient details may also be downloaded from the optional SPCS software.

Once selected, the patient's name will appear at the bottom of the screen.

Use the help button to obtain further information.

New Patient	
Patient Details	
ID	
Last Name	
First Name	
Sex	Male
Origin	Caucasian
Height (cm)	
1 2 3 4 5 6 7 8 9 0	
q w e r t y u i o p	
a s d f g h j k l	
z x c v b n m	
Cancel	Finish

10:16


To add a patient to the database use the on screen keyboard to type a unique patient ID and then touch the enter key.

You will then be prompted for Last Name, First Name, Sex, Ethnic Origin, Height, Weight and date of Birth.

Once all the patient details are added the patient is added to the database and the main menu is displayed with the patient name displayed at the bottom of the screen.

From the main menu select the required test, by touching the icon with the stylus. If the displayed patient is not required for testing then touch the patients name and options to change or remove the current patient will become available

If Relaxed Spirometry is selected then a volume/time graph will be displayed. Note that the unit may be customized to perform a relaxed Vital Capacity with tidal breathing or from a single expiration or single inspiration.

When a maneuver has been obtained touch 'Results' to view the indices, 'Again' to repeat the maneuver, 'Reject' to delete the maneuver or 'Done' to end the test.

Relaxed Result

Base
1*


Indice	Value	%Pred	[Min	Pred	↑
EVC	4.70				
IC	3.61				
TV	0.36				
ERV	1.09				
IRV	3.25				
FR	27.0				
Ti	1.16				
Te	0.97				
Ti/Ttot	54				

Graph
Done

Christopher Lawson 13:28

All the active indices are displayed for any of the maneuvers selected together with an option to review the volume/ time curves. The active indices listed can be changed by using the customization option.


Select 'Done' to proceed to the Spirometry Main Menu.

From this menu the results of the test may be viewed, saved, or printed and notes may be added.

It is also possible to proceed to a forced baseline spirometry test, or a post medication relaxed spirometry test.

Select 'Exit' when all the required functions have been used.


If forced spirometry is selected the default graph will be displayed. This may be changed by touching the arrows at the top of the screen. Flow/Volume, Volume/time or child incentive default displays may be selected using the customize option from the main menu.

When the spirometry maneuver has been completed options to repeat the test, reject the test, and view results will be available.


At the end of the test options to view results, save results, print results, and to add notes will be available from the spirometry main menu.

'Touch' the MVV icon to select this mode of testing and the display will instruct the patient to start breathing hard to commence the test.


It is recommended that the patient perform 3 tidal breathing maneuvers prior to performing hard and fast rapid breathing (required for the MVV maneuver).

The patient should be instructed to tidally breathe. The tidal breaths are automatically detected prior to commencing the MVV maneuver.

Once tidal breathing is complete the display will change and an audible beep heard to instruct the patient to start rapid fast breathing. The start button should be touched using the stylus to start registering the MVV maneuver.

The current maneuver will be displayed in black, during the maneuver the breath rate (BR) will be displayed in green if the breath rate is acceptable (> 65 breaths per min) if the breath rate falls below this level it will be displayed in red to show the operator that the patient needs to be instructed to breath harder and faster during the maneuver. After 12 seconds of hard, fast and rapid breathing the display will show a green line indicating 12 seconds of the maneuver have elapsed – the patient should be encouraged to continue until the display changes to signify the end of the test. The MVV rate, the % variation between maneuvers the breath rate and an ATS quality warning for the maneuver will be displayed.

Note: Patient efforts is acceptable when:

Patient made a maximum effort. Indicated to the user by the breath rate being displayed in green (>65 breaths per minute).


Lasted a full 12 seconds maneuver. Indicated to the user by a green line being displayed – the patient should ideally continue until the test is automatically terminated at 15s.

Interruption i.e. Did not cough.

Once the test has finished the display will show current test (shown in black – if more than one maneuver has been performed then the best maneuver will also be displayed in blue) the MVV rate, the % variation between maneuvers the breath rate and an ATS quality warning for the test session.

'Touch' Again to repeat the maneuver, Reject to reject the current maneuver, Results to display a list of indices, the values obtained, % predicted where applicable and also a quality statement concerning the test session.

To meet the ATS quality criteria for a Good blow the maneuver should last 15 seconds with a breath rate greater than 65 breaths per minute. The ATS reproducibility criterion is two maneuvers with a Good Blow and the MVV variability between maneuvers should not exceed 20%.

Note: The MVV test is an exhausting test. It should not be repeated without a rest period. Some elderly or ill people cannot repeat this test even after the rest period.

MVV Results

Base

1 2 3*

Indice	value	%pred.	(min	pred.	n
MVV	225	139	114	162	20
BR	75				
B.T.	15				
VT	0.00				
Te	0.00				
Ti	0.00				
TTOT	0.00				
Ti/Te	0.00				
Te/TTOT	0.00				
Ti/TTOT	0.00				

Good blow. ATS Criterion Met

Back

Quick Patient 14:45

'Touch' Back to return to testing and the current maneuver.

Note: If the breathing rate is insufficient (less than 65 breaths per minute) then the BR value will be displayed in red – an MVV value will be calculated and a message displayed that the MVV result was extrapolated from a maneuver with a poor breath rate.

Once the number of maneuvers has been completed and the test session has finished select Done and the results with selected indices will be displayed. Each maneuver will be numbered and the best maneuver highlighted with an asterix (*). Touch Graphs to view the graphs of the currently selected maneuver and best maneuver. Select Set best to manually select the best maneuver. Select Done to return to the main MVV menu.


Once testing is complete the MVV main menu will be displayed 'Touch' the appropriate icon to allow a Post 1 MVV test to be performed, View Results, Print Results, Add notes for the patients examination, Save the tests or Exit to return to the main spirometer menu.

Calibration Check

The spirometer is calibrated to read in liters at body temperature, Barometric Pressure Saturated with water vapor (BTPS).

The calibration should remain stable indefinitely, unless the transducer is physically damaged, and the unit should not require re-calibration. However, to ensure the correct functioning of the unit it is recommended that a calibration check is performed after the transducer has been removed for cleaning.


Connect a 3L syringe to the transducer with the minimum of adapters and empty it by pushing the handle fully in.

Note: It is recommended that the transducer is disinfected prior to a calibration check or a MicroGard filter is used during the procedure.

Select 'Calibration Check' from the main menu and then select 'check Calibration'.

Fill the syringe by pulling the handle at a constant rate until the end stop is reached and then immediately empty the syringe completely.


Press 'Reject' to retry the calibration check at the required flow rate.

Press 'Again' to repeat the calibration check at a low flow rate.

Press 'Again' to repeat the calibration check at a high flow rate.

When a calibration check at all three flow rates has been completed press 'Done' to view the calibration check report screen.

The calibration error for expiration and inspiration at each flow rate are displayed. The calibration error should be less than 3%. If a greater error is shown, repeat the procedure ensuring that the syringe is emptied and filled in a smooth manner without jerking the handle. If an error greater than 3% is still shown, inspect the turbine transducer and clean if necessary.

Customization

The 'Customize' option from the main menu may be used to configure many of the features of your MicroLab and are divided into system, spirometry and MVV options.

System options allow you to configure the following:

- Language.
- Height and weight units.
- Date format.
- Date separator.
- Color or monochrome printing (on external printer).
- Personalized printout heading.

Spirometry options allow you to configure the following:

- Relaxed spirometry mode (with or without tidal breathing).
- Predicted value sets.
- Predicted area or line display.
- Display default.
- Incentive display type.
- Printed graphs.
- Best test criteria.
- Interpretation and Lung Age indication.
- Dyspnoea score and smoking status.
- Calibration reminder and manual temperature adjustment.
- Indices selection.

MVV options allow you to configure the following:

- Choice of predicted values
- Display ambient temperature during MVV test
- Include graph of MVV maneuver in the final printout

Note: when the language is selected, the height and weight units, date format, and date separator will be automatically changed. However this automatic selection may be overridden manually.

Administration Mode


Administration mode allows the administrator to restrict the availability of functions to the user by disabling icons on the main menu. For example, after the unit has been configured to the administrator's requirements, disabling of the 'Customize' icon will prevent any further adjustment by the user. Similarly, disabling of the 'Database Management' icon will prevent the user from deleting any patient details or test results.

To enter administration mode turn the unit on whilst holding down the help key. The default access code is 0000. Type this number in using the on-screen keyboard. A number of functions are now available. Press the help button to obtain a full description of the functions.

Paper Loading

To load a new roll of thermal paper lift the paper cover using the side levers, place the paper into the compartment as shown and close the cover firmly. It is recommended that only CareFusion thermal printer paper is used with the MicroLab to avoid damage to the thermal printer head.

To tear off the paper pull the paper towards you and to the right as shown below:

External Printing

When printing, the MicroLab will default to an external USB printer, if connected.

Using the USB cable provided connect the mini A/B socket on the right hand side of the MicroLab to the input socket on the printer. For a list of compatible printers contact your local dealer.

It is recommended that during printing the batteries are on charge with the unit connected to the mains adapter.

Note: Keep the printer out of reach of the patient at all times.

Note: Disconnect the printer during live measurements.

CAUTION: Connect the MicroLab only to printers and computers that comply with EN 60601-1 – ‘Medical Electrical Equipment’ Part 1: General requirements for basic safety and essential performance.

Switching Off

The unit is switched off by pressing the On/Off button.

The unit can be disconnected from the mains by unplugging the charger from the mains socket or unplugging the USB cable if no mains supply is connected.

Charging Procedure

The MicroLab’s internal batteries are discharged when shipped from the factory and should be fully charged on first use. Plug the AC adapter into the mains supply and plug the adapter output plug into the power input socket on the right hand side of the instrument. The orange charging light next to the power input socket will flash to indicate charging and will turn on constantly to indicate full charge.

The batteries will take approximately 4 hours to become fully charged.

Note: Use only the AC adapter supplied. Use of any other type may cause permanent damage to the MicroLab and cause a fire or electrical hazard. Do not plug in and remove the power lead from the AC adapter repeatedly.

Please Note: Dispose of the waste battery in accordance with EU Waste Battery Regulations.

PC connection using SPCS

SPCS is an easy to use PC based windows application that interfaces to the MicroLab via a USB or serial port. It incorporates a database into which patient details can be entered and downloaded to the MicroLab or test results may be uploaded from the MicroLab to the PC.

Using SPCS and the MicroLab, live blows can be performed with the PC directly controlling the operation of the MicroLab.

The results and graphs produced are displayed directly on the PC screen.

The spirometer is connected from the serial or USB port on the PC, to the port on the right hand side of the instrument using the USB or serial cable provided with SPCS.

Note: The MicroLab should only be connected to a computer that is manufactured in accordance with EN60950 1992/1993 – ‘Safety of Information Technology Equipment including Electrical Business Equipment’

Note: Keep the PC and Monitor out of reach of the patient at all times.

It is recommended that whilst the unit is connected to a computer the mains adapter is used.

Looking after your Spirometer

Please observe the following precautions:

- Do not touch the screen with fingers. Use only the stylus provided.
- Use only a damp, lint free, cloth to clean the screen.
- Do not keep the spirometer in a damp place or expose it to extremes of temperature.
- Do not direct the transducer holder towards a strong light source whilst operating the spirometer.
- Check the AC charger for compatibility with local power rating.

Cleaning Instructions

Disinfection of contaminated parts is only effective after having them carefully pre-cleaned. CareFusion recommend the tested solution of PeraSafe sterilizing powder (36-SSC5000A) for pre-cleaning and disinfection. If a different solution is used, please follow the given manufacturer's instructions.

The device must not be wiped with any aqueous solutions and must not be exposed to solvents e.g. alcohol, chloride solutions as there are electronic components inside that will be permanently damaged.

CAUTION: Switch off the device and always unplug your MicroLab before cleaning.

External Surfaces of the Spirometer

CAUTION: Do not attempt to wash or immerse the MicroLab transducer housing in water or cleaning fluid, as there are electronic components inside that will be permanently damaged

The external housing of the spirometer may be wiped with sterile wipes or a damp cloth that has been immersed in a cold sterilizing solution, when required.

CAUTION: Do not wipe the touch screen.

Cleaning accessories

With the use of a filter MicroGard® IIB (V-892381), MicroGard® IIC (V-892384) or a One-way Safety mouthpiece (36-SST1250) for each patient, cleaning of the transducer is recommended once a month.

When using the disposable cardboard mouthpiece (adult: 36-PSA1000, pediatric: 36-PSA1200) without a filter and under the prerequisite that the patient was instructed only to exhale into the transducer, the following parts have to be cleaned once a day: adult/pediatric adapter, transducer.

With any other use as described above all contaminated parts must be disinfected between patients.


Important note: Used single patient nose clip pads, mouthpieces, MicroGard® IIB filter and MicroGard® IIC filter must be disposed off immediately after the use.

If there are changes on the material surfaces (cracks, brittleness) the respective parts must be disposed off.

Cleaning the Transducer

The transducer requires no routine maintenance or servicing. However, if you wish to sterilize or clean the transducer it may be removed by means of the following procedure:

1. Rotating the turbine transducer anti-clockwise until the locating pipe lines up with the small rectangular cut-out in the housing as shown below.
2. Gently pull the transducer away from the housing.
3. The transducer may now be immersed in warm soapy water for routine cleaning or immersed in cold disinfecting solutions e.g. Perasafe (Cat No. 36-SSC5000A) for a maximum of 10 minutes (Alcohol and chloride solutions should be avoided).
4. After cleaning/disinfecting, the transducer should be immersed briefly in distilled water and left to dry.
5. Re-assemble the mouthpiece holder.

CAUTION: Do not attempt to wash or immerse the MicroLab **transducer housing** in water or cleaning fluid, as there are electronic components inside that will be permanently damaged.

Servicing

There is no routine maintenance required for the MicroLab and there are no user serviceable parts in this spirometer. Please return the unit to CareFusion or an authorized agent, if servicing is required.

If your spirometer requires service or repair please see page 96 for contact details.

Trouble Shooting Information

Should you encounter problems operating your MicroLab meter, please consult the table below:

Problem	Possible Cause	Solution
display 'freezes' and the unit does not respond to any key presses	Multiple icons have been selected or accidentally pressed	hold the on/off button down for 10 seconds until the unit switches off and then turn on again
Does not register a blow	Head assembly or cable broken	Replacement of head assembly or return the unit for servicing
No display present	Charger not connected or battery is exhausted	Connect charger to the mains and leave unit to fully charge, or return the unit for servicing
Blows are inverted on the display	Head assembly or cable broken	Replacement of head assembly or return the unit for servicing
Blows tracing ends abruptly although patient still exhaling	Turbine sticking	Clean turbine in warm soapy water or sterilizing solution, if problem persists a replacement turbine may be required
Battery does not hold charge	Exhausted battery	Return the unit for servicing
	Mains charger faulty	Replace the mains charger

Paper not printing	Check paper is housed in printer compartment correctly	Refer to section “paper loading” of this manual, and follow the instructions.
	Incorrect thermal paper being used	Ensure you are using CareFusion recommended thermal paper (See Consumables / Accessories section)
Stylus does not register icons on the display	Touch screen display requires calibration	Select the calibration check icon and choose touch screen, and follow the instructions.
Icons missing from the display	Icon has been de-selected	Hold stylus on to the blue area of the display, a list will appear, ensure required icon is selected.
Calibration procedure failed or cannot be completed	Turbine may be faulty	Repeat calibration procedure, if problem still persists, replace turbine or return unit for servicing
	Turbine not fitted tightly to calibration syringe	Ensure the syringe is fitted to the turbine using the adapter supplied part # 861427. (supplied with the syringe)
	Calibration syringe does not have an inspiratory seal or seal is leaking.	Ensure you are using a CareFusion syringe.
	Shaft of the syringe being pushed down	The syringe should be emptied and filled with one smooth stroke, avoid pushing down on the shaft or banging at the end of each maneuver.

Safety Designation per IEC 60601-1

Type of protection against electrical shock	Internally powered Equipment and Class I
Degree of protection against electrical shock	Type B applied part
Power Equipment	Battery type: NiMH battery pack, 8.4V, 1100mAh
Degree of Electrical connection between equipment and Patient	Equipment designed as non-electrical connection to the patient.
Degree of mobility	Transportable
Mode of operation	Continuous

Classifications according to IEC 60601-1

MicroLab	Applied part, type B
Volume Transducer	Applied part, type B
Option SPO2 Finger Probe	Applied part, type BF
(Nonin IPOD® Integrated Oximetry Device)	

WARNING: No modification of this equipment is allowed.

Note: When you connect other equipment to the unit, always make sure that the whole combination complies with the international safety standard IEC 60601-1 for medical electrical systems. During measurements, connect the MicroLab only to printers and computers that comply with IEC/EN 60601-1/UL 60601-1.

WARNING: The user must not touch any voltage-carrying parts and the patient at the same time.

During database upload the MicroLab may be connected to a computer that complies with EN 60950 – ‘Information technology equipment - Safety - Part 1: General requirements’.

IMPORTANT: Only use the mains adapter supplied (36-PSU1012 12V DC 2.5A) The adapter contains a transformer. Do not cut off the adapter to replace it with another plug as this causes a hazardous situation.

- The adapter transforms the mains voltage (100-240 Volts) to a safe voltage (12V DC)
- Make sure the adapter does not get wet
- Do not use a damaged adapter
- Always unplug your MicroLab before cleaning

WARNING: Do not connect devices that are not specified as part of the system.

Note: If an MPSO (Multiple Portable Socket Outlet) is used with the system, the maximum permitted load should not be exceeded. Do not connect electrical equipment that has not been supplied as part of the system.

Electromagnetic Compatibility (EMC) to EN60601-1:2007

WARNING: use of portable phones or other radio frequency (RF) emitting equipment near the system may cause unexpected or adverse operation

The MicroLab has been tested to EN60601-1-2:2007, regarding its ability to operate in an environment containing other electrical/electronic equipment (including other medical devices).

The purpose of this testing is to ensure that the MicroLab is not likely to adversely affect the normal operation of other such equipment and that other such equipment is not likely to adversely affect the normal operation of the MicroLab.

Despite the testing of the MicroLab that has been undertaken, normal operation of the MicroLab can be affected by other electrical/electronic equipment and portable and mobile RF communications equipment.

As the MicroLab is medical equipment, special precautions are needed regarding EMC (electromagnetic compatibility).

It is important that the MicroLab is configured and installed/put into service, in accordance with the instructions/guidance provided herein and is used only in the configuration as supplied.

Changes or modifications to the MicroLab may result in increased emissions or decreased immunity of the MicroLab in relation to EMC performance.

The MicroLab should be used only with the accessories (USB cables, mains adapter and turbine transducer) supplied (which are referenced in the accessories section of this manual). None of the MicroLab cables should be extended in length by the user.

If any cables are extended by the user or non approved accessories are used, this may result in an increased level of emissions or decreased level of immunity, in relation to the MicroLab's EMC. None of the MicroLab's accessories should be used with other devices, as this may result in an increased level of emissions or decreased level of immunity, in relation to the other devices' EMC.

The MicroLab has an essential performance – the product should continue to operate correctly. In the unlikely event of a Fast Transient / ESD event occurring, the device should be reset and located away from the source of interference.

WARNING: The MicroLab should not be used adjacent to or stacked with other equipment. If adjacent or stacked use with other equipment is necessary, the MicroLab and the other equipment should be observed / monitored, to verify normal operation in the configuration in which it will be used.

Guidance and Manufacturer's Declaration – Electromagnetic Emissions The MicroLab is intended for use in the electromagnetic environment specified below. The customer or the user of the MicroLab should assure that it is used in such an environment.		
Emission Test	Compliance	Electromagnetic Environment - Guidance
RF emissions CISPR 11	Group 1	The MicroLab uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Group B	The MicroLab is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes
Harmonic emissions IEC61000-3-2	Class A	
Voltage fluctuations / flicker emissions IEC61000-3-3	Complies	

Guidance and Manufacturer's Declaration – Electromagnetic Immunity The MicroLab is intended for use in the electromagnetic environment specified below. The customer or the user of the MicroLab should assure that it is used in such an environment.			
Immunity Test	IEC 60601 test level	Compliance Level	Electromagnetic Environment - Guidance
Electrostatic discharge (ESD) IEC61000-4-2	± 6 kV contact ± 8 kV air	± 6 kV contact ± 8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient / burst IEC61000-4-4	± 2 kV for power supply lines ± 1 kV for input / output lines	± 2 kV for power supply lines Input/output line tests not applicable (<3m)	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC61000-4-5	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	Mains power quality should be that of a typical commercial or hospital environment.

Voltage dips, short interruptions and voltage variations on power supply input lines IEC61000-4-11	$<5\% U_T$ $(>95\% \text{ dip in } U_T)$ For 0.5 cycle $40\% U_T$ $(60\% \text{ dip in } U_T)$ for 5 cycles $70\% U_T$ $(30\% \text{ dip in } U_T)$ for 25 cycles $<5\% U_T$ $(>95\% \text{ dip in } U_T)$ For 5 s	$<5\% U_T$ $(>95\% \text{ dip in } U_T)$ For 0.5 cycle $40\% U_T$ $(60\% \text{ dip in } U_T)$ for 5 cycles $70\% U_T$ $(30\% \text{ dip in } U_T)$ for 25 cycles $<5\% U_T$ $(>95\% \text{ dip in } U_T)$ For 5 s	Mains power quality should be that of a typical commercial or hospital environment. If the user of the MicroLab requires continued operation during power mains interruptions, it is recommended that the MicroLab be powered from an uninterruptable power supply or a battery
Power frequency (50/60Hz) Magnetic field IEC61000-4-8	3 A/m	3 A/m	If incorrect operation occurs, it may be necessary to position the MicroLab further from sources of power frequency magnetic fields or to install magnetic shielding. The power frequency magnetic field should be measured in the intended installation location to assure that it is sufficiently low.
NOTE U_T is the a.c. mains voltage prior to application of the test level.			

Guidance and Manufacturer's Declaration – Electromagnetic Immunity The MicroLab is intended for use in the electromagnetic environment specified below. The customer or the user of the MicroLab should assure that it is used in such an environment.			
Immunity Test	IEC 60601 test level	Compliance Level	Electromagnetic Environment - Guidance
Conducted RF IEC61000-4-6	3 Vrms 150 kHz to 80 MHz	3 Vrms	Portable and mobile RF communications equipment should be used no closer to any part of the MicroLab, including any cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance (d) $d = 1.2\sqrt{P}$
Radiated RF IEC61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m	$d = 1.2\sqrt{P}$ 80 MHz to 800 MHz $d = 2.3\sqrt{P}$ 800 MHz to 2.5 GHz Where <i>P</i> is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Fields strengths from fixed RF transmitters, as determined by an electromagnetic site survey,^a should be less than the compliance level in each frequency range.^b Interference may occur in the vicinity of equipment marked with the following symbol: 
NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies. NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			

- ^a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the MicroLab is used exceeds the applicable RF compliance level above, the MicroLab should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orientating or relocating the MicroLab.
- ^b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

Recommended separation distances between portable and mobile RF communications equipment and the MicroLab

The MicroLab is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the MicroLab can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the MicroLab as recommended below, according to the maximum output power of the communications equipment.

Rated Maximum Output Power of Transmitter in Watts (W)	Separation Distance in Meters (m) according to Frequency of Transmitter		
	150 KHz to 80 MHz	80 MHz to 800 MHz	800 MHz to 2.5 GHz
	$d = 1.2\sqrt{P}$	$d = 1.2\sqrt{P}$	$d = 2.3\sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23.3

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures objects and people.

Symbols

Type B device

0086

In accordance with Directive 93/42/EEC

Disposal in compliance with WEEE

Consult the instructions for use

Caution: Consult the accompanying documents

Manufacturer

Date of Manufacture

Serial number

Direct current

Single patient use

Federal U.S. law restricts this device to sale by or on the order of a physician. (Rx only)

Nationally Recognized Test Laboratory (NRTL) officially recognized by the American Occupational Safety and Health Administration (OSHA) for OSHA Electrical safety and compliance.

Specification of the MicroLab

General:

Storage:	>2000 tests including Flow/Volume loops and Volume/Time curves
Printer Output:	Compatible Hewlett Packard USB printers (Specified models only).
Display:	Color 1/4VGA LCD.
Power supply:	Input 100 to 240V, 47 to 63Hz. Output 12V 2.5A (Class 1) Type: MENB1030A1200F03
Battery Pack:	Rechargeable NiMH 8.4V 1A-hours
Battery life:	Approximately 30 hours with a fully charged new battery
Dimensions:	35x255x120 mm. Transducer 50x60x90mm.
Weight:	630g.
Operating Temperature:	0 to +40 °C
Operating Humidity:	30% to 90% RH

Transport and Storage -20 to + 70 °C

Temperature:

Transport and Storage 10% to 90% RH

Humidity:

Spirometry Measurements

Relaxed Expiratory Vital Capacity (VC)

Forced Expired Volume in 0.75 seconds (FEV_{0.75})

Forced Expired Volume in 1 second (FEV₁)

Forced Expired Volume in 3 second (FEV₃)

Forced Expired Volume in 6 seconds (FEV₆)

Forced Vital Capacity (FVC)

Peak Expiratory Flow Rate (PEF)

FEV_{0.75} as a percentage of VC (FEV_{0.75}/VC)

FEV_{0.75} as a percentage of FVC (FEV_{0.75}/FVC)

FEV₁ as a percentage of VC (FEV₁/VC)

FEV₁ as a percentage of FVC (FEV₁/FVC)

FEV₃ as a percentage of VC (FEV₃/VC)

FEV₃ as a percentage of FVC (FEV₃/FVC)

FEV_{0.75} as a percentage of FEV₆ (FEV_{0.75}/FEV₆)

FEV₁ as a percentage of FEV₆ (FEV₁/FEV₆)

Maximum Expired Flow at 75% of FVC remaining (MEF75)

Maximum Expired Flow at 50% of FVC remaining (MEF50)

Maximum Expired Flow at 25% of FVC remaining (MEF25)

Mean Mid-Expiratory Flow Rate (MMEF)

Forced expiratory flow at 50% of volume as a percentage of VC (FEF50/VC)

Forced expiratory flow at 50% of volume as a percentage of FVC (FEF50/FVC)

Maximal voluntary ventilation indicated (MVV (ind))

Forced inspired volume in 1 second (FIV₁)

Forced inspiratory Vital Capacity (FIVC)

Peak Inspiratory Flow Rate (PIF)

FIV₁ as a percentage of FIVC (FIV₁/FIVC)

Forced inspiratory flow at 25% of inhaled volume (FIF25)

Forced inspiratory flow at 50% of inhaled volume (FIF50)	
Forced inspiratory flow at 75% of inhaled volume (FIF75)	
Forced expiratory flow at 50% of volume as a percentage of FIF50 (FEF50/FIF50)	
The time taken between 25% and 75% of the forced expired volume (MET2575)	
Forced Expiratory Time (FET)	
Tidal Volume (TV)	
Expiratory reserve volume (ERV)	
Inspiratory reserve volume (IRV)	
Inspiratory capacity (IC)	
Expiratory Relaxed vital capacity (EVC)	
Inspiratory vital capacity (IVC)	
Breathing frequency rate (FR)	
Inspiratory time (Ti)	
Expiratory time (Te)	
Ti as a % of total breath timely (Ti/Ttot)	
Tidal volume as a % of Ti (TV/Ti)	
Breath Rate	BR
Breathing Time	B.T
Volume Tidal	VT
Expiratory Time - average time of expiration per breath in seconds	Te
Inspiratory Time - average time of inspiration per breath in seconds	Ti
Total Tidal Breath Time in Seconds	TTOT = Ti + Te
Ratio of Average Expiratory and Inspiratory Breaths	Ti/Te
Average Time of Expiration per Breath as a ratio to the Total Tidal Breath Time	Ti/TTOT

Tests per subject: 5 VC maneuver

Predicted Values:

8 FVC maneuvers

Various – depends upon national preference

Transducer:

CareFusion Bi-Directional Digital Volume.

Resolution:

10ml volume 0.03l/s flow

Accuracy:

+/-3%. To ATS recommendations – Standardization of spirometry 1994 update for flows and volumes.

Consumables / Supporting Products

Description

Thermal Printer Paper (pack of 5 rolls)

Archive Printer Paper ~ 25 Year Anti-Fade (pack of 5 rolls)

Adult Disposable Mouthpieces (500 per box)

One-way Safety Mouthpieces (500 per box)

One-way Safety Mouthpieces (250 per box)

Pediatric Disposable Mouthpieces (250 per box)

Paediatric Adaptor

MicroGard Pulmonary filters

MicroGard Pulmonary filters

PeraSafe Sterilizing Powder 81g (to make up 5 liters of solution)

Nose Clips (pack of 5)

3 Liter Calibration Syringe

Silicone adapter "oval" for MicroGard® II B

Calibration syringe adapter

SPCS Software

For further information or to place an order for Disposables/
Supporting products please contact CAREFUSION, your local
distributor or view our website www.carefusion.com/micromedical

NOTE: USE ONLY CAREFUSION ACCESSORIES
--

Note: Information in this manual is subject to change without notice and does not represent commitment on the part of CareFusion 232 UK Ltd. The software may be used or copied only in accordance with the terms of that agreement. No part of this manual may be reproduced or transmitted in any form or by any means, electronic or mechanical, including photocopying and recording for any purpose without the written permission of CareFusion 232 UK Ltd.

Customer contact information

UK Customers only

For all Sales Order processing for products, training and spare parts, Service and Technical Support enquiries, please contact the following:

CareFusion UK 232 Ltd
UK Customer Service & Support
The Crescent
Jays Close
Basingstoke
RG22 4BS

Customer Service Sales Enquiries:

Telephone: 01256 388550

Email: micro.uksales@carefusion.com

Factory Repair and Administration Enquiries:

Telephone: 01256 388552

Email: micro.ukservice@carefusion.com

Technical Support Enquiries:

Telephone: 01256 388551

Email: support.rt.eu@carefusion.com

International customers only

For all Sales Order processing for products and Spare parts, Service and Technical Support enquiries, please contact the following:

CareFusion Germany 234 GmbH
Customer Service & Support International
Leibnizstrasse 7
97204 Hoechberg
Germany

Customer Service Sales Enquiries:

Telephone: 0049 931 4972 670

Email: micro.internationalsales@carefusion.com

Factory Repair and Administration Enquiries:

Telephone: 0049 931 4972 867

Email: support.admin.eu@carefusion.com

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© CareFusion 2012
CareFusion Duitsland 234 GmbH