

This guide is a reference for using the SARS-CoV-2 Antigen Self Test Nasal It is essential that you read the Instructions for Use for patients before using this test.

IVD CE

1.

2.

۲

1 Preparing for a test

1. Carefully read the Instructions for Use for patients for the SARS-CoV-2 Antigen Self Test Nasal.

2 Collecting and preparing nasal sample

Open the foil pouch **2** by tearing along the tear-

and one nozzle cap and place them on the table

spilling the liquid inside the tube. Place the tube

line and take out one of the tubes with liquid

Open the seal of the tube carefully without

2. Wash your hands with soap and water or use a hand sanitizer before performing the test.



in the tube holder.

pouches. Do not use the test if the expiry date has passed

3. Check the expiry date \square on the back of the foil

Image / REF No. 9901-INCOV-06G Image / LOT No. // LOT No.	SARS-CoV	-2 Antigen Sel
	REF No.	9901-NCOV-06G

Blow your nose once using a tissue. 3.

Important

Warning!

test

safety information

· Wash your hands with soap and water or use a hand sanitizer before and after performing the

· Keep the swab clean. Avoid touching the tip of the swab and ensure that it does not come into

contact with any surfaces prior to use.

· Avoid contact with skin and eyes. In case of

thoroughly using a suitable disinfectant.

Do not swallow any of the components.

accidental contact, rinse well in order to avoid skin


irritations. In case of concerns, consult your doctor.

In the event of a spillage, ensure that it is cleaned


Do not use any of the test components in the body with the exception of the swab included in the kit.

Use the test kit only once.

- Remove the swab from the packaging. Ensure 4. that you only touch the handle of the swab and not the soft pad at the tip.
- Rotate the swab 4 times (for a total of approx. 7. 15 seconds) against the lining of the nasal wall before removing it from the nostril.
- 8. Repeat steps 6 and 7 in your right nostril using the same swab Note: Samples must be collected from both nostrils using the same swab.



- Insert the swab into the tube until the soft pad is in 9. the liquid. Squeeze the tube at the bottom and hold it tight. Stir the swab more than 10 times to transfer the biological material from the swab to the liquid.
- 10. Remove the swab while squeezing the sides of the tube to extract the liquid from the swab. Dispose the swab and seal the tube securely with the nozzle cap.



WARNING! Failure to squeeze the tube can lead to incorrect results due to excess buffer in the swab.

3 Performing a test

- 1. Place the test device on a flat surface.
- 2. Hold the tube upright above the circular well marked on the picture below (not over the rectangular result window).



3. Drop exactly 4 drops onto the circular well. Gently squeeze the sides of the tube together if necessary. Note: You can continue with the test even if you accidentally drop 5 drops onto the test device.

your left nostril. Slowly slide the swab approx.

2 cm forward (parallel to the roof of your mouth - not upwards) until you encounter resistance.

Slightly tilt your head backwards.

Do not apply any pressure.

5. 6.

> 4. Set the timer and read the test result after 15 to 30 minutes

WARNING! Test results that are read before 15 minutes or after 30 minutes may be incorrect.

15-30 min

- 5. Wash your hands with soap and water or use a hand sanitizer after performing the test.

Components of the test kit

Pouch 1 containing a test device and a desiccant package (x5)




Sterile swab

desiccant package.










4. Open one of the foil pouches **1** by tearing along

the tear-line and take out the test device and the




Tube holder





- Timer
- Tissue


5. Ensure that the test device is intact and that there are no green beads in the desiccant package. Do not open the desiccant package.








4 Interpreting results





1. If a control line (C) is not visible, the result must be considered invalid.

The test is not working correctly and you should perform another test using a different test kit.

You may have performed the test incorrectly. Carefully read the Instructions for Use and repeat the test. If your test result is still invalid, please contact your doctor or a COVID-19 test center.


2. If a test line (T) is visible together with a control line (C), this means that the result is positive.

Look carefully at the result: The test should be considered positive if two lines are visible - even if they are faint.

A positive test result means it is very likely that you have COVID-19. Please contact your doctor/primary care physician or your local health authority immediately and adhere to the local auidelines regarding self-isolation. Your doctor may require you to undergo a PCR test to confirm the result.


SD Biosensor, Inc.

C-4th&5th, 16, Deogyeong-daero 1556beon-gil Yeongtong-gu, Suwon-si, Gyeonggi-do, 16690, REPUBLIC OF KOREA Made in Korea



Distributed by: Roche Diagnostics GmbH Sandhofer Strasse 116 D-68305 Mannheim, Germany www.roche.com

۲



3. If a control line (C) is visible (regardless of how faint it is) and a test line (T) is not visible, this means that the result is negative. It is unlikely that you have COVID-19.

However, even if your test is negative, continue to observe all hygiene and safety measures.

If you suspect that you have an infection (i.e., if you have prolonged symptoms or if your symptoms are worsening), contact your doctor/primary care physician. You may have another infection, or your test result may be false. You may repeat the test after 1-2 days, as COVID-19 cannot be detected with complete accuracy during all stages of an infection.

Document version: Initial version

Revision date: April 2021 © 2021. All rights reserved.

SARS-CoV-2-AGST-NASAL_1-0_QRG_09436090001-01_EN_pst-ROW.indd 1

07/05/2021 15:34:44