

# SAM IO® Stabilizer


**SAM**  
MEDICAL




## EN INSTRUCTIONS FOR USE




### APPLICATION STEPS




With catheter in position, consider use of Stabilizer.



Place sterile Stabilizer carefully over catheter and against skin.




Attach primed extension set to catheter.  
**KEY POINT:** Extension set design permits single-handed tightening to catheter.




Optional Stabilizer Application Method:


With attached extension set in position gently feed extension set through Stabilizer cylinder. Guide Stabilizer over catheter, and against patient skin.




Ensure catheter is centered within Stabilizer.



Hold Stabilizer in position and gently pull Tab 1. By pulling Tab 1 you are removing adhesive backing. This will expose adhesive to patient skin. Gently press exposed adhesive against patient skin.



Hold Stabilizer in position and gently pull Tab 2. By pulling Tab 2 you are removing adhesive backing. This will expose adhesive to patient skin. Gently press exposed adhesive against patient skin.




























With Stabilizer in position (and adhesive backing removed) gently smooth all exposed adhesive to patient skin.

Notice:

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the regulatory authority/competent authority of the Member State or region in which the user and/or patient is established.

SYMBOL GLOSSARY:

Symbol	Title (Reference)	Description
	Manufacturer (5.1.1 <sup>[1]</sup> )	Indicates the medical device manufacturer.
	Authorized representative in the European Community (5.1.2 <sup>[1]</sup> )	Indicates the authorized representative in the European Community/ European Union.
	Date of Manufacture (5.1.3 <sup>[1]</sup> )	Indicates the date when the medical device was manufactured.
	Use-by date (5.1.4 <sup>[1]</sup> )	Indicates the date after which the medical device is not to be used.
	Batch Code (5.1.5 <sup>[1]</sup> )	Indicates the manufacturer's batch code so that the batch or lot can be identified.
	Catalogue Number (5.1.6 <sup>[1]</sup> )	Indicates the manufacturer's catalog number so that the medical device can be identified.
	Importer (5.1.8 <sup>[1]</sup> )	Indicates the entity importing the medical device into the locale.
	Country of Manufacture (5.1.11 <sup>[1]</sup> )	To identify the country of manufacture of products.
	Sterilized using ethylene oxide (5.2.3 <sup>[1]</sup> )	Indicates a medical device that has been sterilized using ethylene oxide.
	Do not re-sterilize (5.2.6 <sup>[1]</sup> )	Indicates a medical device that is not to be re-sterilized.
	Do not use if package is damaged and consult instructions for use (5.2.8 <sup>[1]</sup> )	Indicates a medical device that should not be used if the package has been damaged or opened and that the user should consult the instructions for use for additional information.
	Single sterile barrier system (5.2.11 <sup>[1]</sup> )	Indicates a single sterile barrier system.
	Single sterile barrier system with protective packaging outside (5.2.14 <sup>[1]</sup> )	Single sterile barrier system with protective packaging outside.
	Keep away from sunlight (5.3.2 <sup>[1]</sup> )	Indicates a medical device that needs protection from light sources.
	Keep dry (5.3.4 <sup>[1]</sup> )	Indicates a medical device that needs to be protected from moisture.
	Do not re-use (5.4.2 <sup>[1]</sup> )	Indicates a medical device that is intended for one single use only.
	Consult instructions for use or consult electronic instructions for use (5.4.3 <sup>[1]</sup> )	Indicates the need for the user to consult the instructions for use.
	Medical Device (5.7.7 <sup>[1]</sup> )	Indicates the item is a medical device.
	Unique Device Identifier (5.7.10 <sup>[1]</sup> )	Indicates a carrier that contains unique device identifier information.
	CE Marking (Regulation (EU) 2017/745 Article 20)	Indicates device is approved for sale in the European Community/European Union.
	Prescription Only (21 CFR 801.109)	Caution: Federal law restricts this device to sale by or on the order of a physician.
	Swiss Authorized Representative (Section 3 <sup>[2]</sup> )	Indicates the authorized representative in Switzerland.
	Package Quantity (No Applicable Reference)	Indicates a quantity of 30 single devices.
	Package Quantity (No Applicable Reference)	Indicates a quantity of 6 five packs
	Not For Sternal Use (No Applicable Reference)	Indicates device is not to be used on the sternum

[1] ISO 15223-1:2021 Medical devices - Symbols to be used with information to be supplied by the manufacturer - Part 1: General requirements  
[2] Swissmedic, Obligations Economic Operators CH, VM-ID: MU600\_00\_016e/V4.0/mea/pmi/06.04.2023

REF IO740-1P

REF IO740-5P



CH REP

MedEnvoy Switzerland  
Gotthardstrasse 28  
6302 Zug  
Switzerland



2862

EC REP

CS Life Sciences Europe Ltd  
The Black Church, St. Mary's Place  
Dublin 7, Dublin  
D07P4AX, Ireland  
eurep@cslifesciences.com



ZIEN Medical Technologies, Inc.  
2490 South 300 West  
Salt Lake City, UT 84115 USA



MedEnvoy  
Prinses Margrietplantsoen 33 - Suite 123  
2595 AM The Hague  
The Netherlands



Customer Service:

Tel: +1.503.639.5474 (USA)

technical.support@sammedical.com

sammedical.com

© 2025 all rights reserved. SAM® Medical and SAM IO® Intraosseous Access System® are trademarks of SAM® Medical.

SAM04109-EN REV A  
ISSUE DATE: 2025-05-08