













Digital thermometer Model Eco Temp Smart

- Digitaltas Fieberthermometer
- Thermomètre digital
- Termometro digitale
- Termómetro digital
- Digitalt termometer
- Digitaalinen kuumemittari
- Digitalt termometer
- Digitale koortsthermometer
- Ψηφιακό θερμομέτρο
- Digitalt klinik termometresi
- Цифровой термометр

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Digital Thermometer English

- A. Sensor Tip
- B. Display Window
- C. I/O Button
- D. Storage Case

Thank you for purchasing the OMRON Digital Thermometer Eco Temp Smart. The thermometer you purchased offers a safe, accurate and quick temperature reading. You can measure your temperature either in the anus (rectally), in the mouth (orally) or in the armpit (axillary). This thermometer eliminates any worries about broken glass or mercury hazards. Note that your temperature is affected by smoking, eating or drinking. Contact your doctor if you have any questions regarding specific temperatures.

Important Safety Notes

To assure the correct use of the product basic safety measures should always be followed including the precautions listed below.

Warning

- A high or prolonged fever requires medical attention, especially for small children. Please contact your doctor.
- Do not attempt to take temperatures at other sites, such as in the ear, as it may result in false readings and may lead to injury.
- Store the thermometer out of the reach of children. Do not allow children to take their temperatures unattended. Children may injure themselves when attempting to take temperatures without supervision.
- Do not leave the battery, battery cover or probe cover where children can get to them. Children may swallow them. Should a child swallow the battery, battery cover or probe cover, contact a doctor immediately.
- Do not attempt measurements when the thermometer is wet as inaccurate readings may result.

3. Displaying the last temperature

Press the I/O button to display the temperature that was measured most recently. For details, please refer to the "2. Use of the thermometer" section.

4. Correct measurement

Measurement accuracy cannot be assured when the method used to measure the temperature is incorrect.

Rectal use:
 • Commonly used for small children when it is difficult to take an oral or armpit temperature.
 • Lubricate with non-irritating skin-friendly gel.
 • Gently insert the probe tip no more than 1.3 cm into the rectum. Do not force the tip into the rectum if resistance is encountered.
 • Disinfect the thermometer after use.

Oral use (for ages 4 and older):
 The mouth should remain closed for 5 minutes before starting a reading.

• Place the thermometer in the mouth under the tongue so that it rests to the left or right of the root of the tongue.
 • Use downward tongue pressure to hold the thermometer in place.
 • Hold the thermometer to keep it from sliding around in the mouth. (Fig. 9)

Armpit (axillary) use:
 • Place the sensor tip in the centre of the armpit. (Fig. 10)
 • Lock the sensor tip under the arm, using the arm to slightly apply pressure inwards.

The following may result in incorrect measurement:
 • Taking the temperature immediately after exercising, bathing, eating or drinking will result in an incorrect measurement.
 • Do not use the thermometer in places where strong static electricity or electromagnetic fields are present. Doing so may lead to inaccurate readings and may contribute to instrument failure.
 • Do not attempt rectal measurements on persons with rectal disorders. Doing so may aggravate or worsen the disorder.
 • Do not step on the unit or the hard case.

General safety precautions
 • Do not use the thermometer to measure anything other than human body temperatures.
 • Do not drop the thermometer and storage case or expose them to shock or vibration.

1. Switching between Fahrenheit/Centigrade (Celsius)
 Press the I/O button for 7 seconds until a beep is emitted. The display will change to the other measurement system at the time when the beep sounds. (Fig. 1, 2, 3)

2. Use of the thermometer
 1. Press the I/O button. (Fig. 4)
 The result of the previous measurement is displayed for a few seconds. Then 37.0°C (98.6°F) is displayed as an internal test-function. (Fig. 5)
Note: When a temperature other than 37.0°C (98.6°F) is displayed, please refer to the "6. Display Indications and Troubleshooting" section for details to correct the problem.

2. Apply the thermometer to the measurement site.
Note: Even if the "L" indicator fails to display (when only the temperature is displayed), this is not an error. The temperature may be taken as is. (Fig. 6)

3. Buzzer notification timing
Rectal: approx. 10 seconds
Oral: approx. 20 seconds
Armpit: approx. 25 seconds*
 * 25 seconds is the earliest estimated temperature indication. To complete the measurement, it is required to measure the temperature in the armpit between 2-5 minutes (depending on the ambient temperature and condition).

Notes:
 • The measurement continues even after the buzzer notification.
 • The buzzer notification timing is based on an environmental temperature of 23 degrees.
 • Since oral and armpit may give lower measurement results it is recommended to apply the device rectally to get the most accurate result.
 • Turn off the thermometer. Clean and disinfect before returning it to its protective case.
Notes:
 • The reading will automatically be stored in memory.
 • Turn off the thermometer manually to conserve the battery. It will automatically switch off 30 minutes after use or 3 minutes when turned on but not used.

5. Battery replacement
Battery: Alkaline-Manganese Button Battery LR41 (available commercially)
Battery Replacement Indicator ()
 Replace the battery when the battery replacement indicator appears when the thermometer is switched on.

1. Use a coin to take off the cover. (Fig. 12)
2. Remove the battery. (Fig. 13)
3. Insert the new battery with the "+" pole facing up as shown in the diagram. Make sure the battery is put in with the right polarity. (Fig. 14)
4. Use a coin to replace the battery cover. (Fig. 15)
Notes:
 • The thermometer comes fitted with a battery for testing purposes. This battery may not provide the same amount of battery life as a new battery.
 • Do not loose the O-ring around the screw. Water-resistant will fail and can lead to misreadings and failure of the unit.

6. Display Indications and Troubleshooting

Indicator	Cause	Remedy
L (low)	Displayed when the temperature of the sensing section is less than 32°C (89.6°F).	This is not an error so it may be ignored and the temperature taken.
H (high)	Displayed when the temperature of the sensing section is higher than 42°C (107.6°F).	Cool the sensing section and then take the temperature.
Air Temperature Displayed	Displayed when the surrounding temperature exceeds 32°C (89.6°F).	Turn the thermometer on and off several times. If a temperature other than that in the range from 36.9°C (98.4°F) to 37.1°C (98.8°F) is displayed as a test display, please contact your nearest authorized OMRON dealer.

7. Storage
Store the thermometer in its storage case.
 • Do not store the thermometer in the following types of places. Doing so may damage the thermometer:
 - Wet locations.
 - Locations with high heat and humidity or those that are exposed to direct sunlight. Areas close to heating equipment, dusty locations, or environments where there are high salt concentrations in the air.
 - Locations where the unit will be subjected to leaning, shock or vibration.
 - Pharmaceutical storage areas or locations where corrosive gases are present.

8. Cleaning and Disinfecting
 OMRON recommends to wipe the device with damp cloth. For disinfection, 70% Ethanol or Isopropyl alcohol can be used. Do not attempt to disinfect the sensing section of the thermometer by immersing in alcohol or in hot water (water over 50°C). Do not submerge deeper than 15 cm and longer than 30 minutes.

9. Specifications
Product Type and Number Digital Thermometer MC-341-E/MC-341-EB
Product Name Eco Temp Smart
Power Consumption: 0.1 mW
Battery Life: Approx. 2 years or more (3 times per day)
Sensing Unit: Thermistor
Measurement Method: Actual measurement (non-predictive)
Temperature Display: 3 digits (°C), or 4 digits (°F) display in 0.1 degree increments
Measurement Accuracy: ±0.1°C (32.0 to 42.0°C) ±0.2°F (89.6 to 107.6°F) (when measured at a standard room temperature of 23°C (73.4°F) in a test tank held at constant temperature) 32.0 to 42.0°C (89.6 to 107.6°F)

Usage Environment Temp and Humidity: Surrounding temperature +10 to +40°C (+50 to +104°F), Relative Humidity 30 ~ 85%RH
Storage Environment Temp and Humidity: Surrounding temperature -20 to +60°C (-4 to +140°F), Relative Humidity 10 ~ 95%RH
External Dimensions: Approx. 112 × 115 mm (with battery installed)
Package Content: Test Battery (Alkaline-Manganese Button Battery LR41), Storage Case, Instruction Sheet.

Measurement Range: 32.0 to 42.0°C (89.6 to 107.6°F)
Use Environment Temp and Humidity: Surrounding temperature +10 to +40°C (+50 to +104°F), Relative Humidity 30 ~ 85%RH
Storage Environment Temp and Humidity: Surrounding temperature -20 to +60°C (-4 to +140°F), Relative Humidity 10 ~ 95%RH
External Dimensions: 19.4 mm (w) × 132.5 mm (l) × 10.0 mm (d)
Package Content: Test Battery (Alkaline-Manganese Button Battery LR41), Storage Case, Instruction Sheet.

Applied part: Type B

Notes:
 • The specifications may be changed without prior notice.
 • This unit is water-resistant.
 • The top of the thermometer contains Nickel.
 • This thermometer does not contain any serviceable part. OMRON recommends that the accuracy is verified every two years by OMRON's representative.
 • This OMRON product is produced under the strict quality system of OMRON Healthcare Co. Ltd., Japan.

Warranty
 This product is guaranteed by OMRON for 3 years after date of purchase. Any damage caused by improper handling shall not be covered by warranty.
 Batteries and packaging are also excluded from the warranty. All other damage claims excluded. A warranty claim must be submitted with the purchase receipt.

Read the Instructions in this manual carefully.

CE 0197

This device fulfills the provisions of the EC directive 93/42/EEC (Medical Device Directive) and the European Standard EN12470:2000, Clinical thermometers - Part 3: Performance of compact electrical thermometers (non-predictive and predictive) with maximum device.