
# **1-PAD** Instructions for Use

# Rescue Steps

- 1. Connect the defibrillator pads to the i-PAD and then place on patient.
- 2. Automatic shock delivery if needed.
- 3. Administer CPR



Quick Reference Guide

# **i-PAD** Instructions for Use

The information in these Instructions for Use applies to the i-PAD NF1201. This information is subject to change.

Please contact CU Medical Systems, Inc. or its authorized representatives for information on revisions.

Printed in the Republic of Korea

Publication Date: March 2013 Part No.: OMP-NF2-E-02

Published by CU Medical Systems, Inc.

#### Copyright

© 2013 CU Medical Systems, Inc.

No part of these Instructions for Use may be reproduced without the permission of CU Medical Systems, Inc.

#### Medical Device Directive

The i-PAD NF1201 complies with the requirements of the Medical Device Directive 93/42/EEC and its revisions.



0470

The i-PAD NF1201 is manufactured by:

CU Medical Systems, Inc.

130-1, Dongwhagongdan-ro, Munmak-eup, Wonju-si,

Gangwon-do, Korea

EU Authorized Representative of CU Medical Systems, Inc.

Medical Device Safety Service

Schiffgraben 41, 30175 Hannover, Germany

#### Product and Order Inquiries:

Overseas Sales Team

CU Medical Systems, Inc.

5F Cheoncaye Bldg., 221, Anyangpangyo-ro,

Uiwang-si, Gyeonggi-do, Korea

Tel: +82 31 421 9700

Fax: +82 31 421 9911

E-mail addresses : sales@cu911.com

Service Request and Technical Support

Customer Service Team

CU Medical Systems, Inc.

130-1, Dongwhagongdan-ro, Munmak-eup, Wonju-si, Gangwon-do, Korea

Tel: +82 33 747 7657

Fax: +82 33 747 7659

E-mail addresses: service@cu911.com

Branch Office: CU Medical Germany GmbH

Ernst-Augustin-Str.5-11, 12489, Berlin, Germany

Tel: +49 30 6781 7804 Fax: +49 30 6782 0901 CU Medical Systems, Inc



CU Medical Systems, Inc.

# CU Medical Systems, Inc.

# Table of Contents

Ins	Instructions for Use Conventions · · · · · · · · · · · · · · · · · · ·	•	•	•	•	•	•	8
Ov	Overview · · · · · · · · · · · · · · · · · · ·	•		•		•		9
1.								10
	Product Description · · · · · · · · · · · · · · · · · · ·						•	10
	Indications for Use						•	10
	Required User Training • • • • • • • • • • • • • • • • • • •		\	•	•	•	•	10
2.	2. Device Orientation · · · · · · · · · · · · · · · · · · ·			•		•	•	12
3.	3. Setting up the i-PAD · · · · · · · · · · · · · · · · · · ·							14
	Package Contents · · · · · · · · · · · · · · · · · · ·						•	14
	Setting up the i-PAD · · · · · · · · · · · · · · · · · · ·	•		•	•	•	•	15
4.	4. Using the i-PAD · · · · · · · · · · · · · · · · · · ·							17
	Preparation for Defibrillation · · · · · · · · · · · · · · · · · · ·	•		•	•	•	•	18
5.	5. After Using the i-PAD · · · · · · · · · · · · · · · · · · ·							
	After Each Use							
	Data Storage and Transfer · · · · · · · · · · · · · · · · · · ·							
	Last Use Data	•		•	•	•	•	27
	Data Transmission • • • • • • • • • • • • • • • • • • •	•		•	•	•	•	28
	Device Configuration · · · · · · · · · · · · · · · · · · ·	•		•	•	•	•	29
6.	J 1							30
	General Operating Guidelines · · · · · · · · · ·							30
	Poutine Maintenance · · · · · · · · · · · · · · · · · · ·							
	Cleaning the i-PAD · · · · · · · · · · · · · · · · · · ·	•		•	•	•	•	35

# CU Medical Systems, Inc.

# Table of Contents

/.	Troubleshooting	
	Automated Tests · · · · · · · · · · · · · · · · · ·	
	Device Status · · · · · · · · · · · · · · · · · · ·	
	Troubleshooting the i-PAD $\cdots \cdots 39$	
	Troubleshooting During Emergencies · · · · · · · · · · · · 39	
	Troubleshooting Outside of Emergencies 41	
8.	Servicing the i-PAD · · · · · · · · · · · · · · · · · · ·	
	Warranty • • • • • • • • • • • • • • • • • • •	
	Warranty Disclaimer • • • • • • • • • • • • • • • • • • •	
	Service • • • • • • • • • • • • • • • • • • •	
	How to Contact Us · · · · · · · · · · · · · · · · · ·	
	APPENDICES	
Δ	Accessories · · · · · · · · · · · · · · · · · · ·	
Λ. /	A.1 Standard Access ories	
	A.2 Optional Accessories	
	A.3 Service Center · · · · · · · · · · · · · · · · · · ·	
	A.5 Service Center 11 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1	
D	Equipment Symbols · · · · · · · · · · · · · · · · · · ·	
D.	B.1 i-PAD Defibrillator · · · · · · · · · · · · · · · · · · ·	
	B.2 i-PAD Packaging · · · · · · · · · · · · · · · · · · ·	
	B.3 Symbols on Accessories · · · · · · · · · · · · · · · · · · ·	
	B.3.1 Battery [CUSA0601F] · · · · · · · · · · · · · · · · · · ·	
	B.3.2 Pads [CUA0903PF]	
	B.3.2 Pads [CUAU903PF]	
_	Glossary · · · · · · · · · · · · · · · · · · ·	
C.	Glossary • • • • • • • • • • • • • • • • • • •	
_	: DAD NE1201 Considerations	
D.	i-PAD NF1201 Specifications · · · · · · · · · · · · · · · · · · 54	
г	Electromagnetic Compatibility 61	
⊏.	ciectioniagnetic compatibility	

# Instructions for Use Conventions

These Instructions for Use the following conventions:



Conditions, hazards, or unsafe practices that can result in serious personal injury or loss of life.



Conditions, hazards, or unsafe practices that can result in minor or moderate personal injury, damage to the device, or loss of data stored in the device, particularly if precautionary steps are not taken.

#### NOTICE

Used to denote items that are important during installation, operation, or maintenance of the device.

Thank you for purchasing i-PAD. This product can be successfully and safely used for a long period if you familiarize yourself with the instructions and precautions by reading the user's manual prior to its use.

This Product is Fully Automatic External Defibrillator (AED: Automated External Defibrillator) that is used for sudden cardiac arrest patients. It can be easily used by laypersons or those who are not trained how to use the



product.

✓ A defibrillator uses the electric shock of the high voltage and current. Therefore, the user must be well-acquainted with the instructions and precautions by thoroughly reading this user's manual so as to understand safety cautions.

The user must follow the instructions below.

- The user must follow the instructions in the user's manual when using this Product.
- We are not responsible for any issues regarding the Product that are caused by the user's negligent operation or management.
- The repair services for this Product are only provided by us or authorized distributors.
- If the Product is intended to be connected to equipment other than those stated in this user's manual, please contact us prior to the use.
- If this Product does not operate properly, please contact us or a service center authorized by us.

# CU Medical Systems, Inc

# 1. Product Information

#### 1.1 Product Description

This product model is the NF 1201 from the i-PAD product family. This product is an easy-to-use Fully Automatic External Defibrillator (AED) that is a small, light, portable equipment with a battery.

The AED is a defibrillator that delivers electric shock to a patient with arrhythmia in order to recover a normal heart rhythm. Sudden cardiac arrest can occur anytime to anyone at any place and may threaten the patient's life if the appropriate CPR and/or electric shock with a defibrillator are not applied within a

#### 1.2 Indications for use

few minutes.

This Product should be used for sudden cardiac arrest patients (ventricular fibrillation, ventricular tachycardia). A sudden cardiac arrest patient refers to a person with

No response

No normal breathing

#### 1.3 User Credentials


This product can be easily used by laypersons or those who are not trained how to use the Fully-Automated External Defibrillator (AED: Automated External Defibrillator). For those who have little or no understanding of the product use, education and training are recommended for proper use of the product and quick response to an emergency.

#### Note


✓ When delivering shocks to children who are older than 1 year and younger than 8 years as well as lighter than 25 kg, use pediatric pads that are sold separately.



# 2. Device Orientation



Top view of the i-PAD NF1201 with its cover closed



Perspective view of the i-PAD NF1201 with its cover open

**Power Button** 

Turns the equipment on or off.

i-Button

Reports equipment usage (the total hours of the last usage and number of shocks), sets the CPR guide (breath guide, chest compression rate per minute, default guide mode), and checks for errors.

Status Indicator

Indicates the current status of the equipment by lighting or flashing in

different colors.

- Flashing in green: Normal. (In standby mode)

- Green lighting: Normal operation

- Flashing in red: Low battery or equipment error.

- Blue lighting: Running a self-diagnostic test.

- White lighting: Administrator mode.

This is for either information regarding the last use of the equipment or communication with a

PC (personal computer).

- Flashing in yellow: Abnormality in either pads or connection

between pads and equipment.

(In standby mode)

Handle

For carrying the equipment.

Cover

Stores pads.

Defibrillator Pads Insert

Connects the connectors of the pads.

Order Indicator

Notifies of the current step while operating the equipment.

Quick Guide

Provides a brief guide on how to use the equipment during an emergency.

**Battery Pack** 

Supplies power with disposable batteries.

Infrared Data Association

(IrDA) Port


Transmits and receives data between the equipment and a PC.

CU Medical Systems, Inc.

# 3. Setting up the i-PAD

## **Package Contents**

The i-PAD NF1201 packaging box contains the following items.



i-PAD NF1201 Fully-automated External Defibriliator



Instructions for Use



Disposable, non-rechargeable Battery Pack




**Adult Defibrillator Pads** 

These items are the minimum requirements for the operation of the i-PAD . For optional accessories, see Appendix A.

## Setting up the i-PAD

- ① Upon opening the packaging box, verify that it contains all the items listed in the packing list.
- ② Familiarize yourself with the controls and features of the i-PAD. Study the functions of the buttons, switches, indicators, and connection ports.
- 3 Place the defibrillator pads package as shown in the following figure. Do not open the packaging pouch of the pads until you are about to use them.



4 Install the battery pack.





After battery insertion, the i-PAD runs an automatic self-test that verifies its readiness for rescue operations. Passing the battery insertion test is indicated by the State LED blinking in green.

# 3. Setting up the i-PAD

- ⑤ Place the i-PAD inside the carrying case if you have purchased this optional accessory.
- ⑥ Store the i-PAD in accordance with your emergency response protocol. The following must be taken into consideration when the storing the i-PAD.
- The conditions in the area must not exceed the environmental limits of the i-PAD. See the General Operating Guidelines in Chapter 6.
- It must be convenient to check the indicators of the i-PAD to ensure that t is always ready for emergencies.
- A telephone must be located close to the storage area so that you can easily call emergency medical service.



Electromagnetic interference may alter device performance.

During operation, the i-PAD should be placed away from sources of electromagnetic interference such as motors, generators, X-Ray equipment, radio transmitters, cellular mobile telephones and others, as these might interfere with the signals being acquired and analyzed.

See Appendix E. Electromagnetic Compatibility for details.

# **MARNING**

Using accessories and cables other than the ones specified in this manual may result in increased ELECTROMAGNETIC EMISSIONS or may decrease the ELECTROMAGNETIC IMMUNITY of the i-PAD.

Replacement accessories and consumables must be sourced only from CU Medical Systems, Inc. or its authorized representatives.

# A CAUTION

Do not connect the defibrillator pad assembly to the i-PAD during storage. Do not open the sealed container of the pads until ready for use to prevent the AED electrode pads from drying out.

# CU Medical Systems, Inc.

# 4. How to Use the Product

#### 4.1 Product Use Procedures

When a sudden cardiac arrest patient is found, perform the Chain of Survival that is recommended by American Heart Association (AHA).



- 1. Call 911 Request help from an emergency rescue team or a person with an equivalent qualification.
- 2. Early CPR Perform CPR.
- 3. Early defibrillation Use this product (i-PAD).
- 4. Early medical treatment Transfer a patient to a medical institute for treatment.
- 5. Integrated treatment after cardiac arrest Perform treatment after cardiac arrest at a specialized treatment facility.
- 3. Early defibrillation Use this product (i-PAD).

The procedures for the use of this product can be summarized in 3 steps as follows: After pressing the Power Button,

- Step 1: Place pads on the patient.
- Step 2: Automatic shock delivery if needed.
- Step 3: Perform CPR.

# Warnings

- ✓ Do not use this product with infants under 1 year old.
- ✓ When delivering shocks to children who are older than 1 year and younger than 8 years as well as lighter than 25 kg, use pediatric pads that are sold separately. However, do not use pediatric pads on adults.

## **∕ !** Warnings

✓ Do not place the patient in a wet location when performing an electric shock.

## Intelligent Public Access Defibrillator

# 4. How to Use the Product

# / Warnings

✓ Remove all other medical devices, which are not protected from the defibrillator, from the patient when performing an electric shock.

## / Warnings

- ✓ The user or other people near the patient should avoid a contact with the patient as follows:
- Do not touch any parts of the patient such as the patient's torso, head, arms, and/or legs.
- Do not touch conductive liquid such as gel, blood (body fluid) and/or saline solution.
- Do not touch conductive metallic objects such as a stretcher and/or wheelchair.

Such contacts may create an unexpected current path.

## 4.2 Preparation for Defibrillation

1 Turn the equipment on by pressing the Power Button




When the power is on, the equipment works as follows:


- Status Indicator: Blue lights on to indicate a self-diagnostic test
- Beeper: Rings for 1 to 2 seconds
- Status Indicator: Green lights on to indicate normal operation
- Voice Instruction: "Attach pads."

This voice instruction will repeat until the user places pads on the patient and connects pads to the equipment with the pad connectors. However, if pads are not connected within 10 minutes, the equipment will automatically turned off.

② Open the cover.



## ③ Remove the clothing of the patient





Do not waste time in removing the patient's clothing. If necessary, rip off or cut off the patient's clothing.

## 4 Tear open the packaging of the pads




## ⑤ Take the pads out of its packaging



# 4. How to Use the Product

**(6)** If the pad connectors are not connected to the pad insert, connect the pad connectors to the pad insert on the defibrillator.



7 Check where pads will be placed by referring to the pictures on the pads.





® Remove the protective film on the pads.



#### Main Rescue Sequence

### STEP 1: Attach the pads





The i-PAD prompts you with "Do not touch the patient" if it senses a connection to the patient.

# MARNING WARNING

The patient should be still during ECG signal acquisition and analysis to minimize motion artifacts in the signal.

## STEP 2: Automatic shock delivery if needed.

The i-PAD begins rhythm analysis as soon as it is connected to the patient.

The i-PAD activates the following Indicators and Voice prompt during rhythm analysis:

- Voice prompt: "Analyzing heart rhythm"
- Graphical Rescue Guide: lights up the second step of the rescue operation.

After rhythm analysis, the i-PAD decides whether the patient needs a shock or not.

# 4. How to Use the Product

#### If a shock is needed:

Indicators and Voice prompts:

- · Beeper: continuous beep
- · Voice prompt: "Stand by for automatic shock delivery. Warning: No one should be touching the patient."
- · Beeper: beeps continuously until the SHOCK deliver
- · Voice prompt: "Shock will be delivered in (3-2-1)"
- After the countdown (3 2 1), the i-PAD NF1201automatically delivers a shock
- It informs you of •• the shock delivery with the voice prompt: "Shock delivered"

If a shock is not needed:

#### Indicator:

· Voice prompt: "No shock advised"

#### NOTICE

While the i-PAD is charging after a shockable rhythm is detected, the ECG of the patient is continuously acquired and analyzed. If the ECG rhythm changes to a non shockable rhythm, the i-PAD disarms itself.

## MARNING.

Do not let anybody touch the patient during automatic shock delivery. Defibrillation shock can cause operator or by tander injury.

The i-PAD directs you to do CPR after delivering one shock

**Indicators and Voice Prompts** 

· Graphical Rescue Guide: lights up the third step of the rescue operation

 $\cdot$  For other indicators and voice prompts, see CPR Protocol

### **CPR Protocol**

· Voice Instruction: "Begin CPR, now"

· Beeper: Beeps for 1 second

· CPR Guide

	CPR Guide		
Category	Setting Range	Default	
Breathing / Not Breathing	Breathing guide; No breathing guide	Voice instruction mode; Pause mode	
Chest Compression Rate	100, 110, or 120 pushes per minute	110	
Default Guide Settings	Voice instruction mode; Pause mode	Pause mode	

## Note

- ✓ Refer to [CPR Guide Settings] on page 24 for detailed CPR guide settings.
- √Refer to [CPR Guide] on page 25 for detailed CPR guide.

## Intelligent Public Access Defibrillator

# 4. How to Use the Product

## **CPR Guide Settings**

This is the CPR Guide Settings mode. Set the guide as follows:

This is the CPR Guide Settings mode. Set the guide as follows:							
Step	Explanation	User's Action					
1	Enter CPR Guide Settings mode Voice Guide: CPR guidance setup mode.	- Enter administrator's mode by pressing the i-button for longer than 1 second. (Status Indicator: white lighting) - After a voice guide gives information on the version and usage of the equipment, press the i-button for longer than 3 seconds in order to enter the CPR guidance setup mode. (Status Indicator: purple lighting)					
2	Set Breathing Guide Voice Guide: Set the ventilation guidance. Voice Guide: Ventilation guidance on / off.	- Set to enable or disable the breathing guide Press the i-button to apply the setting and wait for 3 seconds. (When applied, the purple status indicator flashes twice)					
3	Set Compression Rate per Minute Voice Guide: Set the compression rate. Voice Guide: 100, 110, or 120 compressions per minute.	- Set a chest compression rate Press the Foutton to apply the setting and wait for 3 seconds. (When applied, the purple status indicator flashes twice)					
4	Set Default Voice Guide Voice Guide: Set to the default guidance. Voice Guide: Voice guidance mode. Pause time countdown mode.	Set to the default guidance for CPR.  Press the i-button to apply the setting and wait for 3 seconds. When applied, the purple status indicator flashes twice)  In voice guidance mode, a voice guide gives the CPR instructions.  When instructing CPR in pause time countdown mode, the guide pauses for 2 minutes and notifies the remaining paused time.(1 minute 30 seconds, 1 minute, 40 seconds, 20 seconds) If the user presses the i-button while it is flashing in blue within 20 seconds after pausing, the voice guide starts according to the settings above.					
5	Voice guide to either meet or save and exit after confirming settings: (Guide to the settings above) The CPR guidance settings will be modified to the following. Setting sentilation guidance. Ventilation guidance off. Setting compression rate. 110 compressions per minute. Setting to the default guidance. Pause time countdown mode.  Voice Guide: To shutdown and save the CPR guidance settings, press the power button. To change the CPR guidance settings, press the i-button	<ul> <li>After the voice guide announces the CPR guide settings, the user may either reset or save and exit at this step.</li> <li>After the settings guide, press i-button to reset settings and to go back to step 2.</li> <li>After the settings guide, press the power button to save the settings; a purple status indicator changes to white and flashes 3 times.</li> <li>After completing all steps, the equipment switches to standby mode.</li> </ul>					

## Note

✓ If no breathing guide is selected on the CPR guide settings, then a beat sound plays for 2 minutes at the set chest compression rate after giving a detailed description of chest compressions.

## CPR Guide

While performing CPR, the equipment gives a guide with a beat sound and voice instructions according to the user's CPR guide settings. The default on the equipment is pause mode. When the user presses the i-button within 20 seconds by following the default guide in pause mode, the voice guide explains how to do chest compressions and guides the user only to compress the chest to the beat for 2 minutes. The typical CPR order is as follows:

Step	Explanation	User's Action
1	Voice Instruction: "Begin CPR, now"	- Place both hands in the middle of the patient's chest between the nipples, and prepare for chest compressions.
2	Beat sound: depending on the setting, a beat sound plays 100, 110, or 120 beats per minute.	- According to the beat sound of the equipment, compress the patient's chest at least 5 on deep.
3	Beeper: briefly beeps once	While swiftly lifting the patient's chin up, tilt the head backward to open the airway.  - Prepare for artificial respiration.
4	Voice Instruction: Breathe~ Breathe~	- When performing artificial respiration, breathe longer than 1 second according to the "Breathe~" sound When performing artificial respiration, give the patient enough breath to make the chest significantly rise Breathe twice in 7 seconds.
5	Cycle: One cycle is composed of steps 1 through 4.	- Perform 5 cycles of CPR.  - When the chest compression rate is set to 120 pushes per minute, perform an additional chest compression after 5 cycles.

## Note

✓ When the default guide is set to pause mode, the equipment announces only the remaining time of CPR. (1 minute 30 seconds, 1 minute, 40 seconds, 20 seconds)

# 4. How to Use the Product

# **A** CAUTION

- The i-PAD temporarily stops analyzing the ECG of the patient during the duration of the CPR.
- It automatically resumes ECG analysis after CPR

# 

If it becomes necessary to use another defibrillator on the patient, do not leave the i-PAD connected to the patient. Disconnect the i-PAD from the patient before using any other defibrillator.



# 5.1 After Each Use

1. Check the i-PAD for signs of damage and contamination.

5. Product Care After Use

- Run the battery insertion test, see section on maintenance in Chapter 6. Verify that the State LED is blinking in green to signify that the i-PAD is ready for a rescue operation. If there is dirt contamination, see section on how to clean the i-PAD in Chapter 6.
- 2. Replace the AED electrode pads, see section on maintenance in Chapter 6. The pads are for single use only. Do not reuse the pads.



- · Use only the AED electrode pads recommended by CU Medical Systems, Inc.
- Do not open the packaging of the pads during storage. Open the packaging only when you are about to use them during rescue operations.

### 5.2 Data Storage and Transfer

#### 5.2.1 Last Use Data

The i-PAD automatically stores the following data during rescue operations.

- · ECG data
- · Rescue event data
- · setup information of the i-PAD

These data are stored in the internal memory of the i-PAD and may be transferred to a personal computer. The internal memory of the i-PAD is nonvolatile, thus, the stored data stay in memory even when the i-PAD is turned OFF

# **A** CAUTION

Do not uninstall the battery pack while the i-PAD is acquiring data. If you do so, you will lose the data in that particular rescue operation. If you must uninstall the battery pack, turn OFF the i-PAD properly by pressing the Power ON/OFF button before removing the battery pack.

# **Cautions**


- ✓ This equipment can save up to 75 minutes of ECG data. If the ECG data are recorded for more than 75 minutes, the ECG data after 75 minutes will not be recorded.
- ✓ This equipment saves the usage history and multiple ECG records. When using the product for less than 5 minutes each time, up to 7 entries can be saved.
- ✓ After saving the usage data for 7 times, new data replaces the oldest data. Hence, we recommend to transfer and save the data from the equipment to a PC after each use.

#### 5.2.2 Data Transfer

The data may be transferred to the PC by using the PC software (CU Expert Ver. 3.60 or higher) that is provided by us. [CU Expert] is a PC software that includes ECG review and printing functions.

#### Data Transfer

- ① Run [CU Expert] on a PC. Refer to the [CU Expert] manual for detailed information regarding how to receive data
- 2) Position the IrDA adapter to face the IrDA port on the equipment as shown in the figure below.
- ③ Press the i-button for longer than 1 second in standby mode (the status indicator flashes in green).



- 4) The white status indicator lights with a guide phrase "Administration mode."
- ⑤ The equipment gives the user a summary (the total hours of the last equipment use and the number of electric shocks).

# 5. Product Care After Use

- **(6)** Wait for 3 minutes when connecting to a PC. If settings on [CU Expert] are proper, the equipment connects to a PC in several seconds.
- 7 After transferring data, the equipment automatically turns off.

## 5.2.3 Equipment Settings

The user may set the volume on the equipment by using PC software (CU Expert Ver. 3.60 or higher). Refer to the [CU Expert] manual for details.

Equipment Settings	Default	Setting Range
Volume	10	0 ~ 10

# 6. Maintaining the i-PAD

### **General Operating Guidelines**

The following table shows the General Operating Guidelines of the i-PAD. Be sure that you do not subject the i-PAD to conditions that are beyond the limits specified below.



Do not operate or store the device in conditions that are beyond the following specified limits.

Standby conditions

Temperature 0  $^{\circ}$ C to 43  $^{\circ}$ C 32  $^{\circ}$ F to 109  $^{\circ}$ F) Humidity 5  $^{\circ}$ 8 to 95  $^{\circ}$ 8 (non-condensing)

Operating Conditions (The equipment contains pads combined with batteries; available immediately in emergency)

Temperature 0  $^{\circ}$  to 40  $^{\circ}$  C (32  $^{\circ}$ F to 104  $^{\circ}$ F) Humidity 5 % to 95 % (non-condensing)

Storage Conditions (The equipment does not contain pads and batteries; only the equipment is kept for an extended period of time or moved)


Temperature -20 °C to 60 °C (-4 °F to 140 °F) Humidity 5% to 95 % (non-condensing)



Do not store the device in areas that are directly exposed to sunlight.



Do not store the device in areas with highly fluctuating temperatures.



Do not store the device near heating equipment.



Do not store the device in areas where there is high vibration (in excess of Category 10 of MIL-STD-810E).



Do not operate or store the device in environments with high concentration of flammable gas or anesthetics.



Do not operate or store the device in areas with high concentration of dust.



Only personnel authorized by the manufacturer may open the device for servicing. There are no user serviceable components inside the device.

# 6. Maintaining the i-PAD

#### **Routine Maintenance**

#### **Device Status Monitoring**

The i-PAD conducts automated self-tests while on standby mode. It is on standby mode if the battery pack is inserted and the State LED is blinking in green. The automated self-tests are conducted daily, weekly, and monthly.

If a fault is detected during the automated self-tests, the i-PAD raises an alarm. See section on trouble-shooting in Chapter 7.

Periodically check the State LED of the i-PAD to ensure that it is always ready for an emergency.

#### Consumables

There are two important consumables that must be monitored while the i-PAD is stored on standby mode. These are: the battery pack and the AED electrode pads.

#### **Battery Pack**

- Replace the battery pack if a low battery level condition is indicated by the i-PAD.
- Use only battery packs that are recommended by the manufacturer.
- Make sure that the replacement battery pack is not yet past its expiration date which is indicated by the "Install by: date" marking.
- The battery pack of the i-PAD is disposable and must not be recharged.

#### **Battery Pack Replacement**

1. Remove the spent battery pack. Disengage its lock by pressing the locking mechanisms towards each other while simultaneously pulling the battery pack out of its compartment. These are shown below.





CU Medical Systems, Inc.

2. Insert the new battery pack with the label facing up and in the direction indicated by the guiding arrow printed on the label.



3. Push all the way in until you hear the locking mechanism click





# 

- · Do not charge the battery pack
- · Do not open the case of the battery pack.
- Do not saw off or break apart the case of the battery pack.
- Do not let the battery pack come into contact with open flames and other hot objects. Do not dispose of in fire
  - Do not short-circuit the terminals of the battery pack.
- Do not subject the battery pack to serious physical impact. Do not hit it with a hammer.
- In case of leakage or strange smell, keep away from fire to prevent ignition of any leaked electrolyte.

# 6. Maintaining the i-PAD

# **!** WARNING

- Keep the battery pack out of children's reach.
- · If the battery pack, leaks and the leaked liquid gets in the eyes, wash them with clean water and consult a physician immediately.
- Do not leave the battery pack in direct sunlight or in high temperature areas.
- · Do not have the battery pack in contact with water.
- Keep the battery pack away from direct sunlight, high temperature, and humidity.
- · Follow local regulations when disposing of the battery pack.
- Do not subject the battery pack to conditions beyond the safe environmental conditions for the i-PAD.

#### **AED Electrode Pads**


- Regularly check the expiration date of the AED electrode pads. Make sure that the set of AED electrode pads stored with the i-PAD is not expired.
- Check the integrity of the packaging of the AED electrode pads.
- · Use only AED electrode pads that are recommended by the manufacturer for use with the i-PAD.

#### **AED Electrode Pads Replacement**

1. Verify that the replacement pads are not expired. The expiration date is indicated using a sticker as shown below.



The sticker is located beside the AED Electrode Pads labeling "Multifunction Defibrillation ADULT PADS"




Details of the expiry date sticker Expiry Date YYYY-MM YYYY - Year

MM - Month

CU Medical Systems, Inc.

2. Place the AED electrode pads on top of the front panel of the i-PAD. Close the cover down on the pads.



## Cleaning the i-PAD

Clean the i-PAD after each use with a soft cloth which has been moistened with any one of these approved solutions:

Dilute soap and water

Dilute Chlorine bleach and water mixture (30 ml bleach/liter of water)

Dilute ammonia-based cleaners

Dilute hydrogen peroxide



- · Do not immerse any part of the i-PAD in fluids.
- If the i-PAD has been submerged in water, call CU Medical Systems, Inc. or an authorized representative for service assistance
- Do allow fluid to enter the case of the device.
- · Do not spill liquids on the case of the device.
- Do not use strong, acetone-based cleaners in cleaning the device.
- Do not use abrasive materials in cleaning the unit, especially on the infrared filter of the IrDA port.

Do not sterilize the i-PAD.

# 7. Troubleshooting

## **Automated Tests**

The i-PAD conducts the following tests to verify its readiness for emergency rescue operations:

Test	Test Details
	Runs when the battery is inserted.  Conduct this test:  When the i-PAD is first placed into service.  After the i-PAD is used in an emergency rescue operation.  When the battery pack is replaced.  When you suspect that your i-PAD is damaged.
	Do not run this test when you are about to conduct a rescue operation because this test takes time.
	If a new battery is inserted at the start of a rescue operation, press the Power ON/OFF button to bypass the test. Wait for the i-PAD to turn OFF. Press the ON/OFF button again to turn the i-PAD ON.
Battery insertion test	The i-PAD tests the i-button during this test; you must press the button when the i-PAD prompts you to do so.
	The following prompts are given:  "Press the flashing blue i-button"  If no fault is detected, the i-PAD goes into standby mode with the
16	State LED blinking in green.
	If a fault is detected, the i-PAD gives the prompt "Press the flashing red i-button". This is aside from activating the State LED and the beeper.
Power On Test	Runs when the i-PAD is turned ON.
Run-time test	Runs when the i-PAD is used in a rescue operation.  Monitors the performance of critical components.
Daily, Weekly, Monthly test	Runs daily, weekly, and monthly to check the functionality of important subsystems.

The blue status indicator lights in all self-diagnostic tests except "real-time self-diagnostic test."

When the equipment successfully completes a "battery pack self-diagnostics test" or "periodic self-diagnostics test", the status indicator flashes in green while the equipment switches to standby mode.

When a self-diagnostic test fails, the status indicator flashes in red. Also, the beeper beeps every minute. The equipment is unavailable for emergency while in this state.

When the i-button flashes in red, the user should press the i-button. After pressing it, the equipment will inform the user on the cause of an error via the voice instruction. Refer to [Section 7.3: Trouble-shooting] of this document for more information.

If the status indicator flashes in yellow after performing the "periodic self-diagnostic test", this indicates a problem with the pads such as a connection problem between the equipment and pads. Since the equipment does not have any problem, check pads and resolve relevant problems; either take out the battery and reinsert or turn on the equipment in order to perform a "battery pack self-diagnostic test" or "power self-diagnostic test." When problems are resolved, the status indicator flashes in green in standby mode.

## Intelligent Public Access Defibrillator

# 7. Troubleshooting

## **Equipment Status**

This equipment has different indicators to inform the user regarding the equipment conditions.

- Status Indicator
- Order Indicator
- Beeper
- i-Button

Indicator	Description	Remark
Status Indicator: Flashing in green	The equipment is available for a normal operation.	
Status Indicator: Green lighting	The equipment is properly functioning.	
Status Indicator: Flashing in red Beeper: Beeping 3 times in a 1 minute interval	The equipment detected an error during the self-diagnostic test; the power is off.	
Status Indicator: Red lighting i-Button: Flashing in red Voice Instruction: "Press the flashing red i-button"	The equipment detected an error while operating.	
Status Indicator: Blue lighting	The equipment is performing a self-diagnostic test.	
Status Indicator: White lighting	The equipment is currently on administrator mode.	
First, the order indicator is flashing	The pads are not connected to either the equipment or the patient.	
Second, the order indicator is flashing	The equipment is analyzing the patient's ECG.	
Third, the order indicator is flashing	The equipment is in CPR mode.	
i-Button: Flashing in red	The equipment detected an error. In order to identify the error, press the i-button.	
Status Indicator: Flashing in yellow	The pads are not properly connected to the equipment, or the pads have a problem.	

#### Troubleshooting the i-PAD

The i-PAD makes troubleshooting easy by providing indicators whenever

- $\cdot$  a fault in any of its subsystems is detected
- $\cdot$  or conditions that prevent the successful administration of a rescue operation occur.

Study carefully the following guides. Apply the recommendations when appropriate.

### **Troubleshooting During Emergencies**

Indication(s)	Possible Cause(s)	Remedial Action(s)
Device Mode; Standby Mode	· Low battery level.	Replace the battery pack of the i-PAD with a new battery pack
State LED: Flashing Red	· An error in the i-PAD is detected.	<ul> <li>Replace the i-PAD with a functioning defibrillator. Bring the malfunctioning device to a service center for repair.</li> </ul>
Device Mode; Operation Mode State LED: Flashing Red	· Low battery level.	<ul> <li>The i-PAD shuts down automatically after 10 minutes or after delivering 10 shocks.</li> <li>Replace the battery pack of the i-PAD with a new battery pack</li> <li>If no new battery pack or no other defibrillator is available, check the patient and begin CPR if needed until the EMS team arrives.</li> </ul>
Device Mode; Operation Mode State LED: Solid Red i-button: Flashing red	<ul> <li>An error in the i-PAD is detected.</li> <li>The i-PAD is unusable for a rescue operation in this condition.</li> </ul>	<ul> <li>Perform CPR on the patient.</li> <li>Use another defibrillator</li> <li>If there is no other defibrillator available, continue performing CPR until the EMS team arrives.</li> </ul>

# 7. Troubleshooting

### **Troubleshooting During Emergencies**

Indication(s)	Possible Cause(s)	Remedial Action(s)		
Voice Prompt: "No shock delivered." "Press the pads firmly to the bare skin"	The pads are not attached firmly to the skin of the patient.	<ul> <li>Press the pads firmly on the patient's skin.</li> <li>If needed, shave hair or wipe excess moisture off the skin of the patient and then reattach the pads.</li> <li>Get another set of pads immediately if the pads you are using do not stick to the bare and dry skin of the patient.</li> </ul>		

### **Troubleshooting Outside of Emergencies**

Indication(s)	Possible Cause(s)	Remedial Action(s)
Device Mode; Standby Mode State LED: Flashing Red	Low battery level.     An error in the i-PAD is detected.	<ul> <li>Press the Power ON/OFF button</li> <li>After turning ON, the i-PAD prompts you to:     "Press the flashing red i-button."</li> <li>After pressing the i-button, the i-PAD delivers any of the following prompts:     If the cause of theerror is low battery level, the prompt is:     "Low Battery Level, Replace the battery with a new one."     If the cause of the error is a fault in one of its subsystems, the prompt is:     "System failure. Error code is XXXX"</li> <li>If the cause of the error is low battery level, replace the battery pack. Use only battery packs supplied by CU Medical Systems, Inc.</li> <li>If the cause of the error is a system failure, bring the i-PAD to an</li> </ul>
Device Mode; Operation Mode State LED: Flashing Red	• Low battery level.	<ul> <li>authorized service center.</li> <li>The i-PAD shuts down automatically after 10 minutes or after delivering 10 shocks.</li> <li>Replace the battery pack of the i-PAD with a new battery pack</li> <li>Use only battery packs supplied by CU Medical Systems, Inc.</li> </ul>

# 7. Troubleshooting

### **Troubleshooting Outside of Emergencies**

Indication(s)	Possible Cause(s)	Remedial Action(s)		
Device Mode; Operation Mode State LED: Solid Red i-button: Flashing red	<ul> <li>An error in the i-PAD is detected.</li> <li>The i-PAD is unusable for a rescue operation in this condition.</li> </ul>	<ul> <li>The i-PAD prompts you to:     "Press the flashing red i-button"</li> <li>After you press the flashing red i-button, the i-PAD delivers the following prompt:     "System failure. Error code is XXXX"</li> <li>Bring the i-PAD to an authorized service center.</li> </ul>		

### For Non-emergency Situations

Symptom	Cause	Resolution
Equipment Mode: Standby Mode Status Indicator: Flashing in red	The battery is low. Or a system failure.	Press the power button. Press the i-Button by following the instruction, "Press the flashing red i-Button."  When the instruction says "Low battery level. Replace the battery with a new one." then replace the battery.  When the instruction says "System failure. Error code is OOOO,"then contact us by referring to [Product Service].
Equipment Mode: peration Mode Status Indicator: Flashing in red	The battery is low.	The equipment is turned off in 10 minutes or after 10 electric shocks.  Immediately replace the battery with a new one.
Equipment Mode: Operation Mode Status Indicator: Red lighting i-Button: Flashing in red	System failure.	Press the i-Button by following the instruction, "Press the flashing red i-Button." When the instruction says "System failure. Error code is OOOO," then contact us by referring to [Product Service].
Equipment Mode: Standby Mode Status Indicator: Flashing in yellow	The pads are not properly connected to the equipment, or the pads have a problem.	Check the connection between the pads and equipment, and connect the pad connectors to the pad insert on the equipment.  Reinsert the battery or turn on the equipment in order to perform a "battery pack self-diagnostic test" or "power self-diagnostic test." If problems are resolved, the status indicator flashes in green in standby mode. If the status indicator continues to flash in yellow in standby mode, replace the pads with new pads.

 $<sup>\</sup>checkmark$  If the problem is not solved with the resolution or if no battery is available to replace, contact us by referring to [Product Service].

# 8. Servicing the i-PAD

#### Warranty

- This device is warranted by CU Medical Systems, Inc. against defects in materials and workmanship for two full years from the date of original purchase. During the warranty period, we will repair or, at our option, replace at no charge a product that proves to be defective, provided you return the product, shipping prepaid, to CU Medical Systems, Inc. or its authorized representative.
- This warranty does not apply if the product has been damaged by accident or misuse or as the result of service or modification by other than CU Medical Systems, Inc. or its authorized representatives. IN NO EVENT SHALL CU MEDICAL SYSTEMS BE LIABLE FOR CONSEQUENTIAL DAMAGES.
- Only products with serial numbers and their accessories are covered under this warranty.
   PHYSICAL DAMAGE CAUSED BY MISUSE OR PHYSICAL ABUSE IS NOT COVERED UNDER THE WARRANTY. Items such as cables and modules without serial numbers are not covered under this warranty.

#### Warranty Disclaimer

- · Servicing by unauthorized personnel renders this warranty null and void.
- If the factory seal is broken without proper authorization from CU Medical Systems, Inc., this warranty becomes null and void.

#### Service

- The i-PAD must be serviced only by authorized personnel. Unauthorized servicing during the warranty period renders the warranty null and void.
- The i-PAD will be serviced free of charge during the warranty period. After the warranty period, the cost of material and service shall be shouldered by the user.
- When the i-PAD is not operating properly, immediately bring it for servicing to an authorized service center.
- Please fill up the following table with the necessary information when requesting for service.

Product cla	Product classification Fully Automated External Defibrillator			rillator
Product	Name	i-PAD	Model Name	NF1201
Serial N	lumber		Date of Purchase	
Sales Repr	esentative			
	Name			
User Information	Address			
	Contact 10.			
Brief descrip	otion of the			

#### How to Contact Us

Medical Systems, Inc.

5F Cheonggye Bldg., 221, Anyangpangyo-ro, Uiwang-si, Gyeonggi-do, Korea

TEL: +82 31 421 9700 / FAX: +82 31 421 9911 / Homepage: www.cu911.com

Sales inquiries: sales@cu911.com

Technical inquiries: techinfo@cu911.com

Service: service@cu911.com

## A. Accessories

#### A.1 Standard Accessories

- -. Adult defibrillator electorde Pads [CUA0903PF]
- -. Disposable Battery Pack [CUSA0601F]
- -. Instructions for Use

#### A.2 Optional Accessories

- -. Pediatric defibrillator electorde Pads [CUA0512P]
- -. Wall mounted cabinet
- -. Carrying case
- -. IrDA Adapter
- -. PC S/W (CU Expert Version 3.60 or higher)
- -. i-PAD Trainer

#### A.3 Service Center

#### **Customer Service**

Homepage: http://www.cu911.com

Address : 5F Cheonggye Bldg., 221, Anyangpangyo-ro,

Uiwang-si, Gyeonggi-do, Korea

Tel : +82-31-421-9700 Fax : +82-31-421-9911

#### **Product Purchase**

E-mail : sales@cu911.com
Tel : +82-31-421-9700

#### Service Request

E-mail : service@cu911.com Tel +82-31-421-9700

# B. Equipment Symbols

### B.1 i-PAD NF1201 Defibrillator

Symbol	Description
Turn Onen	The cover of the i-PAD opens if the yellow lid is lifted.
	Power ON/OFF button
i	i-button
	State LED
<b>┤☆</b> ├	BF type, defibrillation-proof equipment
$\triangle$	Attention: Refer to accompanying documents.
<b>C€</b> 0470	CE Mark
SN	Serial Number
	Date of manufacture

### B.2 i-PAD Packaging

Symbol	Description
	Stack up to 6 cartons high only
<u> </u>	This side up
<b>*</b>	Keep dry
Ţ	Fragile; breakable
<b>○</b> *	Use no hooks
-000	Temperature limits: -20 $^\circ$ C to 60 $^\circ$ C
<b>C</b> € 0470	CE Mark
SN	Serial Number

# B. Equipment Symbols

### **B.3 Symbols on Accessories**

### B.3.1 Disposable Battery Pack [CUSA0601F]

Symbol	Description
LiMnO <sub>2</sub>	Lithium Manganese Dioxide battery
LOT	Lot Number
	Date of manufacture
	Do not mutilate the battery or open the battery case
	Do not expose the battery to high heat or open flames.  Do not incinerate the battery.
	Do not crush the battery
	Do not dispose of the battery in municipal waste. Follow local regulations on battery disposal.
A	Attention: Refer to accompanying documents
<b>€</b> 0470	CE Mark

### B.3.2 Adult defibrillator electorde Pads [CUA0903PF]

Symbol	Description
00	Temperature limits: 0°C to 43°C
LOT	Lot Number
	Expiration date
REF	Order reference number
2	Single use only; do not reuse
	Do not fold or bend.
Contains no Latex	Contains no latex
2010-02	Expiration Date and Lot number sticker
$\triangle$	Attention: Refer to accompanying documents
C€	CE Mark

# C. Glossary

Patient

1 Cycle Refers to 30 chest compressions followed by 2 breaths during CPR.

1 CPR 1 CPR consists of 5 cycles.

American Heart Association (AHA) 2010 CPR Guidelines This product instructs to compress the chest immediately after one electric shock according to the 2010 CPR Guidelines. Also, the user may change

the CPR guide depending on the user's CPR proficiency. Please contact us

for additional questions.

Condensation If you take a cold water bottle out from the refrigerator, dewdrops form

on the surface of the water bottle. Likewise, since moisture has an adverse effect on the equipment when it is formed on the cold equipment surface by condensing moisture in the warm air, the Product should be stored in

a place without condensation.

Sudden Cardiac Arrest A patient with symptoms of a sudden cardiac arrest. This product is

intended to be used on patients with the following symptoms: First, this product should be used on patients with no response and no normal

breathing.

We Refers to CU Medical Systems, Inc.

Standby Mode In this mode, the status indicator flashes in green while the equipment is

turned off. The equipment is available for a normal operation.

Operation Mode In this mode, the green status indicator lights while the equipment is

turned on. The equipment is properly operating.

Battery Pack This equipment uses a disposable battery. Never try to recharge the

battery.

Status Indicator This indicates the current status of the equipment with various

colors and flashing.

Adult The adult referred in this manual is defined as a person who is older than

8 years or heavier than 25 kg.

Pediatric The child referred in this manual is defined as a person who is older than

1 year and younger than 8 years as well as lighter than 25 kg.

#### Intelligent Public Access Defibrillator

Self-Test Automated self diagnostic tests that verify the proper operation of

the subsystems of the i-PAD so that it is always ready for a rescue operation. The i-PAD indicates the error when you press the i-button  $\,$ 

after an error has occurred.

Internal discharge (disarm) While the i-PAD is charging after a shockable rhythm is detected, the ECG of

the patient is continuously acquired and analyzed. If the ECG rhythm

changes to a non shockable rhythm, the i-PAD disarms itself.

IrDA Port Port that is used to connect the i-PAD to a personal computer for data

transfer.

Shock The shock that is delivered to defibrillate the heart of a patient.

Charging Time: less than 12 seconds.

Solid The State LED is always ON.

Blinking The State LED is turned ON and OFF regularly.

i-button Last usage data or error codes or CPR guide prompts are played on

the speaker when this button is pressed.

Product The i-PAD defibrillator of CU Medical Systems, Inc. with model name

i-PAD.

Fully Automatic External

Defibrillator

A device that automatically defibrillating shock after analyzing and

recognizing a shockable rhythm.

Carrying Case A case that is used to store the i-PAD and all of the accessories

needed for a rescue operation.

Communications port The port that is used to connect the i-PAD to a personal computer for

data transfer.

AED electrode pads Electrode pads that are used in the acquisition of ECG and the

delivery of defibrillating shock to a patient.

Pads Connector The connector in the defibrillator pads assembly that is used to

connect the pads to the i-PAD.

## Instructions for Use **i-PAD**

# C. Glossary

**Pads Packaging** The pouch that is used to contain the pads. It prevents the pads from

drying out. Instructions for the usage of the pads are printed on this

packaging.

The coating on the electrodes that facilitates the conduction of Electrode pads conducting gel

electrical signals and energy between the patient's skin and the

electrodes.

Left leg electrode LL

RA Right arm electrode

The data management software that is used for the transfer of data PC S/W CU Expert from the i-PAD to a personal computer. Also used for reviewing, (CU-EX1)

printing, and archiving of data in the personal computer.

The port in the i-PAD that mates with the pads connector. Pads Connector Socket

**Emergency Medical Service EMS** 

# D. i-PAD NF1201 Specifications

Model: NF1201

Physical |

Category Nominal Specifications

Size 3.19 in high X 8.66 in wide X 11.06 in deep

(81 mm high X 220 mm wide X 281 mm deep)

- without handle and battery pack

Weight Approximately 4.85 lbs (2.2 kg) with battery pack installed

Environmental

Category Nominal Specifications

Operating Conditions (The equipment contains pads combined with batteries; available immediately in emergency)

Temperature 0  $^{\circ}$ C to 40  $^{\circ}$ C (32  $^{\circ}$ F to 104  $^{\circ}$ F) Humidity 5  $^{\circ}$ 6 to 95  $^{\circ}$ 6 (non-condensing)

Storage Conditions (The equipment does not contain pads and patteries; only the equipment is kept for an extended period of time or moved)

Temperature  $-20 \,^{\circ}\mathrm{C}$  to 60  $\,^{\circ}\mathrm{C}$  (-4  $\,^{\circ}\mathrm{F}$  to 140  $\,^{\circ}\mathrm{F}$ ) Humidity  $5 \,^{\circ}\mathrm{M}$  to 95  $\,^{\circ}\mathrm{M}$  (non-condensing)

Shock/Drop/Abuse

Tolerance

Meets IEC 60601-1 clause 21 (Mechanical Strength)

Vibration Meets EN1789 random and swept sine, road ambulance

specification in operating and standby states.

CU Medical Systems, Inc

Sealing IEC 60529: IP54

ESD Meets IEC 61000-4-2:2001

EMI (Radiated) Meets IEC 60601-1-2 limits, method EN 55011:1998+ A1:1999

+A2:2002, Group 1, Class B

EMI (Immunity) Meets IEC 60601-1-2 limits, method EN 61000-4-3: 2001 Level

3 (10V/m 80MHz to 2500MHz)

## Instructions for Use i-PAD

Defibrillator

Category Nominal Specifications

Operating Mode Fully automatic

Waveform  $e \sim cube$  biphasic (Truncated exponential type); impedance

compensated

Energy 200 Joules nominal into a  $50\Omega$  load

Charge Control Automatic by Software (Arrhythmia Detection System and

**Charging Control)** 

Charge time from "Stand by, for

automatic shock delivery."

< 10 seconds, typical

Charge complete indicator • Text prompt (SHOCK WILL BE DELIVERED IN, 3-2-1)

· Beep from the beeper

Disarm Once charged, the i-PAD disarms itself if:

• Patient's heart rhythm changes to non-shockable rhythm, or • The ON/OFF button is pressed to turn OFF the i-PAD, or

• The defibrillator pads are removed from the patient or the

pads connector is disconnected from the i-PAD

Automatic Shock Delivery Shock is automatically delivered if a shockable rhythm is detected.


Shock Delivery Vector Via adult defibrillator pads in the anterior-anterior (Lead II)

position or via Reduced-energy pediatric pads in the anterior-pos-

terior position.

Patient solation Type BF

# D. i-PAD NF1201 Specifications



Waveform Specifications (200 Joules)

Patient Impedance (Ohms)	Phase A, Duration (milliseconds)	Phase B, Duration (milliseconds)	Energy Delivered (Joules)	
25	1.9	1.9	200	
50	3.8	3.8	200	
75	5.7	5.7	200	
100	7.3	7.3	199	
125	9.2	9.2	199	
150	11.0	11.0	200	
175	12.8	12.8	200	

**ECG** Acquisition

Category Nominal Specifications

Acquired ECG Lead Lead II

Frequency Response 1 Hz to 30 Hz

**ECG Analysis System** 

Category Nominal Specifications

Function Determines the impedance of the patient and evaluates the ECG of

the patient to determine whether it is shockable or non shockable

Impedance Range  $25\Omega$  to  $175\Omega$ 

Shockable Rhythms Ventricular Fibrillation or Fast Ventricular Tachycardia

Fast Ventricular Tachycardia

Sensitivity & Specificity: Meets AAMI DF39 guidelines

# D. i-PAD NF1201 Specifications

### ECG Analysis System - ECG Database Test

ECG Rhythm Class	Rhythms	Minimum test sample size	Perfor- mance goal	Test sample size	Shock Decision	No Shock Decision	Observed Perfor- mance	90% One Sided Lower Confidence Limit
SHOCKABLE	Coarse VF	200	>90% sensitivity	219	213	6	97.26% (213/219) sensitivity	95%
SHOC	Fast VT	50	>75% sensitivity	137	111	26	81.02% (111/137) sensitivity	76%
	Normal Sinus Rhythm	100 minimum (arbitrary)	> 99% specificity	100	0	100	100% (100/100) specificity	97%
NON SHOCKABLE	AF,SB,SVT, heart block, idioven- tricular PVC'	30 (arbitrary)	> 95% specificity	219	1	218	99.54% (218/219) specificity	98%
	Asystole	100	> 95% specificity	132	5	127	96.21% (127/132) specificity	93%

## Instructions for Use i-PAD

Controls, Indicators, and Prompts

Category Nominal Specifications

Controls Power On/Off Button,

i-Button,

Indicators State LED, Graphical Rescue Guide LED

Audio Speaker Provides voice prompts

Beeper Provides various audible indications

Low Battery Detection Automatic during daily testing and Power ON and runtime

testing

Low Battery Indicator State LED and Voice Prompt

Prompts Voice prompts guide the user throughout a rescue operation

Self-Tests

Automatic Power On Self-Test / Run Time Self-Test

Daily / Weekly / Monthly

User Initiated Battery Insertion Test

Battery Pack [CUSA0601F]

Category Nominal Specifications

Battery Type 12 Volt DC, 4.2 Ah, lithium manganese dioxide, disposable

long-life primary cell.

Capacity Minimum 200 shocks or 4 hours of operating time.

· Operating Conditions
Temperature Range

Temperature: 0  $^{\circ}$ C to 40  $^{\circ}$ C (32  $^{\circ}$ F to 104  $^{\circ}$ F)

Storage Conditions

Temperature: -20 °C to 60 °C (-4 °F to 140 °F)

# D. i-PAD NF1201 Specifications

#### Defibrillator Pads (CUA0903PF)

Category Nominal Specifications

Type self-adhesive, disposable, non-polarized defibrillation pads

Adult Pads Defibrillation pads for patients 8 years of age and older or 55 lbs.

(25 kg) and over.

Surface Area Adult: 110cm each

Cable Length 1.5m

#### Defibrillation Pads for infants/children (CUA0512P)

Category Nominal Specifications

Type self-adhesive, disposable, non-polarized defibrillation pads

For infants/children Use on children up to 8 years old or up to 55lbs (25kg).

Surface Area For infants/children: 80cm per piece

Cable Length 1.5m

Data Recording and Transmission

Category Nominal Specifications

InfraRed Wireless transmission of event data to PC through IrDA port.

Data Stored Maximum 75 minutes of ECG and the entire incident's events and

analysis decisions.

# E. Electromagnetic Compatibility

Guidance and manufacturer's declaration - electromagnetic emissions

The i-PAD is intended for use in the electromagnetic environment specified below. The customer or the user of the i-PAD should assure that it is used in such an environment.

Emissions Test	Compliance	Electromagnetic environment - guidance
RF Emissions CISPR 11	Group 1	The i-PAD uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF Emissions CISPR 11	Class B	
Harmonic Emissions IEC 61000-3-2	Class A	The i-PAD is suitable for use in all establishments, including dome tic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic
Voltage fluctuations / flicker emissions IEC 61000-3-3	Complies	purposes.



The i-PAD should not be used adjacent to or stacked with other equipment.

If adjacent or stacked use is necessary, the NF1201 should be observed to verify normal operation in the configuration in which it will be used.

# E. Electromagnetic Compatibility

### Guidance and manufacturer's declaration - electromagnetic immunity

The i-PAD is intended for use in the electromagnetic environment specified below. The customer or the user of the i-PAD should assure that it is used in such an environment.

Immunity Test	IEC 60601-1 test level	Compliance level	Electromagnetic environment - guidance	
Electrostatic discharge (ESD)	±6 kV contact	±6 kV contact	Floors should be wood, concrete or ce-	
IEC 61000-4-2	±8 kV air	±8 kV air	with synthetic material, the relative humidity should be at least 30%.	
Electrical fast transient/burst	$\pm 2$ kV for power supply lines	$\pm 2$ kV for power supply lines	Mains power quality should be that of a typical commercial or hospital environment.	
IEC 61000-4-4	$\pm 1$ kV for input/output lines	±1 kV for input/output lines		
Surge	$\pm 1$ kV differential mode	±1 kV differential mode	Mains power quality should be that of a typical commercial or hospital environment.	
IEC 61000-4-5	±2 kV common mode	±2 kV common mode		
Voltage dips, short interruptions and voltage variations on power supply input lines	<5 % Uτ (>95% dip in Uτ) for 0,5 cycles  40 % Uτ (60% dip in Uτ) for 5 cycles  70 % Uτ (30% dip in Uτ)	<5 % Uτ (>95% dip in Uτ) for 0,5 cycles 40 % Uτ (60% dip in Uτ) for 5 cycles 70 % Uτ (30% dip in Uτ)	Mains power quality should be that of a typical commercial or hospital environment. If the user of the i-PAD requires continued operation during power mains interruptions, it is recommended that the i-PAD be powered from an uninterruptible power supply or a battery.	
	for 25 cycles <5 % UT (>95% dip in UT) for 0,5 cycles	for 25 cycles <5 % Uτ (>95% dip in Uτ) for 0,5 cycles		
Power frequency (50/60 Hz) magnetic field IEC 6 1000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.	

### Guidance and manufacturer's declaration - electromagnetic immunity

The i-PAD is intended for use in the electromagnetic environment specified below. The customer or the user of the i-PAD should assure that it is used in such an environment.

Immunity Test	IEC 60601-1 test level	Compliance level	Electromagnetic environment - guidance
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz outside ISM bands <sup>a</sup> 10 Vrms 150 kHz to 80 MHz in ISM bands	3 Vrms	Portable and mobile RF communications equipment should be used no closer to any part of the i-PAD, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance $d=1.16\sqrt{P}$ $d=1.2\sqrt{P}$
Radiated RF IEC 61000-4-3	10 V/m 80 MHz to 2,5 GHz	10 V/m	$d=1.2\sqrt{P}$ 80MHz to 800MHz $a=2.3\sqrt{P}$ 800MHz to 2,5GHz
	ic.		where is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and <b>d</b> is the recommended separation distance in metres (m) <sup>b</sup> Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey <sup>c</sup> , should be less than the compliance level in each frequency range <sup>d</sup> .  Interference may occur in the vicinity of equipment marked with the following symbol:

NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

The ISM (industrial, scientific and medical) bands between 150 kHz and 80 MHz are 6,765 MHz to 6,795 MHz; 13,553 MHz to 13.557 MHz; 26,957 MHz to 27,283 MHz; and 40,66 MHz to 40,70 MHz.

The compliance levels in the ISM frequency bands between 150 kHz and 80 MHz and in the frequency range 80 MHz to 2,5 GHz are intended to decrease the likelihood that mobile/portable communications equipment could cause interference if it is inadvertently brought into patient areas. For this reason, an additional factor of 10/3 is used in calculating the recommended separation distance for transmitters in these frequency ranges.

Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the i-PAD is used exceeds the applicable RF compliance level above, the i-PAD should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as or relocating the i-PAD.

Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

# E. Electromagnetic Compatibility

## Recommended separation distances between portable and mobile RF communications equipment and the i-PAD

The i-PAD is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the i-PAD can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the i-PAD as recommended below, according to the maximum output power of the communications equipment.

Rated maximum	Separation distance according to frequency of transmitter  m					
output power	150 kHz to	150 kHz to	80 MHz to	800 MHz to		
of transmitter	80 MHz	80 MHz	800 MHz	2,5 GHz		
W	outside ISM 	in ISM bands				
	bands	_		_		
	$d = 1.16\sqrt{P}$	$d = 1.2\sqrt{P}$	$d = 1.2\sqrt{P}$	$d = 2.3\sqrt{P}$		
0.01	0.116 m	0.12 m	0.12 m	0.23 m		
0.1	0.37 m	0.38 m	0.38 m	0.73 m		
1	1.16 m	1.2 m	1.2 m	2.3 m		
10	3.67 m	3.79 m	3.79 m	7.27 m		
100	11.6 m	12 m	12 m	23 m		

For transmitters rated at a maximum output power not listed above, the recommended separation distance **d** in meters (m) can be determined using the equation applicable to the frequency of the transmitter, where **P** is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2 The ISM (industrial, scientific and medical) bands between 150 kHz and 80 MHz are 6,765 MHz to 6,795 MHz; 13,553 MHz to 13,567 MHz; 26,957 MHz to 27,283 MHz; and 40,66 MHz to 40,70 MHz.

NOTE 3 An additional factor of 10/3 is used in calculating the recommended separation distance for transmitters in the ISM frequency bands between 150 kHz and 80 MHz and in the frequency range 80 MHz to 2,5 GHz to decrease the likelihood that mobile/portable communications equipment could cause interference if it is inadvertently brought into patient areas.

NOTE 4 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.





### CU Medical Systems, Inc.

130-1, Dongwhagongdan-ro, Munmak-eup, Wonju-si, Gangwon-do, Korea

TEL: +82 33 747 7657 / FAX: +82 33 747 7659

Homepage: www.cu911.com

-Ver 1.10-