

Product Specification

Version: 02

Date of issue: 22/02/2017

Date of revision: 26/10/2017

Medipal[®]
Infection prevention. Everywhere.

Medipal Alcohol Disinfectant Wipes

Medipal Alcohol Wipes are intended for cleaning and disinfecting hard surfaces of non-invasive medical devices within healthcare and care environments. For professional use only.

Medipal Alcohol wipes are proven to kill most microorganisms such as bacteria & viruses. The wipes have been independently tested to a range of European testing standards (EN) which includes EN16615.



Product Description

Product Name	150 Medipal Alcohol Wipes
Product Code	W505110MPCE
Languages	EN; FR; AE/SA

Pack Type & Dimensions

Count	150
Dimensions	130 x 150mm
Weight	0.279kg
Colour of wipes	White
Wipes material description	20gsm Tritex Technology

Solution

Solution Description	Medipal Alcohol Solution
Active Ingredients Content	Propan-2-ol
Hazard identification	H225: Highly Flammable Liquid & Vapour For further information please refer to the material safety data sheet

Primary Packaging

Pack Type	0.5 Litre HDPE canister
Lid Description	Red PP dispenser lid
Label	150 Alcohol Wipes label
Barcode	5025254 019214

Product Specification

Version: 02

Date of issue: 22/02/2017

Date of revision: 26/10/2017

Medipal[®]
Infection prevention. Everywhere.

Medipal Alcohol Disinfectant Wipes

Secondary Packaging

Case Description	Plain outer case
Barcode	5025254 042199

Packaging & Logistical Information

Product Name	150 Medipal Alcohol Wipes
Product Code	W505110MPCE
Case Quantity	10 canisters/case
Case Dimensions [OD]	435 x 185 x 180
Case Weight	3.1kg
Case/Layer	12
Layers/Pallet	7
Cases/Pallet	84
Pallet Weight	275kg
Hazardous for Transport	Not applicable
Commodity code for export	3808949000

