

Product Data Sheet



Bulbous Male Urinal 875ml

Disposable, single use bulbous male urinals. Ergonomically shaped traditional design, enabling hygienic, independent toileting in bed for male patients.

Key Benefits:

- Reduces the risk of cross contamination and for increased infection prevention.
- Single use with a minimum of 4 hours water retention.
- Ergonomically shaped. Easy to carry and hold.
- Environmentally friendly recyclable cardboard packaging
- Products stack well for storage purposes and separate easily for dispensing.



Product specification

Manufacturer's Product Code	9060
Unit of Issue	Case of 100 units
Usable Capacity	800 ml
Product Dimensions (L x W x H) mm	260 x 100 x 130
Product Weight (g)	35 (+/- 2g)
Outer Case Dimensions (L x W x H) cm	89.5 x 25.0 x 88.5
Outer Case Weight (kg)	3.50

Product information

Composition	100% recyclable and biodegradable Moulded Pulp Paper, consisting of clean Over-Issued News Print with fixative, waterproofing and detergent proofing additives.
Shelf Life	Minimum 4 years if stored in cool, dry conditions
Sterilisation	Non-sterile
Standards and Guidelines	PAS29:1999 Disposable pulp products for use in healthcare (Maceration test, normal cycle) PAS29:1999 Disposable pulp products for use in healthcare Clause 5 (Water retention test at 35°C, +/- 3°C) ISO9001:2008 BRC/IoP Global Standard
Disposal Method	Macerator or in accordance with local waste policy
Country of Origin	United Kingdom