

3M™ Littmann® Stethoscopes  
For Medical Professionals



# The First Test of Commitment to Quality Patient Care



Choosing a stethoscope is one of the first steps you take on the path to becoming a physician. Select a Littmann stethoscope and you are choosing wisely. Littmann stethoscopes are trusted, precision-crafted patient assessment instruments, renowned worldwide for quality and acoustical performance. Two examples:

**3M™ Littmann® Cardiology III™ Stethoscope** is a versatile double-sided stethoscope: a larger adult side and a smaller pediatric side with a tunable diaphragm on each. You can also remove the diaphragm and use as an open bell.

**3M™ Littmann® Master Cardiology™ Stethoscope** is a high-performance stethoscope featuring rugged durability and tunable diaphragm technology. Alternate between low and high frequency sounds simply by pressing the chestpiece.

You deserve the best and so do your patients. Invest in a Littmann cardiology-grade stethoscope now to enhance your learning ability and the care you provide to your patients.



3M™ Littmann®  
Cardiology III™  
Stethoscope



3M™ Littmann®  
Master Cardiology™  
Stethoscope





**FREE Student Accessory Kit**  
with the purchase of a Cardiology or Electronic Stethoscope

# Uniquely Yours

## Littmann cardiology stethoscopes offer:

- Tunable diaphragm technology, enabling you to conveniently alternate between low- and high-frequency sounds using only the diaphragm side of the chestpiece
- Single tube with a double lumen (two channels) design eliminates frictional tube rubbing noise of old dual tubing designs
- Comfortable fit and excellent acoustic seal with 3M™ Littmann® Snap Tight Soft-Sealing Eartips
- Anatomically angled headset for improved sound
- Patient-friendly, nonchill rim and diaphragm
- Many tubing colors and chestpiece finishes allow you to make your Littmann stethoscope your own



**Because your most vital instrument is knowledge**

Basic content for all.  
Premium content for Littmann stethoscope owners.

Littmann stethoscopes invests in and supports your educational needs by offering valuable auscultation training at [www.Littmann.com/education](http://www.Littmann.com/education).

To learn more about Littmann stethoscopes, visit [www.Littmann.com](http://www.Littmann.com) or call the 3M Health Care Helpline at 1-800-228-3957.



**3M Health Care**  
**Littmann® Brand Products**  
2510 Conway Avenue  
St. Paul, MN 55144-1000  
U.S.A.  
1 800 228 3957  
[www.3m.com/healthcare](http://www.3m.com/healthcare)  
[www.Littmann.com](http://www.Littmann.com)

**3M Health Care**  
**3M Canada**  
Post Office Box 5757  
London, Ontario N6A 4T1  
Canada  
1 800 364 3577

[www.3m.com/patents](http://www.3m.com/patents)  
3M, Cardiology III and Master Cardiology are trademarks of 3M. Littmann is a registered trademark of 3M. Used under license in Canada. Please recycle. Printed in U.S.A. © 3M 2014. All rights reserved. 70-2010-7282-7





3M™ Littmann®  
Cardiology III™  
Stethoscope



3M™ Littmann®  
Master  
Cardiology™  
Stethoscope



|  |  |  |
|--|--|--|
| <b>Performance (1-10 Scale)</b>              | 9  | 10   |
| <b>Chestpiece Design and Technology</b>      | Machined Stainless Steel<br>Two-sided Tunable  | Polished Stainless Steel<br>Single-sided   |
| <b>Headset Material</b>                      | Wide Diameter<br>Aerospace Alloy   | Wide Diameter<br>Aerospace Alloy   |
| <b>Length</b>                                | 22"/56 cm, 27"/69 cm   | 22"/56 cm, 27"/69 cm   |
| <b>Total Weight</b>                          | 165, 175 grams   | 175, 185 grams   |
| <b>Diaphragm Details, Diameter</b>           | Two Tunable Diaphragms<br>1.7"/4,3 cm, 1.3"/3,3 cm   | Tunable<br>1.75"/4,4 cm  |
| <b>Warranty</b>                              | 5 Years  | 7 Years  |
| <b>Available finishes* &amp; Tube Colors</b> |  Black <sup>1</sup>  Brass <sup>2</sup>  Copper <sup>3</sup><br><i>* Special chestpiece finishes only available with certain tube colors, as indicated.</i>  Smoke <sup>4</sup>  Rainbow <sup>5</sup> |  Black <sup>1</sup>  Brass <sup>2</sup>  Smoke <sup>4</sup> |
| <b>Black</b>                                 | 3127, 3128, 3131BE <sup>1</sup> , 3128BRS <sup>2</sup> , 3157SM <sup>4</sup> , 3152RBW <sup>5</sup>  | 2159, 2160, 2161 <sup>1</sup> , 2175 <sup>2</sup> , 2176 <sup>4</sup>  |
| <b>Burgundy</b>                              | 3129   | 2163   |
| <b>Navy Blue</b>                             | 3130   | 2164   |
| <b>Hunter Green</b>                          | 3134   | 2165, 2183 <sup>2</sup>  |
| <b>Plum</b>                                  | 3135   | 2167   |
| <b>Chocolate</b>                             | 3137, 3137CPR <sup>3</sup>   |  |
| <b>Gray</b>                                  | 3136   |  |
| <b>Caribbean Blue</b>                        | 3138   | 2178   |
| <b>Ceil Blue</b>                             | 3146   |  |
| <b>Red</b>                                   | 3140   |  |
| <b>Raspberry</b>                             | 3148   |  |
| <b>Coral</b>                                 | 3149   |  |
| <b>Sand</b>                                  | 3165 <sup>4</sup>  |  |
| <b>Dark Olive Green</b>                      | 3166 <sup>4</sup>  | 2182 <sup>4</sup>  |
| <b>Lavender</b>                              | 3158 <sup>5</sup>  |  |
| <b>Turquoise</b>                             | 3159 <sup>4</sup>  |  |
| <b>Lemon-Lime</b>                            | 3161 <sup>4</sup>  |  |

<sup>1</sup> Black Edition – black finish chestpiece and headset

<sup>2</sup> Brass Edition – brass finish chestpiece and headset

<sup>3</sup> Copper Edition – copper finish chestpiece and headset

<sup>4</sup> Smoke Edition – smoke finish chestpiece and headset

<sup>5</sup> Rainbow Edition – rainbow finish chestpiece and sunshine gold headset