

FERNO Combi Head Immobiliser (CHI)

Developed for all Scoop® stretchers and
Spineboards this Head Immobiliser
is now indeed UNIVERSAL



**FOR MORE DETAILED
INFORMATION SEE OVER**

Specifications

Length	Width	Height	Weight
250 mm	320 mm	170 mm	2.5 kg/5.5lbs

FERNO – is a world leading innovator and manufacturer of patient transportation solutions.

Our product range is born of a closer understanding of our customers' needs and the practical demands of their work. **Ferno's Combi Head Immobiliser** typifies our passion for top quality, highly functional design... patient transfer products that don't just do the job, they do it **BETTER**.

call +44 (0) 1274 851 999
email sales@ferno.co.uk
www ferno.co.uk

 **FERNO**

FERNO Combi Head Immobiliser (CHI)

Specialty designed for optimum fitting on Scoop® stretchers and Spineboards

Ferno's classic **Universal Head Immobiliser** design led the way in providing secure immobilisation, preventing further injury to a patient's head and neck. This tried and tested product has been enhanced so it can be used on all versions of Spineboard and Scoop® or breakaway stretcher. The Ferno **Combi Head Immobiliser** has a uniquely designed basepad and strapping system providing optimum fitting in all cases, therefore making it the first choice for all medical professionals.

A further advantage of the Ferno **CHI** is the unique latex-free Polyurethane moulded head-blocks, which are impervious to liquids and also contain anti-bacterial additives. Accidental surface damage will not affect the integrity of the product, which means enhanced durability and many years of effective use.

The **CHI** remains structurally and functionally similar to the original design, which ensures there is little or no staff retraining required.


- **Highly versatile** – compatible with all versions of Scoops and Spineboards, by altering orientation of head-blocks and using a specially designed strapping system
- **Combat infection** – fully infection-control compliant, using anti-microbial Velcro and PU moulded head-blocks that contain anti-bacterial additives throughout
- **Reusable** – accidental surface damage to head blocks does not affect integrity of product, and can still be reused, unlike dipped or disposable head blocks
- **Easy storage** - compact in design, can be stored whilst located on a spineboard or separately in a cupboard
- **Modular design** – allows replacement of individual parts, not the complete product
- **Use in X-ray and MRI scanners** – scans of the patient can be taken with the **CHI** in position
- **Quick to deploy** – easy to apply the **CHI** to a Scoop or Spineboard in minutes

Part Nos:

031271300	Ferno CHI Combi Head Immobiliser
031271600	Set of replacement Grey AB fixing straps for CHI (2 * headblock H-straps plus board fixing straps)
031271500	Replacement Orange head/chin pad (each)
0451-1013	Replacement headblock for CHI (each)
031271400	Spare/replacement baseplate for CHI



How straps work with
Millennia board



CHI FPC004 08.22

call +44 (0) 1274 851 999
email sales@ferno.co.uk
www ferno.co.uk



Ferno (UK) Limited, Stubs Beck Lane, Cleckheaton, West Yorkshire, BD19 4TZ, England