

# SAFETY DATA SHEET

## Clinell Universal Spray

According to Regulation (EU) No 453/2010

Issue Date: 25<sup>th</sup> April 2016

Version Number: 7

### SECTION 1: Identification of the substance/mixture and company/undertaking

#### 1.1 Product Identifier

Product Name Clinell Universal Spray

#### 1.2 Relevant identified uses of the substance or mixture and uses advised against

Identified Use Disinfectant spray for surface disinfection and cleaning of non-invasive medical devices

#### 1.3 Details of the supplier of the safety data sheet

Supplier GAMA Healthcare Ltd  
2 Regal Way  
Watford  
Hertfordshire WD24 4YJ  
United Kingdom  
Tel: +44 (0) 207 993 0030  
Email: [info@gamahealthcare.com](mailto:info@gamahealthcare.com)

#### 1.4 Emergency telephone number

Tel: +44 (0) 207 9930 035

### SECTION 2: Hazards identification

#### 2.1 Classification of the substance or mixture

Classification according to Regulation (EC) No 1272/2008 Mixture not classified as hazardous

#### 2.2 Label Elements

Contains PHMB. May produce an allergic reaction.

### SECTION 3: Composition/information on ingredients

#### 3.1 Mixtures

Declarable components	Conc.(%)	EC No.	CAS No.	Classification of individual components under Regulation EC No1272/2008
Benzalkonium chloride	≤0.5	270-325-2	68424-85-1	Skin Corr 1B (H314) Acute Tox 4 (H302, H312) Aquatic Acute 1 (H400)
Didecyl dimethyl ammonium chloride	≤0.5	230-525-2	7173-51-5	Acute Tox 4 (H302) Skin Corr 1B (H314)
Polyhexamethylene biguanide (PHMB)	≤0.10	NA	27083-27-8	Acute tox 4 (H302) Skin sens 1B (H317) Eye dam 1 (H318) Carc. 2 (H351) STOT RE 1 (H372) Aquatic acute 1 (H400) Aquatic chronic (H410)

Other components:

Water >75  
Additives Each <1

# SAFETY DATA SHEET

## Clinell Universal Spray

According to Regulation (EU) No 453/2010

Issue Date: 25<sup>th</sup> April 2016

Version Number: 7

### SECTION 4: First aid measures

#### 4.1 Description of first aid measures

##### Inhalation

Acute effects following exposure to this product via the inhalation route are not anticipated during normal handling and use.

##### Skin

This product is not intended for skin use. However, it has been dermatologically tested and approved safe for contact with skin. The use of gloves is recommended for prolonged use. If irritation develops, seek medical advice.

Although this product contains components classified as corrosive and sensitive to skin, due to the high volume of water also present in the formulation, the dilution effect means the classification of the formulation through CLP does not result in the hazards being carried through to the product.

A toxicological risk assessment considers this product unlikely to cause dermal irritation.

##### Eye

This product contains components classified as damaging to eyes. Due to the high volume of water also present in the formulation, the dilution effect means that the classification through CLP does not result in the hazard being carried through to the product.

Nevertheless, should eye irritation be experienced, this effect would likely be transient. But should symptoms persist, seek medical advice.

##### Ingestion

This product is for external use only and should be kept away from children. No adverse effects are anticipated from the formulation via the oral route during normal handling and use.

#### 4.2 Most important symptoms and effects, both acute and delayed

This product contains PHMB which may produce an allergic reaction.

#### 4.3 Indication of any immediate medical attention and special treatment needed

Treat symptoms as they occur

### SECTION 5: Firefighting measures

#### 5.1 Extinguishing media

Water spray, carbon dioxide, dry chemical and foam are compatible with the product. No unsuitable extinguishing media are known.

#### 5.2 Special hazards arising from the substance of mixture

The product is water based, therefore not flammable or explosive.

#### 5.3 Advice for fire fighters

Fire fighters should wear an approved self-contained breathing apparatus and full protective clothing.



# SAFETY DATA SHEET

## Clinell Universal Spray

According to Regulation (EU) No 453/2010

Issue Date: 25<sup>th</sup> April 2016

Version Number: 7

Initial boiling point/range	Ca. 100°C
Flash point	Not determined: water based product
Evaporation rate	Not determined: water-based product
Flammability (solid, gas)	Not determined: water- based product
Flammability or explosive limits	Not determined: water-based product
Vapour pressure	24 mmHg (25°C) (water)
Relative density	Not determined: water-based product
Solubility	Liquid is water soluble
Partition coef	No data available
Auto-ignition temperature	Not determined: water based product
Decomposition temperature	Not determined: water based product
Viscosity	Not determined: water based product
Explosive properties	Not determined: water based product
Oxidising properties	Not determined: water based product
<b>9.2 Other information</b>	Not available

### SECTION 10: Stability and reactivity

#### 10.1 Reactivity

Contact with ionic substances for example oils and dyes, may reduce the effectiveness of the product. Contact with oxidising agents should be avoided.

#### 10.2 Chemical stability

This product is considered stable under normal ambient storage and handling conditions or temperature and pressure.

#### 10.3 Possibility of hazardous reactions

No hazardous reactions anticipated.

#### 10.4 Conditions to avoid

None known

#### 10.5 Incompatible materials

Oxidizing agents and anionic formulations.

#### 10.6 Hazardous decomposition products

None known.

# SAFETY DATA SHEET

## Clinell Universal Spray

According to Regulation (EU) No 453/2010

Issue Date: 25<sup>th</sup> April 2016

Version Number: 7

---

### SECTION 11: Toxicological information

---

This preparation has undergone a toxicology risk assessment.

#### 11.1 Information of toxicological effects

Acute toxicity

Not likely to be acutely toxic.

Irritancy

Not likely to cause significant dermal irritation.

Corrosivity

No risk of dermal corrosivity identified under normal handling and use.

Sensitisation

Not likely to cause significant sensitisation or delayed hypersensitivity.

Repeated dose toxicity

No data available on the repeat dose toxicity of this product.

Carcinogenicity

No data on the carcinogenicity of this product.

Mutagenicity

None of the components in the formulation have exhibited confirmed mutagenic characteristics in the evaluation of their toxicity to date.

Toxicity for reproduction

No data available on the toxicity to reproduction of this product.

---

### SECTION 12: Ecological information

---

Ecotoxicological data have not been determined specifically for this product. Based on the classification of the formulation through CLP, the environmental hazards are not carried through to the product.

#### 12.1 Toxicity

Components are classified as toxic to the environment but are not present in the formulation at sufficient levels. The hazards are not carried through to the product.

#### 12.2 Persistence and degradability

Two components of the formulation (DDAC and BAC) have been found to readily degrade in OECD 301D closed bottle tests. However, PHMB was found not to be readily biodegradable under the same protocol.

#### 12.3 Bioaccumulative potential

Due to the distribution coefficient of n-octanol/water, accumulation in organisms is not expected.

#### 12.4 Mobility soil

No information available on the mobility of active substances in soil.

# SAFETY DATA SHEET

## Clinell Universal Spray

According to Regulation (EU) No 453/2010

Issue Date: 25<sup>th</sup> April 2016

Version Number: 7

### 12.5 Results of PBT and vPVP assessment

The formulation does not contain substances that meet the PBT or vPvB criteria of REACH annex XIII.

### 12.6 Other adverse effects

No information available

---

## SECTION 13: Disposal considerations

---

### 13.1 Waste treatment methods

This product may be disposed of via the drains, by landfill, or by incineration. Disposal must be in accordance with current national and local regulations.

In the Healthcare Industry, chemical residues, biocides and infectious substances generated as a result of medical and nursing care may require classification as hazardous waste.

Waste disposal is regulated in the EC member countries through corresponding laws and regulations. In the UK, we recommend that you consult the List of Wastes available through the environment agency. In other countries, contact either the authorities or approved waste disposal companies for advice on disposal of used waste.

---

## SECTION 14: Transport Information

---

Not classified for transport

---

## SECTION 15: Regulatory information

---

### 15.1 Safety, health and environmental regulations/legislation specific for the mixture

The product is classified under the Classification, Labelling and Packaging of Substances and Mixtures (EC) No 1272/2008, it contains substances which are notified under the Biocidal Products Regulation (EU) No 528/2012.

### 15.2 Chemical safety assessment

Not applicable

---

## SECTION 16: Other Information

---

Revisions

Currently in sixth version to bring in line with new regulations.

Basis of classification

The mixture is self-classified on the basis of available information on the ingredients.

This Safety Data Sheet was compiled using ECHA Guidance on the compilation of Safety Data Sheets, Version 1.1 December 2011.

*Disclaimer: This information is furnished without warranty express or implied, except that it is accurate to the best of our knowledge.*

*It relates to the specific material designated herein, and does not relate to the use in combination with any other material or in any process. We assume no legal responsibility for use or reliance upon this information.*