MARSDEN

M-200

Marsden High Capacity Chair Scale











Premium, Class III Approved Chair Scale

Designed to accommodate bariatric patients

Robust construction with comfortable one-piece seat

Rechargeable battery: 80 hours' use between charges

www.marsden-weighing.co.uk



Marsden M-200: A high capacity, high quality chair scale

The M-200 is Marsden's flagship chair scale. **Robustly built**, with **premium features** and an **extra-wide seat**, this chair scale will weigh **almost any patient** - providing reliable weighings for many years.

High accuracy, high capacity

The M-200 is available as either a **250kg or 300kg** chair scale - meaning it can be used to weigh almost any patient.

Graduations start at just **50g**, meaning you'll also get a precise reading of your patient's weight. This means this **Class III Approved** scale can be used for medical diagnosis - as well as gaining a weight for critical treatment - for patients unable to stand on a floor scale.



The M-200 boasts a powder-coated frame to make it more **hard-wearing**. The **extra-wide seat** is a one-piece moulding and provides comfort for the patient.

To make the scale easier to get into (and out of), the **armrests and footrests are hinged**. Brakes are provided on two of the wheels.

An **optional handle** is also available for the back of the scale, for improved control when pushing the scale.





Easy to use features

The M-200 features the popular Marsden DP-3810 indicator. This means its features and functionality is identical to that of many of our other medical scales, so it's **quicker to get to grips with**.

A **full numeric keypad** allows you to **calculate BMI** quickly and easily. Tare allows you to remove unwanted weight and Hold stabilises the weight reading on the clear display. Clear **dual displays** mean weight, height and BMI data is easy to read, instantly

Key features

Class III Approved
Heavy duty construction
One-piece extra wide seat
Hinged armrests/footrests
Brakes on two wheels
Large, bright display
BMI/Hold/Tare

Use it anywhere

The internal rechargeable battery provides 80 hours of use between charges.

For even greater flexibility, optional **Bluetooth** or **wifi** connectivity is available, allowing you to **wirelessly transmit patient data** from the scale to a device or database. An optional **micro USB port** is also available. Please contact us for more information about wireless data transfer.

Specification

Capacity: 250kg or 300kg

Graduations: 50g <150kg> 100g

Power: Mains adaptor or rechargeable battery

Dimensions: 530mm x 460mm Weight of Scale: Approx. 27kg

Optional Extras
Printer (TP-2100)
Bluetooth/Wifi/USB
Handle

Alternative footrests

Marsden M-200: Key features at-a-glance



Tel: 01709 364296 / 0800 169 2775

Fax: 01709 364293

E-mail: sales@marsdengroup.co.uk

www.marsden-weighing.co.uk

f MarsdenWeighing

▶ @MarsdenWeighing

You Tube www.youtube.com/MarsdenWeighing

http://pinterest.com/marsdenweighing

