

Plinth Mechanism

Electric Plinth

Hydraulic Plinth

Product Code

BA02E

BA02H

Two Section Plinth


Provides great value with strength, elegance and reliability. While offering a robust platform for a range of medical procedures, this couch provides ease of use, whilst maintaining patient comfort.

Standard Features:

- 260kg safe working load for electric mechanism
- 225kg safe working load for hydraulic mechanism
- Paper roll holder fitted as standard
- Retractable rubber castors & Adjustable foot for uneven floors
- MRSA resistant anti bacterial & fire retardant upholstery
- Gas spring assisted back section
- Extra low entry height for wheel chair transfer
- World leading Linak 'High Reliability' actuators
- Upholstery available in 14 standard colours
- Manufactured in the UK and CE Marked
- **Warranty:** 7 Years on Couch Frame & 3 Years on moving parts

Popular Options & Accessories:

- Breathing Hole
- Foot Switch (For Electric Versions)
- Adjustable Arm Rests
- Painted/Chrome Cotsides
- See full range of Options & Accessories on page 30/31



Plinth Mechanism

Electric Plinth

Hydraulic Plinth

Product Code

BA03E

BA03H

Three Section Plinth

Medi-Plinth's best selling 3 section couch enables advance patient positioning for a wide range of examination and treatment procedures providing optimum versatility.

Standard Features:

- 260kg safe working load for electric mechanism
- 225kg safe working load for hydraulic mechanism
- Paper roll holder fitted as standard
- Retractable rubber castors & Adjustable foot for uneven floors
- MRSA resistant anti bacterial & fire retardant upholstery
- Gas spring assisted back section
- Extra low entry height for wheel chair transfer
- World leading Linak 'High Reliability' actuators
- Upholstery available in 14 standard colours
- Manufactured in the UK and CE Marked
- **Warranty:** 7 Years on Couch Frame & 3 Years on moving parts

Popular Options & Accessories:

- Detachable Headrest
- Tilting Seat Option
- Battery Operating System
- Breathing Hole
- See full range of Options & Accessories on page 30/31