

NIO™ Pediatric

Quick, Safe, Easy Vascular Access

- Adjustable according to patient age
- Redundant safety features
- Semi-automatic device
- No exposed needle
- Single-use device
- 5-year shelf life
- Self-contained



Vascular Access in a Flash

When you need immediate vascular access, you can be confident that the NIO™ Pediatric is complete and ready to use. The NIO Pediatric is the quickest, most reliable intraosseous device for gaining immediate vascular access in pediatric patients.

NIO™ Pediatric

Product Highlights

Semi-automatic, self-contained unit facilitates rapid vascular access.

Unexposed needle decreases the risk of accidental needle injury.

Redundant safety features prevent accidental activation of the device.



Specifications

Dimensions	148.55mm (5.9in) L 64.5mm (2.6in) W 43.7mm (1.72in) H
Weight	85.6g
Packaging	Packed individually in a peel-back blister pack
Needle Gauge	18G
Needle Effective Length	38.1mm (1.5in)
Needle Penetration Depth	3-9 years: 14mm (Leave red spacer intact) 9-12 years: 18 mm (Remove red spacer)
Safety Mechanism	1. Rotating safety latch prevents accidental activation 2. Simultaneously compressing the device against the skin while squeezing the trigger wings activates device
Shelf Life*	5 years

Specifications are subject to change without notice.
*From date of manufacture

The NIO Pediatric features an innovative stabilizer hub which serves four functions:

- **Unique location arrows** aid insertion site identification
- **Needle stabilizer** assists in securing the needle to the limb.
- **Adjustable needle penetration depth** means one device can accommodate a wide range of pediatric patients
- **Mechanical stopper** limits needle depth, preventing over-penetration

PerSys Medical
5310 Elm St.
Houston, TX 77081 USA
1.888.737.7978
info@ps-med.com

www.persysmedical.com

