


 **Medisave**
www.medisave.co.uk

seca authorised dealer

go to www.medisave.co.uk for the full
seca range

seca 384

Digital baby scale,
also ideal as floor scale
for children 



- Baby scale and floor scale for children in one.
- Weighing tray easy to detach from base.
- Breast-Milk-Intake-Function.
- Mobile and robust for transport.

seca 384:

Convertible and robust –
compact quality with dual purpose.

At home, in the surgery and around the world wherever life begins: reliable weight control is decisive for a healthy start in life. The seca 384 is an all-rounder which has been developed on the basis of practical experience. It can be used not only as a baby scale but also as a floor scale for small children.



Baby scale and floor scale for children in one.

Babies can be weighed safely and comfortably in the tray of the seca 384. The rounded shape of the tray with its slightly raised sides gives the baby a feeling of security. By detaching the tray from the base, the scale is converted into a fully functional digital floor scale with a capacity of 20 kg.

Very safe due to secure attachment of tray.

In no time at all, the baby scale can be converted into a floor scale: pressing a button unfastens the mechanism attaching the tray to the base. When locked, this secure attachment ensures that the baby scale is easy to lift and safe to transport.

Easy to operate due to self-explanatory button functions.

The seca 384 is easy to operate due to its intuitive operating concept with self-explanatory buttons. The HOLD function means that, after being weighed, the baby can first be taken care of before one has to read off the result. With the TARE function, the weight of nappies or soft mats can be reliably deducted. With the Breast-Milk-Intake-Function, one can determine how much milk the baby has taken when breast fed. Because, at the touch of a button, the weight of the baby determined before feeding is deducted from the weight determined after feeding.

Mobile and robust for transport.

Its compact, robust construction predestined the seca 384 as the scale for all those who have to be mobile in the interest of the health of their youngest patients. The power-saving battery operation – up to 20,000 weighings with just one set of batteries – makes the scale independent of electricity power supply.



seca 413



The weighing tray can be removed quickly and easily at the touch of a button.



seca 384

Technical Data

- Capacity: 20 kg
- Graduation: 10 g < 10 kg > 20 g
- Dimensions (WxHxD):
552 x 162 x 362 mm
- Weight: 2.7 kg
- Power supply: Batteries
- Functions: automatic switch-off, TARE, HOLD, BMIF, weighing range switch-over
- Optional: Carry case seca 413
- Approval class: