

°Control
TRUE TEMPERATURE



Lec
MEDICAL

Experts in
Medical Refrigeration

273L Freestanding Pharmacy Fridges



Key Features

- ✓ Represents your vaccine vial
- ✓ Provides you with the true temperature of your vaccines
- ✓ Ensures you always measure what matters
- ✓ Built-in warning system for safeguarding vaccines
- ✓ Now with 5 year warranty as standard *UK only

5 YEARS
PARTS AND LABOUR
WARRANTY*

*UK Only



Stoney Lane, Prescot, Merseyside, L35 2XW

Tel: 0844 815 3742
Email: sales@gdpa.co.uk

Control
TRUE TEMPERATURE



PGRC273UK Pharmacy fridge – Glass door

SKU Code	444443778
EAN Code	5052263037786
Gross Capacity (Litres) / Net Capacity (Litres)	273 / 269
Unit External Dimensions W / D / H (mm)	595 / 625 (660 with handle) / 1565
Product Weight (kg)	68.5
Internal Dimensions W / D / H (mm)	480 / 445 / 1382 (step:175D x 255H)
Temperature Range	+2°C to +8°C
Secondary Probe	Yes
Bottle Contents / Size	Silicone / 50ml
External Finish	Painted Aluminium Door
Internal Finish	HIPS*
Type of Plastic	SRL800 (HIPS/ABS Mix)
Door Finish / Type	Painted Aluminium / Glass
Door Style	Hinged, Right Handed only, Non-reversible
Door Open Depth (mm)	1180
Replacement Door Seal	Yes
Open Door Sensor	Yes
Lock	Yes
Adjustable Feet (number & position)	2 Adjustable Front Feet, 2 Rear Rollers
Temperature Display on Panel / Controller Type	Yes / Digital, External
Temperature Alarm	Yes
Number of Shelves / Shelf Positions / Shelf Material	5 / 12 / Wire
Shelf Size (mm) (WxDxH)	4: 385 x 500 x 6, 1: 340 x 497 x 6
Number of Baskets / Basket Material	1 / Wire
Interior Light / Fan / Spark free	Yes / Yes / No
Auto or Manual Defrost	Automatic
Climate Class	SN
Refrigerant Type	R600a
CFC Free / HFC Free / HCFC Free / Ammonia Free	Yes
UK 3 Pin Plug / Other Plug Types Available	Yes / Yes
Voltage Supply	230
Aux Output	Yes
Insulation Thickness (mm) / Insulation Type	45 / Cyclopentane
Anti Microbial Handle	Yes
Optional Drawer System Available	Yes