










LARGE CAPACITY FREESTANDING PHARMACY FRIDGE SOLID DOOR 600L WHITE



www.facebook.com/Lec-Medical
www.twitter.com/Lec_Medical

- 
+2°C to +8°C Operating Temperature Range
 Conforms to MHRA guidelines.
- 
NEW Advanced Controller
 Records temperatures accurate to 0.1°C.
- 
USB Data Logging
 Allows user to accurately monitor temperature.
- 
Audible Warning Alarms
 Alerts user should internal temperature deviate from set parameters or in the event of a power failure.
- 
Fitted Lock as Standard
 Maintains security and ensures the unit can only be accessed by approved personnel.
- 
Fitted Castors (Lockable)
 For ease of movement.
- 
Fitted with 5 Adjustable Wire Shelves
 Enables optimum air flow inside the unit.

SKU CODE:	PRODUCT CODE:	DESCRIPTION:
444440367	PSR600UK	Large Capacity Pharmacy Fridge Solid Door 600L White



FULL SPECIFICATION – PSR600UK LARGE CAPACITY FREESTANDING PHARMACY FRIDGE SOLID DOOR 600L WHITE



Description	Freestanding
SKU code	444440367
GDPA code	PSR600UK
Gross capacity (litres)	600
Unit external width (mm)	690
Unit external depth (mm)	840
Unit external height (mm)	1995
Product weight (kg)	135
Packed width (mm)	840
Packed depth (mm)	990
Packed height (mm)	2195
Packed weight (kg)	150
Country of origin	UK
Internal dimensions (mm) (WxDxH)	560 x 680x 1230
Temperature range	+2°C to +8°C
Colour	White/green
External finish	White
Internal finish	Stainless steel
Type of plastic	SRL800 (HIPS/ABS Mix)
Door finish	Painted stainless steel
Door type	Solid
Door style	Hinged
Door open depth	1145
Lock	Yes
Adjustable feet (number & position)	Fitted castors
Temperature display on panel	Yes
Air Temperature alarm	Yes (Audio & Visual)
Mains failure alarm	Yes (Audio & Visual)
Load alarm delay (with no delay)	Yes (Audio & Visual)
Door open alarm	Yes (Audio & Visual)
Compressor alarm	Yes (Visual)
Air/Load/Defrost sensor fail alarm	Yes (Visual)
Battery low alarm	Yes (Visual)
Number of shelves	5
Number of shelf positions	20
Shelf size (WxD)	530 x 650
Shelf material	Plastic coated wire
Interior light	No
Fan	Yes
Spark free	No
Auto or manual defrost	Automatic
Controller type	Digital, external
Refrigerant type	R134a
Climate Class	ST
CFC free	Yes
HFC free	Yes
HCFC free	Yes
Ammonia free	Yes
UK 3 pin plug	Yes
Other plug types available	Yes
Voltage supply	230
Aux output	Yes
Insulation thickness (mm)	60
Insulation type	Cyclopentane
Optional drawer system available	No