




Experts in  
Medical Refrigeration

# 151L Under-Counter Pharmacy Fridge



## Key Features

- ✓ Represents your vaccine vial
- ✓ Provides you with the true temperature of your vaccines
- ✓ Ensures you always measure what matters
- ✓ Built-in warning system for safeguarding vaccines
- ✓ Now with 5 year warranty as standard \*UK only



\*UK Only



**Control**  
TRUE TEMPERATURE



## PSRC151UK Pharmacy Fridge – Solid Door

|  |   |
|--|---|
| SKU Code   | 444443767                               |
| EAN Code   | 5052263037670                           |
| Gross Capacity (Litres) / Net Capacity (Litres)      | 151 / 141                               |
| Unit External Dimensions W / D / H (mm)              | 595 / 600, 635 with handle / 845        |
| Product Weight (kg)                                  | 37                                      |
| Internal Dimensions W / D / H (mm)                   | 500 / 450 / 705 (step:165D x 185H)      |
| Temperature Range                                    | +2°C to +8°C                            |
| Secondary Probe                                      | Yes                                     |
| Bottle Contents / Size                               | Silicone / 50ml                         |
| External Finish                                      | Painted Steel                           |
| Internal Finish                                      | HIPS*                                   |
| Type of Plastic                                      | SRL800 (HIPS/ABS Mix)                   |
| Door Finish / Type                                   | Painted Steel / Solid                   |
| Door Style   | Hinged, Reversible                      |
| Door Open Depth (mm)                                 | 1165                                    |
| Replacement Door Seal                                | No                                      |
| Open Door Sensor                                     | Yes                                     |
| Lock   | Yes                                     |
| Adjustable Feet (number & position)                  | 2 Adjustable Front Feet, 2 Rear Rollers |
| Temperature Display on Panel / Controller Type       | Yes / Digital, External                 |
| Temperature Alarm                                    | Yes                                     |
| Number of Shelves / Shelf Positions / Shelf Material | 3 / 6 / Wire                            |
| Shelf Size W / D / H (mm)                            | 400 / 520 / 6                           |
| Number of Baskets / Basket Material                  | 1 / Wire                                |
| Interior Light / Fan / Spark free                    | No / Yes / No                           |
| Auto or Manual Defrost                               | Automatic                               |
| Climate Class  | SN                                      |
| Refrigerant Type                                     | R600a                                   |
| CFC Free / HFC Free / HCFC Free / Ammonia Free       | Yes                                     |
| UK Plug / Other Plug Types Available                 | Yes / Yes                               |
| Voltage Supply                                       | 230                                     |
| Aux Output   | Yes                                     |
| Insulation Thickness (mm) / Insulation Type          | 45 / Cyclopentane                       |
| Anti Microbial Handle                                | Yes                                     |
| Optional Drawer System Available                     | Yes                                     |