

Product Specification



Version: 03

Date of issue: 22/02/2017

Date of revision: 28/05/2020

Medipal Disinfectant Wipes

Medipal Disinfectant Wipes are intended for cleaning and disinfecting surfaces of non-invasive medical devices within healthcare and care environments. For professional use only.

Medipal Disinfectant wipes are proven to kill most microorganisms such as bacteria, yeasts and viruses. The wipes have been independently tested to a range of European testing standards (EN) which includes EN16615.



Product Description

Product Name	Medipal Disinfectant Wipes (50 wipes)
Product Code	S524110MPCE
Languages	EN, FR, AR

Pack Type & Dimensions

Count	50
Dimensions	200 x 175mm
Weight	0.221kg
Colour of wipes	White
Wipes material description	30gsm Tritex Technology

Solution

Solution Description	Medipal Disinfectant Solution
Active Ingredients Content	Didecyldimethylammonium chloride, 6.4g/kg N-(3-Aminopropyl)-N-dodecylpropane-1, 3-diamine, 5.31g/kg
Hazard identification	H315: Causes skin irritation H319: Causes serious eye irritation H412: Harmful to aquatic life with long lasting effects For further information please refer to the material safety data sheet

Primary Packaging

Pack Type	50 wipes PET/PE laminate film
Lid Description	Medipal Disinfectant Wipes Lens
Label	Not applicable
Barcode	5025254 012734

Product Specification



Version: 03

Date of issue: 22/02/2017

Date of revision: 28/05/2020

Medipal Disinfectant Wipes

Secondary Packaging

Case Description	Plain outer case
Barcode	5025254 042151

Packaging & Logistical Information

Product Name	50 Medipal Disinfectant Wipes
Product Code	S524110MPCE
Case Quantity	24 Packs/case
Case Dimensions [OD]	380 x 290 x 200
Case Weight	5.7kg
Case/Layer	8
Layers/Pallet	5
Cases/Pallet	40
Pallet Weight	243kg
Hazardous for Transport	Not applicable
Commodity code for export	3808941000

