

Technical Information



Version: 03

Date of revision: 30.03.2021

Medipal Alcohol Surface Wipe with 2% chlorhexidine (single sachet)

Description:

Medipal Alcohol Surface wipes with 2% Chlorhexidine are recommended for use on non invasive medical devices such as cannulas, needle free devices, hubs and ports and blood culture bottle tops. This is a Class IIa medical device and regulated under the MDD 93/42/EEC.



Microorganisms

Testing Standard

Contact time

Bactericidal

Pseudomonas aeruginosa
Staphylococcus aureus
Enterococcus hirae
Pseudomonas aeruginosa
Staphylococcus aureus
Enterococcus hirae

EN 16615
EN 16615
EN 16615
EN 14561
EN 14561
EN 14561

2 min
2 min
2 min
30 sec
30 sec
30 sec