

Spencer 406

Transport chair with 2 wheels



Spencer 406 is a transport chair for the evacuation and transport of patient from the rescue place to the ambulance or to another device intended for transport on the vehicle.

Specific features

Easy to clean PVC sheets
Robust and lightweight structure
Large rear wheels to reduce vibration during patient transfer
Quick and intuitive closure system
Safety system against accidental closure

Reference standard

NA

Standard equipment

ST00423A	Copule of belts for chairs
----------	----------------------------

Optional accessories

ST42704E	Wall fastener 10g certified for chair 406
----------	---

Technical data

Length	840 ± 10 mm
Width	410 ± 10 mm
Height opened	940 ± 10 mm
Height folded	770 ± 10 mm
Thickness folded	235 ± 10 mm
Rear wheels	∅ 200 mm
Materials	Steel/Al/PVC
Weight	10 kg ± 0,3
Loading capacity	150 kg

Class I Med.Device according to Dir. 93/42/CEE

406 SPENCER AMBULANCE
CARRYING CHAIR
ST00406A



CND Classification	V9099
Registration number	93547

Rev.1 (18/07/2017)

COPIA NON CONTROLLATA – successive revisioni saranno disponibili al link <http://support.spencer.it>

Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

spencer.it