

MARSDEN

T-210 Automatic Wall Mounted Thermometer



- MDD Approved
- Wall mounted, non-contact thermometer
- Simple and easy to use
- Measures temperature in °C or °F
- Instant audio and visual warning
- 1 year warranty

Marsden T-210:

Automatic Wall Mounted Thermometer

- Accurate to 0.1°C
- MDD Approved
- Affix to a wall or optional stand
- Clear and easy to read digital reading
- Instant visual and audible warning
- Non contact device - responds to motion
- Measures temperature in °C or °F
- Auto-power off to prolong battery life
- Replaceable batteries to aid portability
- Detects signs of virus in seconds
- An alert will signal when higher than 37.8 °C
- Powered by USB Type C lead



What is the T-210?

The T-20 is an MDD Approved wall mounted digital infrared thermometer that makes taking a persons temperature incredibly simple and easy to use.

The device can be mounted directly onto a wall, or attached to our optional stand. An individual simply then walks up to the T-210 thermometer and their temperature will clearly show on the display

What are the benefits?

An automatic thermometer allows you to easily and instantly check the temperature of employees, customers or patients at entrance points.

The T-210 is MDD Approved, meaning dependable accuracy when measuring temperature - peace of mind for you. MDD Approved thermometers can also be used in a medical setting.

How accurate is the T-210?

The MDD Approved T-210 thermometer provides a highly accurate temperature reading within 0.2°C. All readings will appear clearly on the display.

The T-210 has a rechargeable battery fitted as standard and will last for up to 300 measurements. To charge the device simply use the 3m USB-C cable. The device can still be used whilst on charge.

Method:	Wall Mounted or Stand	Power Supply:	Replaceable Batteries
Graduations:	0.2°C	Battery Life:	6 hours of continuous use
Functions:	Unit Switchable, Portable	Approvals:	MDD Approved
Dimensions:	146mm x 37mm x 39mm	Optional Extras:	Stand
Weight of Scale:	1kg	Warranty Period:	1 Year