

FCP-02, 6" Emergency Bandage®

aka The Israeli Bandage®

- Best-in-class multiple use bandage addressing various wounds and wraps
- Coverage for small, medium, and large wounds
- Wrap severe wounds and small area amputations
- Can be used as an improvised tourniquet
- 6" non-adherent pad
- 8-year shelf life



Full Coverage and Multiple Functions in One Simple Package

The FCP-02 Emergency Bandage® provides everything needed to address small, medium, and large wounds in one compact package.

The versatile design is an effective solution for the control of severe bleeding, wound coverage, splinting, slings, head wraps, and creating an improvised tourniquet. The device is intuitive and simple to use for trained and untrained responders.



6" Green Emergency Bandage®

Product Highlights

- Trusted original Emergency Bandage® design
- Multiple use bandage in compact packaging
- Intuitive design facilitates easy, rapid application




This compact, all-in-one bandage is an ideal component for any first aid kit

- Bleeding control, compression, and coverage for a variety of wounds
- Easy, intuitive application
- Compact packaging reduces weight and volume
- Integrated pressure bar applies direct pressure at wound
- Single sterile pad
- Elasticized wrapping leader for compression
- Stop-N-Go™ release for controlled wrapping
- Simple closure bar
- Sterile packaging can be used as an improvised chest seal
- Sleeve for self-application

Specifications

| | |
|--------------------------|---|
| Description | 6" (15cm) wide bandage |
| Part Number | FCP-02 |
| NSN | 6510-01-492-2275 |
| Pad Size | 6" wide x 9" long (15cm x 23cm) |
| Full Length | 65" (165cm) unstretched, including pad |
| Shelf Life | 8 years |
| Colors | OD green elastic leader, white pad |



PerSys Medical
5310 Elm St.
Houston, TX 77081 USA
1.888.737.7978
info@ps-med.com
www.ps-med.com

