



ULTRASPEC®

SINGLE-USE EAR-IRRIGATION TANK

▶ WITH NEW ANTI-SPLASH FEATURE!

Stable

Solid base for improved stability compared to pulp



Single-use

Cost-effective. Reduced risk of cross-contamination and infection



Anti-splash tab

Improved patient experience



Green

Lighter than pulp tank and 50% less bulky to store/transport



Leak-resistant

Smooth, coated medical card forms a clean seal against patients' skin



Space-saving

Thin-wall construction takes up 50% less storage space



New, Improved Noots-Type Ear Tank

Facilitates easier, cleaner treatment of cerumen. The anti-splash tab significantly reduces back-wash during ear-irrigation procedures.

UltraSpec®

products for the people who care™

Instructions for use



1 Simply fold anti-splash tab into the cup

2 With tab towards patient, place the Ear-irrigation tank under the ear and press gently against the skin

3 Dispose of in accordance with Local Health Authority guidelines following single use.

Order Information

Description: Ultraspec Ear-Irrigation Tank

Inner pack: Sleeve of 50

Outer pack: 20 Sleeves

Order code: 3USNWHB



products for the people who care™

Ultraspec - A Splice Cast Ltd Brand

Email: ultraspec@splicecast.co.uk
Web: www.ultraspecmedical.co.uk

© Splice Cast Ltd 2017 – Ver 1-0617

AVAILABLE FROM:

Call Splice Cast Ltd for your nearest distributor

Tell: +44 (0)1686 625 750

Call for a free sample today!