

Product Specification



Version: 03

Date of issue: 22/02/2017

Date of revision: 28/05/2020

Medipal Disinfectant Wipes

Medipal Disinfectant Wipes are intended for cleaning and disinfecting surfaces of non-invasive medical devices within healthcare and care environments. For professional use only.

Medipal Disinfectant wipes are proven to kill most microorganisms such as bacteria, yeasts and viruses. The wipes have been independently tested to a range of European testing standards (EN) which includes EN16615.



Product Description

Product Name	Medipal Disinfectant Wipes (240 wipes)
Product Code	W520110MPCE
Languages	EN, FR, AR

Pack Type & Dimensions

Count	240
Dimensions	195 x 165mm
Weight	0.710kg
Colour of wipes	White
Wipes material description	20gsm Tritex Technology

Solution

Solution Description	Medipal Disinfectant Solution
Active Ingredients Content	Didecyltrimethylammonium chloride, 6.4g/kg N-(3-Aminopropyl)-N-dodecylpropane-1, 3-diamine, 5.31g/kg
Hazard identification	H315: Causes skin irritation H319: Causes serious eye irritation H412: Harmful to aquatic life with long lasting effects For further information please refer to the material safety data sheet

Primary Packaging

Pack Type	2 Litre HDPE canister
Lid Description	Green PP dispenser lid
Label	240 Medipal Disinfectant Wipes label
Barcode	5025254 008553

Product Specification

Version: 03

Date of issue: 22/02/2017

Date of revision: 28/05/2020

Medipal Disinfectant Wipes

Secondary Packaging

Case Description	Plain outer case
Barcode	5025254 042205

Packaging & Logistical Information

Product Name	240 Medipal Disinfectant Wipes
Product Code	W520110MPCE
Case Quantity	10 canisters/case
Case Dimensions [OD]	560 x 225 x 245
Case Weight	7.5kg
Case/Layer	9
Layers/Pallet	5
Cases/Pallet	45
Pallet Weight	420kg
Hazardous for Transport	Not applicable
Commodity code for export	3808941000

