

# Product Specification



Version: 03

Date of issue: 22/02/2017

Date of revision: 28/05/2020

## Medipal Disinfectant Wipes

Medipal Disinfectant Wipes are intended for cleaning and disinfecting surfaces of non-invasive medical devices within healthcare and care environments. For professional use only.

Medipal Disinfectant wipes are proven to kill most microorganisms such as bacteria, yeasts and viruses. The wipes have been independently tested to a range of European testing standards (EN) which includes EN16615.



### Product Description

|              |  |
|--------------|--|
| Product Name | Medipal Disinfectant Wipes (350 wipes) |
| Product Code | W522110MPCE                            |
| Languages    | EN, FR, AR                             |

### Pack Type & Dimensions

|                            |                         |
|----------------------------|-------------------------|
| Count                      | 350                     |
| Dimensions                 | 250 x 225mm             |
| Weight                     | 1.740kg                 |
| Colour of wipes            | White                   |
| Wipes material description | 20gsm Tritex Technology |

### Solution

|                            |  |
|----------------------------|--|
| Solution Description       | Medipal 3in1 Disinfectant Solution   |
| Active Ingredients Content | Didecyltrimethylammonium chloride, 6.4g/kg<br>N-(3-Aminopropyl)-N-dodecylpropane-1, 3-diamine, 5.31g/kg  |
| Hazard identification      | H315: Causes skin irritation<br>H319: Causes serious eye irritation<br>H412: Harmful to aquatic life with long lasting effects<br>For further information please refer to the material safety data sheet |

### Primary Packaging

|                 |                                      |
|-----------------|--------------------------------------|
| Pack Type       | 3 Litre PP bucket                    |
| Lid Description | Green PP dispenser lid               |
| Label           | 350 Medipal Disinfectant Wipes label |
| Barcode         | -                                    |

# Product Specification



Version: 03

Date of issue: 22/02/2017

Date of revision: 28/05/2020

## Medipal Disinfectant Wipes

### Secondary Packaging

Case Description Plain outer case  
Barcode

### Packaging & Logistical Information

|                           |                                |
|---------------------------|--------------------------------|
| Product Name              | 350 Medipal Disinfectant Wipes |
| Product Code              | W522110MPCE                    |
| Case Quantity             | 4 buckets/case                 |
| Case Dimensions [OD]      | 408 x 223 x 348                |
| Case Weight               | 7.5kg                          |
| Case/Layer                | 10                             |
| Layers/Pallet             | 4                              |
| Cases/Pallet              | 40                             |
| Pallet Weight             | 315kg                          |
| Hazardous for Transport   | Not applicable                 |
| Commodity code for export | 3808941000                     |

