



HEAD OFFICE

Owen Mumford Ltd
Brook Hill, Woodstock
Oxfordshire OX20 1TU, United Kingdom
T : +44 (0)1993 812021
E : info@owenmumford.co.uk

UK / IRELAND SALES & MARKETING

Owen Mumford Ltd
Primsdown Industrial Estate
Worcester Road, Chipping Norton
Oxfordshire OX7 5XP, United Kingdom
T : +44 (0)1608 645555
E : info@owenmumford.co.uk

NORTH AMERICA

T : +1 770-977-2226
E : info@owenmumfordinc.com

FRANCE

T : +33 (0)2 32 51 88 70
E : fr.info@owenmumford.com

GERMANY

T : +49 (0)6026/97750
E : de.info@owenmumford.com

CHINA

T : +86 (0)21 6032 3528
E : cn.info@owenmumford.com

MALAYSIA

T : +60 (0)7 553 9617
E : my.info@owenmumford.com

MEXICO

T : +52 (55) 3733 8759
E : mx.info@owenmumford.com

owenmumford.com

OUR ACCREDITATIONS INCLUDE:

- ISO13485 Certification
- Medical Device Directive 93/42/EEC
- CMDCAS 13485 (Canadian Health Authority)
- Environmental Management Systems ISO14001



© Copyright Owen Mumford 2016
AHJUP-GLO-LT/MCC/0116/02

Premium quality designed
for patient comfort
That's the point.



Unifine® Pentips®

Premium quality pen needles



More reason to smile

owenmumford.com

Unifine® Pentips®

Premium quality pen needles

Unifine® Pentips® have been developed for a comfortable injection experience. Three bevelled, thin gauge needles combined with anti-coring treatment and high quality silicone, designed to reduce penetration force for a comfortable injection.

Compatible with all major brands of injection pen.¹

Designed for compatibility so that the benefits of Unifine® Pentips® can be enjoyed by all.

Designed to reduce penetration force.

Unifine® Pentips® have three bevelled needles, designed to reduce penetration force.

Clean piercing of insulin cartridge.

Anti-coring technology designed to ensure a clean piercing of the cartridge.

Product range

4mm	5mm	6mm	8mm	12mm
0.23mm (32G) x 4mm	0.25mm (31G) x 5mm	0.25mm (31G) x 6mm	0.25mm (31G) x 8mm	0.33mm (29G) x 12mm
Box of 100	Box of 100	Box of 100	Box of 100	Box of 100
OM code AN3540	OM code AN3550	OM code AN3590	OM code AN3530	OM code AN3529

Compatible with all major brands of injection pens.²



Including:



Owen Mumford: **Autopen® Classic**



Owen Mumford: **Autopen® 24**



Including:



Novo Nordisk: **NovoPen® 4**



Novo Nordisk: **NovoPen Echo®†**



Novo Nordisk: **NovoPen® FlexTouch†**



Novo Nordisk: **Victoza®**



Including:



Lilly: **Byetta®**



Lilly: **Humalog® KwikPen™**



Lilly: **Humapen® Luxura**



Lilly: **Humapen® Savvio™ (All Colours)**



Including:



Sanofi: **Apidra® Solostar®**



Sanofi: **ClikSTAR®**



Sanofi: **Lantus® Solostar®**

¹. Compatibility data correct as of 01/09/14.

Unifine® Pentips® and Autopen® are registered trademarks of Owen Mumford Ltd. All rights reserved. All other trademarks are the property of their respective owners.