

COMMUNITY SIDE RAILS

SIDHIL

SOLITE SAFE-SIDE MESH SIDE RAIL

SOLITE/SR



The Solite Safe-side full length mesh side rails have been designed specifically for the Sidhil Solite bed range. They are padded along the top rail to protect the user. The mesh panels aid visibility and help air circulation around the user. The side panels are easily removed and washable. The sides can be raised and lowered by the carer with minimal disturbance to the patient. The Safe-Sides come in an easy to assemble flat pack form, with no tools needing to be supplied by the customer.

SPECIFICATIONS		FEATURES AND BENEFITS	
Length	1985mm	<ul style="list-style-type: none"> • Easy assembly and operation • Breathable mesh side panels • Mesh aids visibility for both patient and carer • Washable and waterproof • Folds away to below level of mattress • Reduced risk of entrapment 	
Weight	16kg (boxed pair)		
Height	(From top of mattress platform) max - 500mm min - 140mm		
Bed Suitability	SOLITE LASER /4, 2 AND VH SOLITE 4, 2 AND VH		

www.sidhil.com