



Accessories

Master Key LABLOCK3KEYS

This is a master key for use with Labcold DIGLOCK IntelliCold® pharmacy refrigerators, both solid and glass door models.

Exclusively from Labcold, this key serves as an ideal override solution, providing access to the refrigerator when the code is forgotten. It allows for quick removal of contents or repurposing of the unit, especially useful during power failures if backup batteries are not installed, providing nothing but peace of mind.

Specifications

Specifications	
Product Code	LABLOCK3KEYS

Features

- Suitable for any RLDF/G DIGLOCK units
- Comes as a pair of two
- Override keys

