

PRODUCT INFORMATION

Item Number	SNA
Product Colour	Cream
Size Range	5.5-9
Fingertip Thickness	0.20 mm
Palm Thickness	0.18 mm
Typical Weight	14.00 g
Shelf Life	3 Years

FEATURES

Sterile: Sterilised by irradiation method

Sensitivity: Accelerator free synthetic gloves removing the risk of contact dermatitis caused by these chemicals

Grip: Micro Roughened to provide excellent feel and instrument grip in both wet and dry conditions

Comfort: Manufactured on ergonomically designed formers from an advanced synthetic formulation, offers ultimate comfort and tactile sensitivity

Fit: Hand specific to reduce hand fatigue

Sensitivity: Latex free, reducing the risk of hypersensitivity and allergic reactions

Material Composition: Manufactured from Accelerator Free Polychloroprene.

SECTORS

Medical, Vet and Dental


STANDARDS

AQL 0.65

This product is classified as Class IIa Medical Device in compliance with the requirements of Medical Directive 93/42/EEC on Medical Devices (MDD)

Manufactured in accordance with Personal Protective Equipment (PPE) by the European PPE Regulation (EU) 2016/425

Compliant with European Standards EN455: Medical Gloves for Single Use parts 1,2,3&4

Tested to European Standards EN420 & EN ISO 374-1 & EN ISO 374-5



STERILE R



KMNPST VIRUS



For further technical information, please contact Customer Services on 0333 320 8550 or the Technical team on 0333 320 8555

Neotex PureDerma Polychloroprene Accelerator-Free Surgical Gloves

STORAGE INSTRUCTIONS

Store in dry conditions, away from direct sunlight and extreme heat. It is recommended to retain the packaging for reference. Always inspect gloves prior to use and do not use if damaged or showing signs of wear and tear.

DISPOSAL INSTRUCTIONS

Dispose of the product responsibly in accordance with local regulations.

PACKING

1 pair per pouch
50 pairs per dispenser
4 dispensers per case
30 cases per pallet

INNER PACKAGING



Size	Inner Packaging Barcode
5.5	5060040798621
6	5060040798638
6.5	5060040798645
7	5060040798652
7.5	5060040798669
8	5060040798676
8.5	5060040798683
9	5060040798690

Inner Packaging	
Inner dimensions (LxWxH)	290 x 153 x 225mm

Size	Outer Packaging Barcode
5.5	5060040798539
6	5060040798546
6.5	5060040798553
7	5060040798560
7.5	5060040798577
8	5060040798584
8.5	5060040798591
9	5060040798614

Outer Packaging	
Outer Dimensions (LxWxH)	595 x 315 x 235mm

For further technical information, please contact Customer Services on 0333 320 8550 or the Technical team on 0333 320 8555

**Neotex PureDerma Polychloroprene Accelerator-Free
Surgical Gloves****PLEASE NOTE:**

Weight does not form part of the specification and is subject to change. Palm thickness is approximate and based on an agreed manufacturing tolerance. All thicknesses and weights given are an average for medium (8) gloves. All thicknesses and weights given are correct to the best of Polyco Healthline Ltd.'s knowledge at the time of input, however due to the manufacturing process, both thickness and weight may vary between batches.

For further technical information, please contact Customer Services on 0333 320 8550 or the Technical team on 0333 320 8555