



HEAD OFFICE

Owen Mumford Ltd
Brook Hill, Woodstock
Oxfordshire OX20 1TU, United Kingdom
T : +44 (0)1993 812021
E : info@owenmumford.co.uk

UK / IRELAND SALES & MARKETING

Owen Mumford Ltd
Primsdown Industrial Estate
Worcester Road, Chipping Norton
Oxfordshire OX7 5XP, United Kingdom
T : +44 (0)1608 645555
E : info@owenmumford.co.uk

NORTH AMERICA

T : +1 770-977-2226
E : info@owenmumfordinc.com

FRANCE

T : +33 (0)2 32 51 88 70
E : fr.info@owenmumford.com

GERMANY

T : +49 (0)6026/97750
E : de.info@owenmumford.com

CHINA

T : +86 (0)21 6032 3528
E : cn.info@owenmumford.com

MALAYSIA

T : +60 (0)7 553 9617
E : my.info@owenmumford.com

MEXICO

T : +52 (55) 3733 8759
E : mx.info@owenmumford.com

owenmumford.com

OUR ACCREDITATIONS INCLUDE:

- ISO13485 Certification
- Medical Device Directive 93/42/EEC
- CMDCAS 13485 (Canadian Health Authority)
- Environmental Management Systems ISO14001



© Copyright Owen Mumford 2016
AHJAN/4200/LTMCC/0116/02

Simplifying insulin injections



Autopen[®]

Insulin delivery pens



More reason to smile

owenmumford.com

Autopen®

Insulin delivery pens


Every day, millions of people inject themselves with insulin using an injection pen.¹ **Autopen®** has been designed to make this experience as easy as possible with a side-mounted release button.

Compatible² with all major brands of pen needle.

Autopen® has been designed to be compatible with all major brands of insulin pen needles including **Unifine® Pentips®** and **Unifine® Pentips® Plus**.

Simple, automatic delivery using a side-mounted release button.

Easy to use dose selector with audible clicks.



1. Kellie J. Antinori-Lent. Analysis of Comparison of Patient Preference for Two Insulin Injection Pen Devices in Relation to Patient Dexterity Skills. Journal of Diabetes Science and Technology. Volume 6, Issue 4, (July 2012). 2. Universal ISO 11608 thread on Autopen designed for compatibility with all major brands of pen needles including Unifine. 3. ACS22/V10/1/TF/01 (1998)

Accessories included with Autopen®

Autopen® accessories create an even easier injection experience for people with limited hand movement or grip.



Release Button Extension

Increases area to make it easy to push the release button.



Dose Selector Adaptor

Soft grip and larger area to make dialling a dose easy.

Pens to suit every user.

Autopen® 24

Compatible with **Sanofi** insulin.



AN4210 1-21 units in increments of 1 unit



AN4200 2-42 units in increments of 2 unit

Autopen® Classic

Compatible with **Lilly**³ and **Wockhardt** insulin.



AN3810 1-21 units in increments of 1 unit



AN3800 2-42 units in increments of 2 unit