

TECHNICAL INFORMATION SHEET

BD Vacutainer® Barricor™ Plasma Blood Collection Tube with BD Hemogard™ Safety Closure



Product Catalogue Number: **365055**

Product Description

Single use, evacuated, sterile blood collection tubes containing an inert mechanical separator and Lithium Heparin intended for the primary containment and preservation of specimens for the purposes of in-vitro diagnostic examination. Used to obtain and separate a plasma sample. These products are intended for use by healthcare professionals.

Manufacturing Information

Manufacturer:	Becton, Dickinson and Company, 1 Becton Drive, Franklin Lakes, NJ 07417-1885, USA
Standards & Certificate Numbers:	EN ISO 13485:2012 - MD19.2137
Country of origin:	UK
Certification body:	NSAI
Notified Body:	NSAI (0050)
EU Authorised Representative:	BD Switzerland Sarl, Terre Bonne Park - A4, Route de Crassier 17, 1262 Eysins, Switzerland

Sterilisation

Method:	Gamma Irradiation, Co-60
SAL:	10 ⁻⁶
Standards applied:	EN ISO 11137-1:2006

Product Standards & Guidelines

Standards:	ISO 6710:1995, EN 14820:2004
Guidelines:	Clinical and Laboratory Standards Institute (CLSI; Formerly NCCLS): Tubes and Additives for Venous Blood Specimen Collection; Approved Guideline (6th Edition). Document GP39-A6. Wayne, PA, USA, 2010.

Compliance

Directive:	European In Vitro Diagnostic Medical Devices Directive 98/79/EC
Classification:	Non Annex II / General IVD

Product Specification

Tube material:	Polyethylene Terephthalate (PET)
Tube size (mm):	13 X 75
Draw volume (mL):	3
Additives:	Lithium Heparin (17 IU/mL)
Separator:	Mechanical Separator: Thermoplastic Elastomer (TPE) & High Density Polypropylene
Closure material (cap):	Low Density Polyethylene (LDPE)
Closure material (stopper):	Bromobutyl Elastomer
Closure colour:	Lime Green
Product Storage:	Do not expose to direct sunlight Store product between 4° and 25°C
Label type:	Block
Shelf-life:	16 months
Global medical device nomenclature (GMDN):	61902
Material Safety Data Sheet (MSDS):	Please follow directions on http://regdocs.bd.com/
Fill line indicator:	Yes



Materials

Latex (NRL):	No
Dry Natural Rubber (DNR):	No
Phthalates:	No
Material of animal origin:	Lithium Heparin: Porcine Origin.
Other Materials of Concern:	Low levels of 2-(2H-benzotriazol-2-yl)-4,6-ditertpentylphenol (UV-328), CASN 25973-55-1, Should you have any further questions about BD's REACH process, please contact: BD_Product_Stewardship@bd.com

Packaging Specifications

100 unit pack weight (kg):	0.64	100 unit packaging material:	Expanded Polystyrene (EPS) / Polyolefin film
100 unit pack volume (m ³):	0.00227	100 unit packaging weight (kg):	0.02
100 unit pack dimensions LxHxW (mm):	178 x 85 x 150	100 unit packaging volume (m ³):	0.00068
1000 unit pack weight (kg):	6.71	1000 unit packaging material:	Cardboard
1000 unit pack volume (m ³):	0.025337	1000 unit packaging weight (kg):	0.635
1000 unit pack dimensions LxHxW (mm):	438 x 311 x 186	1000 unit packaging volume (m ³):	0.025337

