



ST10501


4BELL STAIR+

Foldable patient transfer chair with slides

The combination of reinforced but lightweight frame, telescopic front handles and the addition of slides makes patient transport smooth and safe for the operator. The chair's fastener is 10G certified.



Load Capacity



Stair Tracks



Foldable



Lightweight

UNI EN 1865-4
EN 1789



4BELL STAIR+ 10G FASTENER

It is installed in ambulances, special vehicles or airplanes to secure all 4Bell+ chairs in foldable position. The device can be fixed both on wall or floor, according to the vehicle's structure and configuration.

OPTIONAL ACCESSORIES AVAILABLE:

- Armrests **ST10450**
- Rear auxiliary handles **ST10451**
- XL Seat **ST10452**
- Headrest **ST10454**
- 4BELL STAIR+ 10G FASTENER **ST10455**



- **Telescopic front handles** facilitate movement
- **Two rear foldable handles** help saving space and facilitate lifting in any situation.
- **Pivoting front wheels mounted on a reinforced frame** increase stability and optimize movement
- **Foldable footrest** increases patient stability
- **Soft backrest** to increase patient comfort
- **Rear brakes** easy to spot and activate
- **The locking system of the handles** increases safety during stair descent or transport
- **Folded dimensions 980 x 510 x 310 mm**

- Length with handles and slide retracted **750 mm**
- Length sitting position with handles extracted **1170 mm**
- Width **510 mm**
- Height with headrest retracted **980 mm**
- Seat height **480 mm**
- Seat dimensions **390 x 390 mm**
- Backrest height from the seat **450 mm**
- Backrest width **450 mm**
- Front wheels **Ø100 mm**
- Rear wheels **Ø150 mm**
- Materials **Al, Steel, PU, PVC**
- Weight **14,5 kg**