



Product Specification – Reliance Medical Infant Warming Mattress

Product Information

REF:	7521
Description:	Reliance Medical Infant Warming Mattress
Manufacturer:	Danyang Rapid Aid Hot and Cold Therapy Products Co., LTD
Address:	Glasses industrial Park, Danfu Road, Danyang City, Jiangsu Province, 212300, P.R. China
EC Rep:	Lotus NL B.V
EU Importer:	Reliance Medical (Ireland)
UK Responsible Person:	Lotus Global Co, LTD
Characteristics:	Provides warmth during transportation of an infant within or between hospitals. It is recommended for full-term infants. Activated Temperature from 24°C(75°F) start: 40°C(104°F) to 37°C(98°F) after 1 hour. 400mm x 170mm
Shelf life:	5 Years from date of manufacture
Storage Conditions:	Store at room temperature. Protect from freezing and excessive heat.



Barcode Information

UNIT	QUANTITY	BARCODE	DIMENSIONS	WEIGHT
Single	1	5032483001040	230mm x 230mm x 50mm	1.5kg
Case	6	15032483001047	240mm x 240mm x 320mm	10kg

The information contained in this product specification is accurate for the above revision date.
 However, components are subject to change without notice, when another specification will be issued.