



CAT III PPE
(EU) 2016/425

AQL
1.5

EN
374

EN
455



Vinyl Examination Gloves - Blue

High quality medical grade vinyl gloves, providing a barrier against bacteria, fungi & viruses. Suitable for contact with all food types. Textured fingertips for enhanced grip.



S
SMALL

M
MEDIUM

L
LARGE

XL
EXTRA
LARGE



MediPro Vinyl Examination Gloves are a cost-effective alternative to latex and nitrile gloves. Ideal gloves for industries requiring frequent glove changes.

MD

CE
2777

A comfortable user experience

- Beaded cuff prevent roll-down.
- Provides some chemical protection.
- Protection against bacteria & viruses.
- Textured fingertips to enhance grip.
- Supplied in boxes of 100.

High quality with protection and comfort

- Latex and powder free - tested for skin irritation.
- Safe for contact with all food types to EN 1186.
- Chemical grading to EN374 1-5, TYPE B, K, P, T.
- Medical grade to EN 455 1-4
- 1.5 AQL rating.



Vinyl Examination Gloves



Recommended for:

- Health & primary care settings.
- Janitorial & industrial cleaning.
- Dental & beauty treatment settings.
- Food preparation areas inc canteens.

Ordering information

Please quote size when ordering.

| SKU | BLUE-VINYL-(SIZE) | Description | Blue Vinyl Examination Gloves |
|-----|-------------------|-------------|-------------------------------|
|-----|-------------------|-------------|-------------------------------|

Standards

Designed and manufactured in conformity with Medical Devices Directive (EU) 2017/745, EN 455 and EN 374 PPE Regulations (EU) 2016/425 for Category III PPE. Certified Class I Medical Device. Notified body: SATRA Technology Europe Limited. CE 2777. Non-sterile gloves.

Specifications

| | |
|-------------------------------|--------------------------|
| Material | Polyvinyl Chloride (PVC) |
| Grade | Medical Grade |
| Colour | Blue |
| Cuff | Beaded cuff |
| Powder content | Power free |
| Surface finish | Textured fingertips |
| Freedom from holes (level 1) | AQL 1.5 |
| Gloves length (mm/inches) min | 240mm / 9.5inch |
| Application temperature | ≤ 40°C |
| Shelf life | 5 years |

| Glove size | Palm (thickness) | Finger (thickness) | Length (mm) | Width (mm) | Weight (gram) |
|-------------|------------------|--------------------|-------------|------------|---------------|
| Small | 0.07mm ± 0.02 | 0.08mm ± 0.02 | ≥ 240mm | 84 ± 5 | 4 ± 0.3 |
| Medium | 0.07mm ± 0.02 | 0.08mm ± 0.02 | ≥ 240mm | 95 ± 5 | 4.5 ± 0.3 |
| Large | 0.07mm ± 0.02 | 0.08mm ± 0.02 | ≥ 240mm | 105 ± 5 | 5 ± 0.3 |
| Extra large | 0.07mm ± 0.02 | 0.08mm ± 0.02 | ≥ 240mm | ≥ 115 ± 5 | 5.5 ± 0.3 |

SKUs

- BLUE-VINYL-SMALL.
- BLUE-VINYL-MEDIUM.
- BLUE-VINYL-LARGE.
- BLUE-VINYL-EX-LARGE.

Key details:

- Ideal for industries where frequent glove changes are needed.
- Blue vinyl gloves are perfect for the food production sector.
- Cost-effective alternative to latex and nitrile gloves
- Latex & powder-free gloves
- Tested against skin irritation.
- Cat III PPE.
- Medical Graded gloves.
- Protects from bacteria, fungi & viruses.



www.medipro.eu