


Product Code: 85568

## H-F Antidote First Aid Kit



# In View



## Do you require HF Gel? Have you completed a risk assessment?

Due to rising sales of hybrid and electric cars increasing by approximately 60% in the last 12 months, the risk associated with HF burns has also increased. Therefore preventative measures must be taken to reduce this risk.

It has been concluded via safety assessments, that electrolyte liquid spills from Lithium-Ion batteries react with air and water to produce hydrofluoric acid, an extremely hazardous liquid. If the skin is exposed to hydrofluoric acid, it can penetrate tissue and bones, removing calcium from them and vital organs which can be lethal. Normal sterile wash products WILL NOT neutralise hydrofluoric acid, therefore HF Antidote Gel and only HF Antidote Gel, should be applied within 1 minute of exposure.

The HF first aid kits are supplied with 4 x 25g tubes of HF Gel that can be conveniently carried upon your person. The kit also includes 2 pairs of gloves for applying the gel to the patient.



### Occupational Health and Safety

These guidelines should form part of a site specific approach to the prevention and treatment of HF burns.



- All persons using HF should have appropriate personal protective equipment (PPE).
- There should be immediate access to high flow water for decontamination and easy access to calcium gluconate gel in areas where HF exposures may occur.
- All persons working on a HF plant should receive education on the spectrum of HF burns from life threatening to delayed burns.
- Workers should be provided with calcium gluconate gel to take home in the event of delayed symptoms from low concentration burns. They should present to hospital if they are concerned or have used the gel at home.

<http://www.npis.org/HFGuidance.pdf>

#### HF Antidote Gel Product Range

Product Code	Product Description	Size	Type	UoM	Case Qty
85568	HF Antidote First Aid Kit With Wall Bracket	-	Kit	Each	3
90889	25g HF Antidote Gel	25g	Tube	Each	3
86164	25g HF Antidote Gel	25g	Tube	Box of 12	1

#### Contact

T: 0845 230 2092  
E: [sales@crestmedical.co.uk](mailto:sales@crestmedical.co.uk)  
Crest Medical, 3 Chesford Grange,  
Woolston, Warrington, WA1 4RQ

Twitter: @CrestMedicalLTD

**To discuss HF Antidote gel or the full Burn Care range  
please call your dedicated sales account manager  
or the customer services team on 0845 230 2092.**