



CAT III PPE
(EU) 2016/425

AQL
1.5

EN
374

EN
455



Vinyl Examination Gloves - Clear

High quality medical grade vinyl gloves, providing a barrier against bacteria, fungi & viruses. Suitable for contact with all food types. Textured fingertips for enhanced grip.



MediPro Vinyl Examination Gloves are a cost-effective alternative to latex and nitrile gloves. Ideal gloves for industries requiring frequent glove changes.



A comfortable user experience

- Beaded cuff prevent roll-down.
- Provides some chemical protection.
- Protection against bacteria & viruses.
- Textured fingertips to enhance grip.
- Supplied in boxes of 100.

High quality with protection and comfort

- Latex and powder free - tested for skin irritation.
- Safe for contact with all food types to EN 1186.
- Chemical grading to EN374 1-5, TYPE B, K, P, T.
- Medical grade to EN 455 1-4
- 1.5 AQL rating.



Vinyl Examination Gloves



Recommended for:

- Health & primary care settings.
- Janitorial & industrial cleaning.
- Dental & beauty treatment settings.
- Food preparation areas inc canteens.

Ordering information

Please quote size when ordering.

SKU	CLEAR-VINYL-(SIZE)	Description	Clear Vinyl Examination Gloves
-----	--------------------	-------------	--------------------------------

Standards

Designed and manufactured in conformity with Medical Devices Directive (EU) 2017/745, EN 455 and EN 374 PPE Regulations (EU) 2016/425 for Category III PPE. Certified Class I Medical Device. Notified body: SATRA Technology Europe Limited. CE 2777. Non-sterile gloves.

Specifications

Material	Polyvinyl Chloride (PVC)
Grade	Medical Grade
Colour	Clear
Cuff	Beaded cuff
Powder content	Power free
Surface finish	Textured fingertips
Freedom from holes (level 1)	AQL 1.5
Gloves length (mm/inches) min	240mm / 9.5inch
Application temperature	≤ 40°C
Shelf life	5 years

Glove size	Palm (thickness)	Finger (thickness)	Length (mm)	Width (mm)	Weight (gram)
Small	0.07mm ± 0.02	0.08mm ± 0.02	≥ 240mm	84 ± 5	4 ± 0.3
Medium	0.07mm ± 0.02	0.08mm ± 0.02	≥ 240mm	95 ± 5	4.5 ± 0.3
Large	0.07mm ± 0.02	0.08mm ± 0.02	≥ 240mm	105 ± 5	5 ± 0.3
Extra large	0.07mm ± 0.02	0.08mm ± 0.02	≥ 240mm	≥ 115 ± 5	5.5 ± 0.3

SKUs

- CLEAR-VINYL-SMALL.
- CLEAR-VINYL-MEDIUM.
- CLEAR-VINYL-LARGE.
- CLEAR-VINYL-EX-LARGE.

Key details:

- Ideal for industries where frequent glove changes are needed.
- Clear vinyl gloves are perfect for the cleaning & hygiene sector.
- Cost-effective alternative to latex and nitrile gloves
- Latex & powder-free gloves
- Tested against skin irritation.
- Cat III PPE.
- Medical Graded gloves.
- Protects from bacteria, fungi & viruses.



www.medipro.eu