

PRODUCT SPECIFICATION

Product name	Nitrile Examination Gloves – Ultrasensitive
Range	Premier Sensitive
Type	Powder free, Non-Sterile
Marking	Product name, Premier Logo, Size, Code Ref, Lot Number, Expiry date, Latex Free symbol, non-sterile symbol, CE Mark, and single use symbols
Shape	Standard cuff length and ambidextrous.
Material	Nitrile
Colour	Blue



Primary use	A single use, non-sterile examination glove, intended to be used for medical purposes to prevent contamination between patient and examiner.
--------------------	--

Ref	Size	Length (mm)	Palm width (mm)	Thickness (mm) (minimum at palm)	Residual Powder (mg/glove)	Force at Break (N)
P2740	S	≥240	80 ± 10	0.07	< 2	≥ 6.0 (during shelf life) ≥ 6.0 (after challenge testing)
P2741	M	≥240	95 ± 10	0.07	< 2	≥ 6.0 (during shelf life) ≥ 6.0 (after challenge testing)
P2742	L	≥240	110 ± 10	0.07	< 2	≥ 6.0 (during shelf life) ≥ 6.0 (after challenge testing)
P2743	XL	≥240	≥ 110	0.07	< 2	≥ 6.0 (during shelf life) ≥ 6.0 (after challenge testing)

Specifications according to EN 455

Packaging	Dispenser box S; M; L; XL.	200 Gloves per dispenser
	Outer Carton	2000 (10 Dispenser Boxes of 200 Gloves)

Quality certification	ISO 13485
Conformity to Regulations	Class I – European Medical Device Regulation 2017/745
Conformity to standards	EN 455 Medical Gloves for single use. Parts 1-4 AQL 1,5
Storage conditions	Store in a cool dry place away from heat and direct sunlight
Shelf life	Refer to product packaging for date of expiry
Additional information	Finger textured to provide excellent reliability.

Medical Device and Reporting System:

Incidents caused by this device must be reported immediately to our offices.

REVISION:	5.0	DATE	12/06/2020	AUTHOR	A Gidley
------------------	-----	-------------	------------	---------------	----------

The information provided within this Data Sheet is correct as of last revision but it is subject to change without notice.

365 Healthcare
T: +44 (0) 1530 278 111
W: www.365healthcare.com

Unit 6, Delta Park Industrial Estate
Millmarsh Lane, London
Enfield, EN3 7QJ

365 Healthcare is a trading division of Bunzl Retail & Healthcare Supplies
Limited Registered Office: York House, 45 Seymour Street, London W1H 7JT. Registered in England No. 6253