



Test Report

No.: SHHG1512047994MD

Date: DEC. 09, 2015

Page: 1 of 3

SHANDONG INTCO MEDICAL PRODUCTS CO., LTD

QIWANG ROAD, NAOSHAN INDUSTRY PARK, QINGZHOU, SHANDONG, CHINA

The following sample(s) was/were submitted and identified by the client as:

Sample Description : NITRILE GLOVES
SGS Ref. No. : QDHG1512005952OT
Sample Receiving Date : DEC. 03, 2015
Testing Period : DEC. 03, 2015 TO DEC. 09, 2015
Test Performed : SELECTED TEST(S) AS REQUESTED BY APPLICANT
Test Requested : POWDER (EN 455-3-2006 MEDICAL GLOVES FOR
SINGLE USE—PART 3:REQUIREMENTS AND TESTING
FOR BIOLOGICAL EVALUATION CLAUSE 4.4)
Test Result(s) : FOR FURTHER DETAILS, PLEASE REFER TO THE
FOLLOWING PAGE(S)
Conclusion : THE SUBMITTED SAMPLE MET THE TEST
REQUIREMENT.

Signed for and on behalf of
SGS-CSTC Standards Technical Services (Shanghai) Co., Ltd.

Vincent Feng
Technical Manager



This document is issued by the Company subject to its General Conditions of Service printed overleaf, available on request or accessible at <http://www.sgs.com/en/Terms-and-Conditions.aspx> and, for electronic format documents, subject to Terms and Conditions for Electronic Documents at <http://www.sgs.com/en/Terms-and-Conditions/Terms-e-Documents.aspx>. Attention is drawn to the limitation of liability, indemnification and jurisdiction issues defined therein. Any holder of this document is advised that information contained hereon reflects the Company's findings at the time of its intervention only and within the limits of Client's instructions, if any. The Company's sole responsibility is to its Client and this document does not exonerate parties to a transaction from exercising all their rights and obligations under the transaction documents. This document cannot be reproduced except in full, without prior written approval of the Company. Any unauthorized alteration, forgery or falsification of the content or appearance of this document is unlawful and offenders may be prosecuted to the fullest extent of the law. Unless otherwise stated the results shown in this test report refer only to the sample(s) tested and such sample(s) are retained for 30 days only.

Attention: To check the authenticity of testing / inspection report & certificate, please contact us at telephone: (86-755) 8307 1443, or email: CN.Doccheck@sgs.com

3rd Building, No.889, Yishan Road, Xuhui District Shanghai, China 200233

t (86-21) 61402666*2013

f (86-21) 54500353

www.sgs.com

中国·上海·徐汇区宜山路889号3号楼 邮编: 200233

t (86-21) 61402666*2013

f (86-21) 54500353

e sgs.china@sgs.com

Member of the SGS Group (SGS SA)



Test Report

No.: SHHG1512047994MD

Date: DEC. 09, 2015

Page: 2 of 3

Test Conducted:

EN 455-3-2006 Medical gloves for single use—Part 3:Requirements and testing for biological evaluation

Number of test sample	:	5 Pieces
Finishes of gloves	:	Powdered-free gloves other than surgeon's gloves
Defects observed before testing	:	No defects
Test Result	:	Pass

Clause	Test Items	Result	Note
4.4	Powder	Pass	#1

Note:

1. Test according to EN ISO 21171-2006.
2. The quantity of powder was 0.2mg<2mg.



This document is issued by the Company subject to its General Conditions of Service printed overleaf, available on request or accessible at <http://www.sgs.com/en/Terms-and-Conditions.aspx> and, for electronic format documents, subject to Terms and Conditions for Electronic Documents at <http://www.sgs.com/en/Terms-and-Conditions/Terms-e-Documents.aspx>. Attention is drawn to the limitation of liability, indemnification and jurisdiction issues defined therein. Any holder of this document is advised that information contained hereon reflects the Company's findings at the time of its intervention only and within the limits of Client's instructions, if any. The Company's sole responsibility is to its Client and this document does not exonerate parties to a transaction from exercising all their rights and obligations under the transaction documents. This document cannot be reproduced except in full, without prior written approval of the Company. Any unauthorized alteration, forgery or falsification of the content or appearance of this document is unlawful and offenders may be prosecuted to the fullest extent of the law. Unless otherwise stated the results shown in this test report refer only to the sample(s) tested and such sample(s) are retained for 30 days only.

Attention: To check the authenticity of testing /inspection report & certificate, please contact us at telephone: (86-755) 8307 1443, or email: CN.Doccheck@sgs.com

3rd Building, No.889, Yishan Road, Xuhui District Shanghai, China 200233

t (86-21) 61402666*2013

f (86-21) 54500353

www.sgsgroup.com.cn

中国·上海·徐汇区宜山路889号3号楼 邮编: 200233

t (86-21) 61402666*2013

f (86-21) 54500353

e sgs.china@sgs.com

Member of the SGS Group (SGS SA)

Test Report

No.: SHHG1512047994MD

Date: DEC. 09, 2015

Page: 3 of 3

Sample Photo:

Received sample



SGS authenticate the photo on original report only

End of Report