



Single Use Dermal Biopsy Punch

**Ergonomically** designed  
with **seamless, precision**  
ground **stainless steel** blade

- ✔ Single use instrument  
- always clean, sharp and sterile
- ✔ Seamless precision ground  
high quality stainless steel blade
- ✔ Robust comfortable handle  
to improve grip and control
- ✔ Depth stop feature ensures optimum  
diagnostic sample depth
- ✔ Biocompatibility tested handle  
polymer to ISO 10993
- ✔ Available in a range of sizes - Ø 3/4/5mm
- ✔ Product and packaging marked with the  
blade size for ease of identification



RHC Code	Description	Unit of Sale
8180	Single Use Dermal Biopsy Punch - Ø 3mm	10 Packs of 10
8181	Single Use Dermal Biopsy Punch - Ø 4mm	10 Packs of 10
8182	Single Use Dermal Biopsy Punch - Ø 5mm	10 Packs of 10



To request samples, visit [www.robinsonhealthcare.com](http://www.robinsonhealthcare.com)