



SKU: MP-BP-BPM

Automatic Upper Arm Blood Pressure Monitor

Take your blood pressure quickly and comfortably in your own home. Accurate blood pressure measurements at the push of a button.



The MediPro home use blood pressure monitor offers the best value for your money, providing simple to understand instructions and accurate blood pressure readings.



Supporting a comfortable user experience

- Obtain blood pressure measurements with speed.
- Comfortably check blood pressure with accuracy.
- Automatic cuff inflation at the push of a button.
- Supplied with a comprehensive guidance leaflet.
- Large LCD display for easy to read results.

Clinical accuracy for your home

- Detects for Irregular Heart Beat (IHB).
- WHO classification for high blood pressure.
- Class IIa Medical Device, CE 0123.
- British & Irish Hypertension Society (BIHS) approved.



supplied with a wide ranging cuff 22-42cm

Automatic Upper Arm Blood Pressure Monitor



Ideal for home monitoring

- High blood pressure.
- Low blood pressure.
- Irregular Heart Beat (IHB).
- Temporary home monitoring.

Ordering information

SKU	MP-BP-BPM	Description	Automatic Upper Arm Blood Pressure Monitor
-----	-----------	-------------	--

Standards

Designed and manufactured in conformity with Medical Devices Directive 93/42/EEC and ISO 13486 Medical Devices Quality Management Systems. Class IIa Medical Device. Notified body: TUV SUD Products Service GmbH. CE 0123.

Specifications

Size	105 L x 120 W x 55 H mm
Weight	290g (without batteries)
Measuring method	Oscillometric measurement
Cuff pressure	290 mmHg
Measuring range	40 to 180 mmHg (diastolic pressure) 60 to 260 mmHg (systolic pressure) 40 to 160 beats/minute (pulse rate)
Measuring accuracy	±3 mmHg for static pressure ±5% of the reading for the pulse rate
Inflation	Automatic by the pump
Rapid deflation	Automatic electronic valve
Batteries	x4 alkaline .5v batteries (AA, R6P or LR6)
Adapter	6v 600mA (optional part)
Memory	2 users with 90 sets of memory each
Supplied cuff size	Wide range cuff 22-42cm
Number of measurements	Approx. 700 times (LR6 alkaline batteries) Approx. 200 times (R6P manganese batteries)

What's included

- x1 blood pressure monitor.
- x1 tube and connector.
- x1 wide ranging cuff.
- x1 instructions for use.

Key facts

- Perfect for home use.
- BIHS approved upper arm device.
- 90 sets of memory to record readings.
- Takes an average value from three blood pressure readings.
- Screens for Irregular Heart Beat (IHB).
- Takes blood pressure readings against WHO blood pressure classification.
- Wide ranging cuff - 22-42cm.
- Powered by x4 AA batteries (not included).
- Automatic cuff inflation.



www.medipro.eu