

Otoscope & Ophthalmoscope Desk Set - Rechargeable

A compact, modern and robust design ideal for busy, fast paced environments such as wards, surgeries, practices and consulting rooms.



This set offers the best value for your money, providing all diagnostic needs in one handy package giving up to 20 hours of operational use.

Supporting a safe and comfortable patient experience

- Quick and simple to operate with minimal patient skin contact.
- Obtain patient observations with speed.
- Easy to clean handles to help minimise cross-contamination.
- Over 20 hours of continuous operational use.

Enhancing staff satisfaction and productivity

- Combined ophthalmoscope and otoscope desk set.
- Supplied with a double charging base, allowing both handles to be charged at the same time by usb or power adapter.
- Powered by Optic SMT LED technology – 2000 hours of use.
- Value adding accessories and consumables available.



Otoscope & Ophthalmoscope Desk Set



Ideal care environments

- Emergency care departments
- GP surgeries & practices
- School health services
- Healthcare clinics

Ordering information

SKU	MP-DESK-SET	Description	Otoscope & Ophthalmoscope Desk Set
-----	-------------	-------------	------------------------------------

Standards

Designed and manufactured in conformity with Medical Devices Regulation (EU) 2017/745 and ISO 13486 Medical Devices Quality Management Systems. Class 1 Medical Device.
Notified body: TUV SUD Products Service GmbH.

Specifications

PART	OTOSCOPE	OPHTHALMOSCOPE
Illumination	Fibre optic illumination	Direct illumination
Magnification	x3	x15
Bulb model	B109L	B206L
Bulb electrical parameter	LED 3.0v DC 60mA	
Luminance	24 lumen	14 lumen
Colour temp	3200K or 5700K +5%	3500K +5%
Tip sizes	2.4mm to 6mm	N/A
Battery	B116 Ni-Mh 2.4v 1600mAh rechargeable batteries	
Charging dock	Input: DC 5.0v 1A	Output: DC 3.0v 2x250mA

What's included

- x1 otoscope.
- x1 ophthalmoscope.
- x1 double desk recharging point.
- x1 power adapter (DC 5V).
- x1 inflation bulb.
- x8 ear specula, sizes: 2.4mm, 3mm, 4mm, 6mm and 8mm
- x1 instructions for use.

Key facts

- 5 different apertures – wide, micro, fixation, hemispot & red-free.
- 15x magnification on ophthalmoscope.
- Dioptic range: -20d to +20d.
- Six levels of brightness.
- 3x magnification on otoscope.
- Up to 20 hours of operational time.