

Pocket LED Otoscope Set

A compact, modern and robust design ideal for busy, fast paced environments such as wards, surgeries, practices and consulting rooms.



The pocket otoscope set offers the best value for your money, providing all diagnostic needs in one handy package giving up to 25 hours of operational use.

Supporting a safe and comfortable patient experience

- Quick and simple to operate with minimal patient skin contact.
- Obtain patient ear observations with speed & accuracy.
- Easy to clean handles to help minimise cross-contamination.
- Over 25 hours of continuous operational use.

Enhancing staff satisfaction and productivity

- The perfect pocket set for practitioners on the move.
- Powered by AA 2.4v batteries or rechargeable 116 batteries.
- Powered by Optic SMT LED technology – 2000 hours of use.
- Value adding accessories and consumables available.



Pocket LED Otoscope Set



Ideal care environments

- Emergency care departments
- GP surgeries & practices
- School health services
- Healthcare clinics

Ordering information

SKU	MP-POCKET-OT	Description	Pocket LED Otoscope Set
-----	--------------	-------------	-------------------------

Standards

Designed and manufactured in conformity with Medical Devices Regulation (EU) 2017/745 and ISO 13486 Medical Devices Quality Management Systems. Class 1 Medical Device. Notified body: TUV SUD Products Service GmbH.

Specifications

PART	OTOSCOPE	OPHTHALMOSCOPE
Illumination	Fibre optic illumination	Direct illumination
Magnification	x3	x15
Bulb model	B109N	B206
Bulb electrical parameter	LED 3.0v DC 60mA	
Luminance	24 lumen	14 lumen
Colour temp	3200K or 5700K +5%	3500K +5%
Tip sizes	2.4mm to 6mm	N/A
Battery	B116 Ni-Mh 2.4v 1600mAh rechargeable batteries	
Charging dock	Input: DC 5.0v 1A	Output: DC 3.0v 2x250mA

What's included

- x1 otoscope head.
- x1 handle.
- x1 heavy duty carry case.
- x8 ear specula, sizes: 2.4mm, 3mm, 4mm, 6mm and 8mm
- x1 instructions for use.

Key facts

- Built-in insufflation bulb socket for pneumatic testing.
- Advanced fibre optic illumination
- 3x magnification viewing window.
- Warm white 3200k light temperature.
- Six levels of brightness.
- Three-minute automatic power off.
- Up to 25 hours of operational time.
- Powered by AA 2.4v batteries (not supplied).