

## Otoscope & Ophthalmoscope Pocket Set

A compact, modern and robust design ideal for busy, fast paced environments such as wards, surgeries, practices and consulting rooms.



**This pocket set offers the best value for your money, providing all diagnostic needs in one handy package giving up to 15 hours of operational use.**

### Supporting a safe and comfortable patient experience

- Quick and simple to operate with minimal patient skin contact.
- Obtain patient observations with speed.
- Easy to clean handles to help minimise cross-contamination.
- Over 15 hours of continuous operational use.

### Enhancing staff satisfaction and productivity

- Combined ophthalmoscope and otoscope pocket set.
- Powered by AA 2.4v batteries or rechargeable 116 batteries.
- Powered by Optic SMT LED technology – 2000 hours of use.
- Value adding accessories and consumables available.



## Otoscope & Ophthalmoscope Pocket Set



### Ideal care environments

- Emergency care departments
- GP surgeries & practices
- School health services
- Healthcare clinics

### Ordering information

SKU	MP-POCKET-SET	Description	Otoscope & Ophthalmoscope Pocket Set
-----	---------------	-------------	--------------------------------------

### Standards

Designed and manufactured in conformity with Medical Devices Regulation (EU) 2017/745 and ISO 13486 Medical Devices Quality Management Systems. Class 1 Medical Device. Notified body: TUV SUD Products Service GmbH.

### Specifications

PART	OTOSCOPE	OPHTHALMOSCOPE
Illumination	Fibre optic illumination	Direct illumination
Magnification	x3	x15
Bulb model	B109N	B206
Bulb electrical parameter	LED 3.0v DC 60mA	
Luminance	24 lumen	14 lumen
Colour temp	3200K or 5700K +5%	3500K +5%
Tip sizes	2.4mm to 6mm	N/A
Battery	B116 Ni-Mh 2.4v 1600mAh rechargeable batteries	
Charging dock	Input: DC 5.0v 1A	Output: DC 3.0v 2x250mA

### What's included

- x1 otoscope head.
- x1 ophthalmoscope head.
- x1 handle.
- x1 heavy duty carry case.
- x8 ear specula, sizes: 2.4mm, 3mm, 4mm, 6mm and 8mm
- x1 instructions for use.

### Key facts

- 5 different apertures – wide, micro, fixation, hemispot & red-free.
- 15x magnification on ophthalmoscope.
- Dioptic range: -20d to +20d.
- Six levels of brightness.
- 3x magnification on otoscope.
- Up to 15 hours of operational time.
- Powered by AA 2.4v batteries (not supplied).