

MARSDEN

M-605

Marsden Hoist Weighing Attachment

Lightweight and portable for care homes and the community

Designed for larger hoists

Class III and MDD approved

Accurate to 100g

www.marsden-weighing.co.uk



Marsden M-605:

A versatile hoist attachment for flexible patient weighing

The Marsden M-605 is a lightweight and portable hoist weighing attachment designed for use within care homes and the community.

Class III and MDD approved hoist attachment

Accurate weighing for bed-bound patients

The M-605 is an easy to use Class III hoist weighing attachment designed for use on larger hoists with coat hanger shaped spreader bars. It is wider than our standard M-600 hoist weighing attachment.

This hoist weighing scale is designed to clip onto your coat-hanger shaped spreader bar easily. The clear, digital display shows the weight of the patient to the nearest 100g, up to a capacity of 200kg.

Hoist weighing attachments are now required by MHRA to be approved under the Medical Devices Directive (93/42/EEC). The Marsden M-600 and M-605 are now MDD approved making them the only hoist weighing attachments that are legal for use in the UK.

Easy to use function keys





Optional carry case



Suits most 'coathanger' sling supports



MARSDEN
WEIGHING

Record weight data easily

As well as Hold and Tare functions for clear and accurate weight readings, this scale is available with Bluetooth.

Bluetooth connectivity allows the weight data to be instantly, wirelessly transferred from the scale to a central database, PC or handheld device.

A printer can also be connected to this scale to provide a printed record of the patient's weight, height and BMI reading along with the date and time the readings were taken.

Built to last

Like other scales in the Marsden range the M-605 is built to last. It comes with a two year warranty, and for extra peace of mind it can be kept accurate and reliable with a Marsden service contract.

Product Specifications

Capacity: 200kg
Graduations: 100g
Weight of scale: 6.7kg
Power Supply: Mains/Rechargeable Battery
Dimensions: 730 x 275 x 160mm

Optional Extras:

Carry Case
Thermal Printer



Accuracy Assured

Tel: 01709 364296 / 0800 169 2775

Fax: 01709 364293

E-mail: sales@marsdengroup.co.uk

www.marsden-weighing.co.uk

 MarsdenWeighing

 @MarsdenWeighing

 www.youtube.com/MarsdenWeighing

 <http://pinterest.com/marsdenweighing>