



SKU: MP-PULSE-OXI

MediPro Pulse Oximeter

A non-invasive medical device for the accurate spot monitoring of blood oxygen levels (SpO₂) and heart pulse rate for adults, adolescents and children aged 4 and over.



A best value for money pulse oximeter providing up to 30 hours of continuous use. Ideal for use within hospital settings, doctors surgeries, care facilities and home use.



Low power mode



Rotating screen



8 second power off

Supporting a safe and comfortable user experience

- Quick and simple - simply clips to the fingertip.
- Obtains pulse rate and blood oxygen levels with speed.
- Perfect for use with adults, adolescents and children over 4.
- Non-invasive and painless method of measurement.

Enhancing staff satisfaction and productivity

- Large and 360° rotatable OLED display for ease of use.
- Automatic power-off after 8 seconds to preserve battery life.
- Lightweight, durable and portable - perfect for patient rounds.
- IP22 rated against water ingress and Type BF rating for protection against electrical shocks.



MediPro Pulse Oximeter



Ideal care environments

- Emergency care departments
- GP surgeries & practices
- Hospitals & care facilities
- Home use for self-monitoring

Ordering information

SKU	MP-PULSE-OXI	Description	MediPro Pulse Oximeter
-----	--------------	-------------	------------------------

Standards

Designed and manufactured in conformity with Medical Devices Regulation (EU) 2017/745 and ISO 13485 Medical Devices Quality Management Systems. Class IIa Medical Device. Notified body: TUV SUD Products Service GmbH. CE 0123.

Specifications

PART	INFORMATION
Display	360 degree rotating OLED display
Age Range	Adults, adolescents & children aged 4+
Operation	Continuous mode of operation: 30 hours
SpO2 Range	70% to 90%
SpO2 Accuracy	±2% stage 70% to 99% unspecified ≤70%
SpO2 Resolution	±1%
Pulse Range	30 BPM to 240 BPM
Pulse Accuracy	±2 BPM or ±2%
Power Requirements	Two AAA 1.5v Alkaline Batteries (not included)
Compliance & Ratings	IP22 rated against water ingress TYPE BF protection against electric shocks Class IIa Medical Device CE 0123

What's included

- x1 fingertip pulse oximeter.
- x1 instruction manual.
- x1 wrist lanyard.

Key pulse oximetry facts

- Gives insight into how well someone is breathing (SpO2) along with how fast the heart is beating (BPM).
- Indication to any respiratory or cardiac issues like Chronic Obstructive Pulmonary Disorder (COPD).
- Detecting silent Hypoxia and screening for Irregular Heart Beat (IHB).

Handy health check tool for:

- Heart issues
- COPD
- Lung cancer
- Anemia
- Asthma



www.medipro.eu